

SphK2 (N-19): sc-22703

BACKGROUND

Sphingosine kinase (SphK or SphK1) is a key enzyme catalyzing the phosphorylation of sphingosine to form sphingosine 1-phosphate (SPP or S1P). SPP is a bioactive lipid that exerts multiple biological effects in a large variety of cell types, acting as either an intracellular messenger or an extracellular ligand coupled to Edg-family receptors. Competitive inhibitors of SphK1 block formation of SPP and selectively inhibit cellular proliferation induced by a variety of factors. One potent inhibitor of SphK1 activity is DMS (N,N-dimethylsphingosine). SPP/SphK1 has been implicated as a signaling pathway that regulates diverse cellular functions, including cell growth, proliferation and survival. Specifically, SphK1 is involved in the signaling pathway(s) that protects human hepatocytes from the apoptotic action of TNF α . Furthermore, SPP/SphK1 may play an important role in neuronal survival by regulating activation of SAPKs and caspases. SphK1 is widely expressed with highest levels in adult liver, kidney, heart and skeletal muscle; however, activation of SphK1 disengages cells from their liver-specific phenotype. SphK1 is highly homologous with SphK2, another member of a growing class of sphingolipid kinases. Expression of SphK2 mRNA exhibits a markedly different tissue distribution than that of SphK1 and appears at a later stage in embryonic development.

REFERENCES

- Xia, P., et al. 2000. An oncogenic role of sphingosine kinase. *Curr. Biol.* 10: 1527-1530.
- Liu, H., et al. 2000. Molecular cloning and functional characterization of a novel mammalian sphingosine kinase type 2 isoform. *J. Biol. Chem.* 275: 19513-19520.
- Osawa, Y., et al. 2001. Sphingosine kinase regulates hepatoma cell differentiation: roles of hepatocyte nuclear factor and retinoid receptor. *Biochem. Biophys. Res. Commun.* 286: 673-677.
- Edsall, L.C., et al. 2001. Sphingosine kinase expression regulates apoptosis and caspase activation in PC12 cells. *J. Neurochem.* 76: 1573-1584.
- Osawa, Y., et al. 2001. TNF- α -induced sphingosine 1-phosphate inhibits apoptosis through a phosphatidylinositol 3-kinase/Akt pathway in human hepatocytes. *J. Immunol.* 167: 173-180.
- Hayashi, S., et al. 2002. Identification and characterization of RPK118, a novel sphingosine kinase-1-binding protein. *J. Biol. Chem.* 277: 33319-33324.
- Meacci, E., et al. 2002. Sphingosine 1-phosphate evokes calcium signals in C2C12 myoblasts via Edg3 and Edg5 receptors. *Biochem. J.* 362: 349-357.
- SWISS-PROT/TrEMBL (Q9NYA1). World Wide Web URL: <http://www.expasy.ch/sprot/sprot-top.html>

CHROMOSOMAL LOCATION

Genetic locus: SPHK2 (human) mapping to 19q13.33; Sphk2 (mouse) mapping to 7 B4.

SOURCE

SphK2 (N-19) is an affinity purified goat polyclonal antibody raised against a peptide mapping near the N-terminus of SphK2 of human origin.

PRODUCT

Each vial contains 200 μ g IgG in 1.0 ml of PBS with < 0.1% sodium azide and 0.1% gelatin.

Blocking peptide available for competition studies, sc-22703 P, (100 μ g peptide in 0.5 ml PBS containing < 0.1% sodium azide and 0.2% BSA).

APPLICATIONS

SphK2 (N-19) is recommended for detection of SphK2 of mouse, rat and human origin by Western Blotting (starting dilution 1:200, dilution range 1:100-1:1000), immunofluorescence (starting dilution 1:50, dilution range 1:50-1:500) and solid phase ELISA (starting dilution 1:30, dilution range 1:30-1:3000).

SphK2 (N-19) is also recommended for detection of SphK2 in additional species, including equine and porcine.

Suitable for use as control antibody for SphK2 siRNA (h): sc-39225, SphK2 siRNA (m): sc-39226, SphK2 shRNA Plasmid (h): sc-39225-SH, SphK2 shRNA Plasmid (m): sc-39226-SH, SphK2 shRNA (h) Lentiviral Particles: sc-39225-V and SphK2 shRNA (m) Lentiviral Particles: sc-39226-V.

Molecular Weight of SphK2: 70 kDa.

Positive Controls: mouse brain extract: sc-2253.

RECOMMENDED SECONDARY REAGENTS

To ensure optimal results, the following support (secondary) reagents are recommended: 1) Western Blotting: use donkey anti-goat IgG-HRP: sc-2020 (dilution range: 1:2000-1:100,000) or Cruz Marker™ compatible donkey anti-goat IgG-HRP: sc-2033 (dilution range: 1:2000-1:5000), Cruz Marker™ Molecular Weight Standards: sc-2035, TBS Blotto A Blocking Reagent: sc-2333 and Western Blotting Luminol Reagent: sc-2048. 2) Immunofluorescence: use donkey anti-goat IgG-FITC: sc-2024 (dilution range: 1:100-1:400) or donkey anti-goat IgG-TR: sc-2783 (dilution range: 1:100-1:400) with UltraCruz™ Mounting Medium: sc-24941.

STORAGE

Store at 4° C, ****DO NOT FREEZE****. Stable for one year from the date of shipment. Non-hazardous. No MSDS required.

RESEARCH USE

For research use only, not for use in diagnostic procedures.

PROTOCOLS

See our web site at www.scbt.com or our catalog for detailed protocols and support products.



Try **SphK2 (9C5E1): sc-517192**, our highly recommended monoclonal alternative to SphK2 (N-19).