

Formin 2 (K-18): sc-22735

BACKGROUND

The temporal genetic hierarchy influencing normal limb development can deregulate and mediate mammalian developmental syndromes. In mice, the limb deformity (Id) locus influences normal limb development and gives rise to alternative mRNAs that can translate into a family of proteins known as formins. Formins play a crucial role in cytoskeletal reorganization by influencing Actin filament assembly. Formins co-localize with the Actin cytoskeleton and can translocate into the cell cytosol and into the nucleus in an HGF-dependent manner. Vertebrate nuclear formins can control polarizing activity in limb buds through establishment of a Sonic hedgehog/FGF-4 feedback loop. Deficiency mutations at the mammalian Id locus lead to profound developmental defects in limb and kidney formation. The human Formin 1 and 2 genes map to chromosome 15q13.3 and 1q43, respectively.

REFERENCES

1. Maas, R.L., Jepeal, L.I., Elfering, S.L., Holcombe, R.F., Morton, C.C., Eddy, R.L., Byers, M.G., Shows, T.B. and Leder, P. 1991. A human gene homologous to the formin gene residing at the murine limb deformity locus: chromosomal location and RFLPs. *Am. J. Hum. Gen.* 48: 687-695.
2. Wynshaw-Boris, A., Ryan, G., Deng, C.X., Chan, D.C., Jackson-Grusby, L., Larson, D., Dunmore, J.H. and Leder, P. 1997. The role of a single formin isoform in the limb and renal phenotypes of limb deformity. *Mol. Med.* 3: 372-384.
3. Zeller, R., Haramis, A.G., Zuniga, A., McGuigan, C., Dono, R., Davidson, G., Chabanis, S. and Gibson, T. 1999. Formin defines a large family of morphoregulatory genes and functions in establishment of the polarising region. *Cell Tissue Res.* 296: 85-93.
4. Tanaka, K. 2000. Formin family proteins in cytoskeletal control. *Biochem. Biophys. Res. Commun.* 267: 479-481.
5. O'Rourke, D.A., Liu, Z.X., Sellin, L., Spokes, K., Zeller, R. and Cantley, L.G. 2000. Hepatocyte growth factor induces MAPK-dependent Formin 4 translocation in renal epithelial cells. *J. Am. Soc. Nephrol.* 11: 2212-2221.
6. Leader, B. and Leder, P. 2000. Formin 2, a novel formin homology protein of the cappuccino subfamily, is highly expressed in the developing and adult central nervous system. *Mech. Dev.* 93: 221-231.
7. Sawin, K.E. 2002. Cell polarity: following formin function. *Curr. Biol.* 12: R6-R8.
8. LocusLink Report (LocusID: 2325). <http://www.ncbi.nlm.nih.gov/LocusLink/>

CHROMOSOMAL LOCATION

Genetic locus: FMN2 (human) mapping to 1q43.

SOURCE

Formin 2 (K-18) is an affinity purified goat polyclonal antibody raised against a peptide mapping within an internal region of Formin 2 of human origin.

RESEARCH USE

For research use only, not for use in diagnostic procedures.

PRODUCT

Each vial contains 200 µg IgG in 1.0 ml of PBS with < 0.1% sodium azide and 0.1% gelatin.

Blocking peptide available for competition studies, sc-22735 P, (100 µg peptide in 0.5 ml PBS containing < 0.1% sodium azide and 0.2% BSA).

APPLICATIONS

Formin 2 (K-18) is recommended for detection of Formin 2 of human origin by Western Blotting (starting dilution 1:200, dilution range 1:100-1:1000), immunofluorescence (starting dilution 1:50, dilution range 1:50-1:500) and solid phase ELISA (starting dilution 1:30, dilution range 1:30-1:3000).

Formin 2 (K-18) is also recommended for detection of Formin 2 in additional species, including equine, canine, porcine and avian.

Suitable for use as control antibody for Formin 2 siRNA (h): sc-43765, Formin 2 shRNA Plasmid (h): sc-43765-SH and Formin 2 shRNA (h) Lentiviral Particles: sc-43765-V.

Molecular Weight of Formin 2: 195 kDa.

RECOMMENDED SECONDARY REAGENTS

To ensure optimal results, the following support (secondary) reagents are recommended: 1) Western Blotting: use donkey anti-goat IgG-HRP: sc-2020 (dilution range: 1:2000-1:100,000) or Cruz Marker™ compatible donkey anti-goat IgG-HRP: sc-2033 (dilution range: 1:2000-1:5000), Cruz Marker™ Molecular Weight Standards: sc-2035, TBS Blotto A Blocking Reagent: sc-2333 and Western Blotting Luminol Reagent: sc-2048. 2) Immunofluorescence: use donkey anti-goat IgG-FITC: sc-2024 (dilution range: 1:100-1:400) or donkey anti-goat IgG-TR: sc-2783 (dilution range: 1:100-1:400) with UltraCruz™ Mounting Medium: sc-24941.

SELECT PRODUCT CITATIONS

1. Shieh, D.B., Li, R.Y., Liao, J.M., Chen, G.D. and Liou, Y.M. 2010. Effects of Genistein on β -catenin signaling and subcellular distribution of Actin-binding proteins in human umbilical CD105-positive stromal cells. *J. Cell. Physiol.* 223: 423-434.

STORAGE

Store at 4° C, ****DO NOT FREEZE****. Stable for one year from the date of shipment. Non-hazardous. No MSDS required.

PROTOCOLS

See our web site at www.scbt.com or our catalog for detailed protocols and support products.



Try **Formin 2 (C-3): sc-376787**, our highly recommended monoclonal alternative to Formin 2 (K-18).