

Rpp38 (D-18): sc-23035

BACKGROUND

Ribonuclease P (RNase P) and Ribonuclease MRP (RNase MRP) are small nuclear ribonucleoproteins (snRNPs) that act on RNA substrates *in vitro*. RNase P and RNase MRP, which accumulate in the nucleolus, have a similar RNA component and also have several protein subunits in common. RNase P, which consists of a complex of an RNA species (H1 RNA), POP1 (processing of precursor 1), POP5 (processing of precursor 5), and at least 7 Rpps (including Rpp14, Rpp29, Rpp30 and Rpp38), removes the 5' leader sequences from precursor tRNA molecules. In particular, the nucleolar-localizing RNase P catalyzes the hydrolysis of a specific phosphodiester bond in precursor tRNA to form the mature 5' end of tRNA. The structurally related RNase MRP plays an essential role in the formation of the 5' end of 5.8S rRNA. Both RNase P and RNase MRP are associated with Th/To ribonucleoproteins; Rpp30 and Rpp38 have specifically been implicated as Th autoantigens which contribute to the autoimmune disease systemic sclerosis.

REFERENCES

1. Karwan, R. 1993. RNase MRP/RNase P: a structure-function relation conserved in evolution? FEBS Lett. 319: 1-4.
2. Jarrous, N., Eder, P.S., Guerrier-Takada, C., Hoog, C. and Altman, S. 1998. Autoantigenic properties of some protein subunits of catalytically active complexes of human ribonuclease P. RNA 4: 407-417.

CHROMOSOMAL LOCATION

Genetic locus: RPP38 (human) mapping to 10p13.

SOURCE

Rpp38 (D-18) is an affinity purified goat polyclonal antibody raised against a peptide mapping within an internal region of Rpp38 of human origin.

PRODUCT

Each vial contains 200 µg IgG in 1.0 ml of PBS with < 0.1% sodium azide and 0.1% gelatin.

Blocking peptide available for competition studies, sc-23035 P, (100 µg peptide in 0.5 ml PBS containing < 0.1% sodium azide and 0.2% BSA).

APPLICATIONS

Rpp38 (D-18) is recommended for detection of Rpp38 of human origin by Western Blotting (starting dilution 1:200, dilution range 1:100-1:1000), immunoprecipitation [1-2 µg per 100-500 µg of total protein (1 ml of cell lysate)], immunofluorescence (starting dilution 1:50, dilution range 1:50-1:500) and solid phase ELISA (starting dilution 1:30, dilution range 1:30-1:3000).

Suitable for use as control antibody for Rpp38 siRNA (h): sc-106729, Rpp38 shRNA Plasmid (h): sc-106729-SH and Rpp38 shRNA (h) Lentiviral Particles: sc-106729-V.


Molecular Weight of Rpp38: 38-40 kDa.

Positive Controls: Hep G2 cell lysate: sc-2227 or PC-3 cell lysate: sc-2220.


RECOMMENDED SECONDARY REAGENTS

To ensure optimal results, the following support (secondary) reagents are recommended: 1) Western Blotting: use donkey anti-goat IgG-HRP: sc-2020 (dilution range: 1:2000-1:100,000) or Cruz Marker™ compatible donkey anti-goat IgG-HRP: sc-2033 (dilution range: 1:2000-1:5000), Cruz Marker™ Molecular Weight Standards: sc-2035, TBS Blotto A Blocking Reagent: sc-2333 and Western Blotting Luminol Reagent: sc-2048. 2) Immunoprecipitation: use Protein A/G PLUS-Agarose: sc-2003 (0.5 ml agarose/2.0 ml). 3) Immunofluorescence: use donkey anti-goat IgG-FITC: sc-2024 (dilution range: 1:100-1:400) or donkey anti-goat IgG-TR: sc-2783 (dilution range: 1:100-1:400) with UltraCruz™ Mounting Medium: sc-24941.

DATA



Rpp38 (D-18): sc-23035. Western blot analysis of Rpp38 expression in Hep G2 (A) and PC-3 (B) whole cell lysates.



Rpp38 (D-18): sc-23035. Immunofluorescence staining of formalin-fixed Hep G2 cells showing nucleolar and nuclear localization.

STORAGE

Store at 4° C, **DO NOT FREEZE**. Stable for one year from the date of shipment. Non-hazardous. No MSDS required.

RESEARCH USE

For research use only, not for use in diagnostic procedures.

PROTOCOLS

See our web site at www.scbt.com or our catalog for detailed protocols and support products.


 MONOS
Satisfaction
Guaranteed

Try **Rpp38 (A-2): sc-398112** or **Rpp38 (B-6): sc-398113**, our highly recommended monoclonal alternatives to Rpp38 (D-18).