

ACLP (N-14): sc-23051

BACKGROUND

Adipocyte enhancer binding protein 1 (AEBP 1), of which aortic carboxypeptidase-like protein (ACLP) is an isoform, is a transcriptional repressor with carboxypeptidase activity that is expressed in vascular smooth muscle cells, and at lower levels in adipose and osteoblastic cells. The ACLP isoform of AEBP 1 contains a signal peptide sequence, a lysine- and proline-rich repeating motif, a discoidin-like protein and a carboxypeptidase-like domain. AEBP 1 is secreted into the extracellular matrix and may play a role in abdominal wall development and dermal wound healing. In differentiating murine osteoblastic cells, AEBP 1 is expressed but shut off in the final calcification phase, indicating a transcriptional repressive effect on genes for bone formation. The expression of AEBP 1 is downregulated during adipogenesis and upregulated during vascular smooth muscle cell differentiation.

REFERENCES

- Ohno, I., et al. 1996. A cDNA cloning of human AEBP1 from primary cultured osteoblasts and its expression in a differentiating osteoblastic cell line. *Biochem. Biophys. Res. Commun.* 228: 411-414.
- Layne, M.D., et al. 1998. Aortic carboxypeptidase-like protein, a novel protein with discoidin and carboxypeptidase-like domains, is upregulated during vascular smooth muscle cell differentiation. *J. Biol. Chem.* 273: 15654-15660.
- Ro, H.S., et al. 2001. Gene structure and expression of the mouse adipocyte enhancer-binding protein. *Gene* 280: 123-33.
- Gagnon, A., et al. 2002. Downregulation of aortic carboxypeptidase-like protein during the early phase of 3T3-L1 adipogenesis. *Endocrinology* 143: 2478-2485.
- Entrez-Protein (NP_001120). World Wide Web URL: <http://www.ncbi.nlm.nih.gov:80/entrez>

CHROMOSOMAL LOCATION

Genetic locus: AEBP1 (human) mapping to 7p13; Aebp1 (mouse) mapping to 11 A1.

SOURCE

ACLP (N-14) is an affinity purified goat polyclonal antibody raised against a peptide mapping near the N-terminus of ACLP of human origin.

PRODUCT

Each vial contains 200 µg IgG in 1.0 ml of PBS with < 0.1% sodium azide and 0.1% gelatin.

Blocking peptide available for competition studies, sc-23051 P, (100 µg peptide in 0.5 ml PBS containing < 0.1% sodium azide and 0.2% BSA).

STORAGE

Store at 4° C, **DO NOT FREEZE**. Stable for one year from the date of shipment. Non-hazardous. No MSDS required.

APPLICATIONS

ACLP (N-14) is recommended for detection of ACLP of human origin by Western Blotting (starting dilution 1:200, dilution range 1:100-1:1000), immunofluorescence (starting dilution 1:50, dilution range 1:50-1:500) and solid phase ELISA (starting dilution 1:30, dilution range 1:30-1:3000).

Suitable for use as control antibody for ACLP/AEBP1 siRNA (h): sc-40327, ACLP/AEBP1 shRNA Plasmid (h): sc-40327-SH and ACLP/AEBP1 shRNA (h) Lentiviral Particles: sc-40327-V.

Molecular Weight of ACLP: 79 kDa.

RECOMMENDED SECONDARY REAGENTS

To ensure optimal results, the following support (secondary) reagents are recommended: 1) Western Blotting: use donkey anti-goat IgG-HRP: sc-2020 (dilution range: 1:2000-1:100,000) or Cruz Marker™ compatible donkey anti-goat IgG-HRP: sc-2033 (dilution range: 1:2000-1:5000), Cruz Marker™ Molecular Weight Standards: sc-2035, TBS Blotto A Blocking Reagent: sc-2333 and Western Blotting Luminol Reagent: sc-2048. 2) Immunoprecipitation: use Protein A/G PLUS-Agarose: sc-2003 (0.5 ml agarose/2.0 ml). 3) Immunofluorescence: use donkey anti-goat IgG-FITC: sc-2024 (dilution range: 1:100-1:400) or donkey anti-goat IgG-TR: sc-2783 (dilution range: 1:100-1:400) with UltraCruz™ Mounting Medium: sc-24941.

RESEARCH USE

For research use only, not for use in diagnostic procedures.

PROTOCOLS

See our web site at www.scbt.com or our catalog for detailed protocols and support products.



MONOS
Satisfation
Guaranteed

Try **ACLP/AEBP1 (G-1): sc-271374**, our highly recommended monoclonal alternative to ACLP (N-14).