

donkey anti-rabbit IgG-HRP: sc-2313

BACKGROUND

Santa Cruz Biotechnology's secondary antibodies are available conjugated to either an enzyme, biotin or fluorophore for use in a variety of antibody-based applications including Western Blot, immunostaining, flow cytometry and ELISA. Secondary antibodies are commonly affinity purified against immobilized whole IgG isotypes, including IgG₁, IgG_{2a}, IgG_{2b}, IgG₃ and IgG₄. Santa Cruz Biotechnology offers a wide selection of secondary antibodies, which are used in conjunction with our Cruz Marker™ molecular weight standards. We also provide specialized secondaries, such as pre-adsorbed secondary antibodies, which are pre-adsorbed with human IgG and mouse IgG for immunoglobulin-rich tissues and cells, F(ab')₂ fragment secondary antibodies that reduce non-specific secondary antibody binding to Fc receptors on the cell surface, and isotype-specific secondary antibodies against IgM, IgA and IgY.

SOURCE

donkey anti-rabbit IgG-HRP is an affinity purified secondary antibody raised in donkey against rabbit IgG and conjugated to HRP (horseradish peroxidase).

PRODUCT

Each vial contains 200 µg IgG in 0.5 ml of 1X PBS containing 40% glycerol.

APPLICATIONS

donkey anti-rabbit IgG-HRP is recommended for detection of rabbit IgG by Western Blotting (starting dilution: 1:5000, dilution range 1:5000-1:10000; optimal dilution to be determined by titration).

RECOMMENDED SUPPORT PRODUCTS

- Western Blotting Luminol Reagent, for 2,000 cm² membrane area: sc-2048
- RIPA Lysis Buffer, 50 ml, cell lysis buffer with protease inhibitors: sc-24948
- Electrophoresis Sample Buffer, 2X, 25 ml, reducing buffer: sc-24945
- Complete™ Protease Inhibitor Cocktail Tablet, 20 tablets: sc-29130
- Running Buffer, 10X, 1 L, TRIS-Glycine WB running buffer, pH 8.3: sc-24949
- Towbin, with SDS, 10X, 1 L, WB transfer buffer pH 8.3: sc-24954
- Bovine Serum Albumin (BSA), 100 g, blocking/incubation agent: sc-2323
- TBS Blotto A, lyophilized powder in single-use bottle: sc-2333
- UltraCruz™ PVDF Transfer Membrane, 0.45 µm, 30 cm x 3 m roll: sc-3723
- UltraCruz™ Nitrocellulose Pure Transfer Membrane, 0.22 µm, 30 cm x 3 m roll: sc-3718
- UltraCruz™ Tissue Culture Dish, 100 mm polystyrene dish: sc-200286
- UltraCruz™ Cell Scrapers, 25 cm, sterile, 100 per case: sc-213229
- UltraCruz™ Electrophoresis Cell: sc-201625: runs up to 10 or 15 sample by SDS – PAGE protein electrophoresis
- UltraCruz™ Autoradiography Film, Blue, 8 x 1, 100 sheets: sc-201697
- UltraCruz™ Gel Incubation Trays, 100 per pack: sc-201755 (blue), sc-201756 (green), sc-201757 (pink), sc-201758 (yellow), sc-201759 (orange)


RESEARCH USE

For research use only, not for use in diagnostic procedures.

STORAGE

Store at 4° C, ****DO NOT FREEZE****. Stable for one year from the date of shipment. Non-hazardous. No MSDS required.

DATA



donkey anti-rabbit IgG-HRP: sc-2313. Western blot analysis of p53 expression in Lac-Z whole cell lysate. Antibody tested: p53 (FL-393): sc-6243.

SELECT PRODUCT CITATIONS

- Wu, J. 2004. β -Amyloid directly inhibits human α 4 β 2-nicotinic acetylcholine receptors heterologously expressed in human SH-EP1 cells. *J. Biol. Chem.* 279: 37842-37851.
- Ducut Sigala, J.L., et al. 2004. Activation of transcription factor NF κ B requires ELKS, an I κ B kinase regulatory subunit. *Science* 304: 1963-1967.
- Wang, C.T., et al. 2009. The effects of green tea (-)-epigallocatechin-3-gallate on reactive oxygen species in 3T3-L1 preadipocytes and adipocytes depend on the glutathione and 67 kDa laminin receptor pathways. *Mol. Nutr. Food Res.* 53: 349-360.
- Li, Y., et al. 2009. Modularity of the oncoprotein-like properties of platelet glycoprotein Ib α . *J. Biol. Chem.* 284: 1410-1418.
- Borbély, G., et al. 2010. Small-molecule inhibitors of NADPH oxidase 4. *J. Med. Chem.* 53: 6758-6762.
- Alexaki, V.I., et al. 2010. GLI2-mediated melanoma invasion and metastasis. *J. Natl. Cancer Inst.* 102: 1148-1159.
- Muehlbauer, S.M., et al. 2010. Proteasome inhibitors prevent caspase-1-mediated disease in rodents challenged with anthrax lethal toxin. *Am. J. Pathol.* 177: 735-743.
- Wang, H., et al. 2010. Oridonin induces G₂/M cell cycle arrest and apoptosis through MAPK and p53 signaling pathways in HepG2 cells. *Oncol. Rep.* 24: 647-651.
- Wang, H., et al. 2010. Proteomic and functional analyses reveal the potential involvement of endoplasmic reticulum stress and α -CP1 in the anti-cancer activities of oridonin in HepG2 cells. *Integr. Cancer Ther.* 10: 160-167.
- Lau, T., et al. 2010. Somatodendritic serotonin release and re-uptake in mouse embryonic stem cell-derived serotonergic neurons. *Neurochem. Int.* 57: 969-978.
- Basco, D., et al. 2011. Absence of aquaporin-4 in skeletal muscle alters proteins involved in bioenergetic pathways and calcium handling. *PLoS ONE* 6: e19225.