

# JIP-3 (V-18): sc-23401

## BACKGROUND

The JNK-interacting proteins (JIPs) are a family of scaffold proteins that mediate JNK signaling by organizing specific components of the MAPK cascade together to form a functional JNK signaling molecule. JIP-3 (JNK-interacting protein 3), also known as JSAP1 or MAPK8IP3 (Mitogen-activated protein kinase 8-interacting protein 3), is a 1,336 amino acid protein that localizes to the cytoplasm and belongs to the JIP family. Expressed in a variety of tissues, including brain and heart, JIP-3 forms homo- or heterooligomeric complexes that can interact with several components of the JNK signaling pathway, thereby functioning as a regulator of kinesin-dependent axonal transport that may also play a role in scaffold formation within neuronal cells. Human JIP-3, which may be phosphorylated upon DNA damage, shares 69% similarity with its mouse counterpart, suggesting a conserved role between species. Multiple isoforms of JIP-3 exist due to alternative splicing events.

## REFERENCES

1. Ito, M., Yoshioka, K., Akechi, M., Yamashita, S., Takamatsu, N., Sugiyama, K., Hibi, M., Nakabeppu, Y., Shiba, T. and Yamamoto, K.I. 1999. JSAP1, a novel Jun N-terminal protein kinase (JNK)-binding protein that functions as a scaffold factor in the JNK signaling pathway. *Mol. Cell. Biol.* 19: 7539-7548.
2. Kelkar, N., Gupta, S., Dickens, M. and Davis, R.J. 2000. Interaction of a mitogen-activated protein kinase signaling module with the neuronal protein JIP3. *Mol. Cell. Biol.* 20: 1030-1043.
3. Ferrer, I., Friguls, B., Dalfo, E. and Planas, A.M. 2003. Early modifications in the expression of mitogen-activated protein kinase (MAPK/ERK), stress-activated kinases SAPK/JNK and p38, and their phosphorylated substrates following focal cerebral ischemia. *Acta Neuropathol.* 105: 425-37.
4. Murr, M.M., Yang, J., Fier, A., Gallagher, S.F., Carter, G., Gower, W.R. and Norman, J.G. 2003. Regulation of Kupffer cell TNF gene expression during experimental acute pancreatitis. The role of p38-MAPK, ERK 1/2, SAPK/JNK, and NFκB. *J. Gastrointest. Surg.* 7: 20-25.
5. Smith, W.E., Kane, A.V., Campbell, S.T., Acheson, D.W., Cochran, B.H. and Thorpe, C.M. 2003. Shiga toxin 1 triggers a ribotoxic stress response leading to p38 and JNK activation and induction of apoptosis in intestinal epithelial cells. *Infect. Immun.* 71: 1497-1504.

## CHROMOSOMAL LOCATION

Genetic locus: MAPK8IP3 (human) mapping to 16p13.3; Mapk8ip3 (mouse) mapping to 17 A3.3.

## SOURCE

JIP-3 (V-18) is an affinity purified goat polyclonal antibody raised against a peptide mapping near the C-terminus of JIP-3 of human origin.

## STORAGE

Store at 4° C, **\*\*DO NOT FREEZE\*\***. Stable for one year from the date of shipment. Non-hazardous. No MSDS required.

## PRODUCT

Each vial contains 200 µg IgG in 1.0 ml of PBS with < 0.1% sodium azide and 0.1% gelatin.

Blocking peptide available for competition studies, sc-23401 P, (100 µg peptide in 0.5 ml PBS containing < 0.1% sodium azide and 0.2% BSA).

## APPLICATIONS

JIP-3 (V-18) is recommended for detection of JIP-3 of mouse, rat and human origin by Western Blotting (starting dilution 1:200, dilution range 1:100-1:1000), immunofluorescence (starting dilution 1:50, dilution range 1:50-1:500) and solid phase ELISA (starting dilution 1:30, dilution range 1:30-1:3000); may cross-react with JIP-4.

JIP-3 (V-18) is also recommended for detection of JIP-3 in additional species, including equine, canine, bovine, porcine and avian.

Suitable for use as control antibody for JIP-3 siRNA (h): sc-37123, JIP-3 siRNA (m): sc-37124, JIP-3 shRNA Plasmid (h): sc-37123-SH, JIP-3 shRNA Plasmid (m): sc-37124-SH, JIP-3 shRNA (h) Lentiviral Particles: sc-37123-V and JIP-3 shRNA (m) Lentiviral Particles: sc-37124-V.

Molecular Weight of JIP-3: 147 kDa.

Positive Controls: IMR-32 cell lysate: sc-2409, rat brain extract: sc-2392 or PC-12 cell lysate: sc-2250.

## RECOMMENDED SECONDARY REAGENTS

To ensure optimal results, the following support (secondary) reagents are recommended: 1) Western Blotting: use donkey anti-goat IgG-HRP: sc-2020 (dilution range: 1:2000-1:100,000) or Cruz Marker™ compatible donkey anti-goat IgG-HRP: sc-2033 (dilution range: 1:2000-1:5000), Cruz Marker™ Molecular Weight Standards: sc-2035, TBS Blotto A Blocking Reagent: sc-2333 and Western Blotting Luminol Reagent: sc-2048. 2) Immunofluorescence: use donkey anti-goat IgG-FITC: sc-2024 (dilution range: 1:100-1:400) or donkey anti-goat IgG-TR: sc-2783 (dilution range: 1:100-1:400) with UltraCruz™ Mounting Medium: sc-24941.

## RESEARCH USE

For research use only, not for use in diagnostic procedures.

## PROTOCOLS

See our web site at [www.scbt.com](http://www.scbt.com) or our catalog for detailed protocols and support products.


 MONOS  
 Satisfation  
 Guaranteed

Try **JIP-3 (F-6): sc-46663** or **JIP-3 (D-3): sc-48392**, our highly recommended monoclonal alternatives to JIP-3 (V-18).