

Neurabin-I (N-20): sc-23494

BACKGROUND

Brain-specific Neurabin-I (neural tissue-specific F-Actin binding protein I) is highly concentrated in the synapse of developed neurons; it localizes in the growth cone lamellipodia during neuronal development. Suppression of endogenous Neurabin in rat hippocampal neurons inhibits neurite formation. Neurabin-I recruits active PP1 via a PP1-docking sequence; mutation of the PP1-binding motif halts filopodia and restores stress fibers in Neurabin-I expressing cells. Neurabin-II (Spinophilin) is ubiquitously expressed but is abundantly expressed in brain. Neurabin-II localizes to neuronal dendritic spines, which are the specialized protrusions from dendritic shafts that receive most of the excitatory input in the CNS. Neurabin-II may regulate dendritic spine properties as Neurabin-III(-) mice have increased spine density during development *in vitro* and exhibit altered filopodial formation in cultured cells. Neurabin may also play a role in glutamatergic transmission as Neurabin-III(-) mice exhibit reduced long-term depression and resistance to kainate-induced seizures and neuronal apoptosis. Neurabin-II complexes with the catalytic subunit of protein phosphatase 1 (PP1) *in vitro*, thus modulating the activity of PP1

REFERENCES

1. Nakanishi, H., Obaishi, H., Satoh, A., Wada, M., Mandai, K., Satoh, K., Nishioka, H., Matsuura, Y., Mizoguchi, A. and Takai, Y. 1997. Neurabin: a novel neural tissue-specific Actin filament-binding protein involved in neurite formation. *J. Cell Biol.* 139: 951-961.
2. Allen, P.B., Ouimet, C.C. and Greengard, P. 1997. Spinophilin, a novel protein phosphatase 1 binding protein localized to dendritic spines. *Proc. Natl. Acad. Sci. USA* 94: 9956-9961.
3. McAvoy, T., Allen, P.B., Obaishi, H., Nakanishi, H., Takai, Y., Greengard, P., Nairn, A.C. and Hemmings, H.C., Jr. 1999. Regulation of Neurabin-I interaction with protein phosphatase 1 by phosphorylation. *Biochemistry* 38: 12943-12949.
4. Feng, J., Yan, Z., Ferreira, A., Tomizawa, K., Liauw, J.A., Zhuo, M., Allen, P.B., Ouimet, C.C. and Greengard, P. 2000. Spinophilin regulates the formation and function of dendritic spines. *Proc. Natl. Acad. Sci. USA* 97: 9287-9292
5. Oliver, C.J., Terry-Lorenzo, R.T., Elliott, E., Bloomer, W.A., Li, S., Brautigan, D.L., Colbran, R.J. and Shenolikar, S. 2002. Targeting protein phosphatase 1 (PP1) to the actin cytoskeleton: the Neurabin I/PP1 complex regulates cell morphology. *Mol. Cell. Biol.* 22: 4690-701.

CHROMOSOMAL LOCATION

Genetic locus: PPP1R9A (human) mapping to 7q21.3; Ppp1r9a (mouse) mapping to 6 A1.

SOURCE

Neurabin-I (N-20) is an affinity purified goat polyclonal antibody raised against a peptide mapping within an internal region of Neurabin-I of human origin.

RESEARCH USE

For research use only, not for use in diagnostic procedures.

PRODUCT

Each vial contains 200 µg IgG in 1.0 ml of PBS with < 0.1% sodium azide and 0.1% gelatin.

Blocking peptide available for competition studies, sc-23494 P, (100 µg peptide in 0.5 ml PBS containing < 0.1% sodium azide and 0.2% BSA).

APPLICATIONS

Neurabin-I (N-20) is recommended for detection of Neurabin-I of mouse, rat and human origin by Western Blotting (starting dilution 1:200, dilution range 1:100-1:1000), immunofluorescence (starting dilution 1:50, dilution range 1:50-1:500) and solid phase ELISA (starting dilution 1:30, dilution range 1:30-1:3000).

Suitable for use as control antibody for Neurabin-I siRNA (h): sc-45982, Neurabin-I siRNA (m): sc-45983, Neurabin-I shRNA Plasmid (h): sc-45982-SH, Neurabin-I shRNA Plasmid (m): sc-45983-SH, Neurabin-I shRNA (h) Lentiviral Particles: sc-45982-V and Neurabin-I shRNA (m) Lentiviral Particles: sc-45983-V.

Molecular Weight of Neurabin-I: 180 kDa.

Positive Controls: mouse brain extract: sc-2253 or PC-12 cell lysate: sc-2250.

RECOMMENDED SECONDARY REAGENTS

To ensure optimal results, the following support (secondary) reagents are recommended: 1) Western Blotting: use donkey anti-goat IgG-HRP: sc-2020 (dilution range: 1:2000-1:100,000) or Cruz Marker™ compatible donkey anti-goat IgG-HRP: sc-2033 (dilution range: 1:2000-1:5000), Cruz Marker™ Molecular Weight Standards: sc-2035, TBS Blotto A Blocking Reagent: sc-2333 and Western Blotting Luminol Reagent: sc-2048. 2) Immunofluorescence: use donkey anti-goat IgG-FITC: sc-2024 (dilution range: 1:100-1:400) or donkey anti-goat IgG-TR: sc-2783 (dilution range: 1:100-1:400) with UltraCruz™ Mounting Medium: sc-24941.

STORAGE

Store at 4° C, ****DO NOT FREEZE****. Stable for one year from the date of shipment. Non-hazardous. No MSDS required.

PROTOCOLS

See our web site at www.scbt.com or our catalog for detailed protocols and support products.

MONOS
Satisfaction
Guaranteed

Try **Neurabin-I (D-4): sc-377407** or **Neurabin-I (52): sc-136327**, our highly recommended monoclonal alternatives to Neurabin-I (N-20).