

# HCS (C-15): sc-23734

## BACKGROUND

Holocarboxylase synthetase (HCS, HLCS) catalyzes the posttranslational covalent attachment of biotin to biotin-dependent carboxylases. These carboxylases include pyruvate carboxylase, propionyl coenzyme A (CoA) carboxylase,  $\beta$ -methylcrotonoyl CoA carboxylase, and acetyl CoA carboxylase. HCS is an obligate participant in biotin-mediated regulation of its own expression and of biotin-dependent carboxylase mRNA levels in human cells. In addition, biotin regulates the genetic expression of HCS and mitochondrial carboxylases. Multiple carboxylase deficiency (MCD) is a rare metabolic life-threatening disease accountable for early onset biotin-responsive multiple carboxylase deficiency. MCD presents with lactic acidosis, tachypnea, temperature instability, and shock in neonates. Interestingly, MCD can be detected prenatally and, if diagnosed, prenatally administered biotin is effectively taken up by the fetus and prevents functional deficiency of the carboxylases in an affected newborn. The HCS gene maps to human chromosome 21q22.13 and encodes a cytoplasmic and mitochondrial protein.

## REFERENCES

1. Squires, L., et al. 1997. Resolution of subependymal cysts in neonatal holocarboxylase synthetase deficiency. *Dev. Med. Child Neurol.* 39: 267-269.
2. Baumgartner, E.R., et al. 1997. Multiple carboxylase deficiency: inherited and acquired disorders of biotin metabolism. *Int. J. Vitam. Nutr. Res.* 67: 377-384.
3. Suormala, T., et al. 1998. Late-onset holocarboxylase synthetase-deficiency: pre- and post-natal diagnosis and evaluation of effectiveness of antenatal biotin therapy. *Eur. J. Pediatr.* 157: 570-575.
4. Hwu, W.L., et al. 2000. Late-onset holocarboxylase synthetase deficiency with homologous R508W mutation. *J. Formos. Med. Assoc.* 99: 174-177.
5. Yang, X., et al. 2001. Structure of human holocarboxylase synthetase gene and mutation spectrum of holocarboxylase synthetase deficiency. *Hum. Genet.* 109: 526-534.
6. Rodriguez-Melendez, R., et al. 2001. Biotin regulates the genetic expression of holocarboxylase synthetase and mitochondrial carboxylases in rats. *J. Nutr.* 131: 1909-1913.
7. Solorzano-Vargas, R.S., et al. 2002. Holocarboxylase synthetase is an obligate participant in biotin-mediated regulation of its own expression and of biotin-dependent carboxylases mRNA levels in human cells. *Proc. Natl. Acad. Sci. USA* 99: 5325-5330.
8. Online Mendelian Inheritance in Man, OMIM<sup>™</sup>. 2002. Johns Hopkins University, Baltimore, MD. MIM Number: 253270. World Wide Web URL: <http://www.ncbi.nlm.nih.gov/omim/>

## CHROMOSOMAL LOCATION

Genetic locus: HLCS (human) mapping to 21q22.13; Hlcs (mouse) mapping to 16 C4.

## RESEARCH USE

For research use only, not for use in diagnostic procedures.

## SOURCE

HCS (C-15) is an affinity purified goat polyclonal antibody raised against a peptide mapping at the C-terminus of HCS of human origin.

## PRODUCT

Each vial contains 200  $\mu$ g IgG in 1.0 ml of PBS with < 0.1% sodium azide and 0.1% gelatin.

Blocking peptide available for competition studies, sc-23734 P, (100  $\mu$ g peptide in 0.5 ml PBS containing < 0.1% sodium azide and 0.2% BSA).

## APPLICATIONS

HCS (C-15) is recommended for detection of HCS of mouse, rat and human origin by Western Blotting (starting dilution 1:200, dilution range 1:100-1:1000), immunofluorescence (starting dilution 1:50, dilution range 1:50-1:500) and solid phase ELISA (starting dilution 1:30, dilution range 1:30-1:3000).

HCS (C-15) is also recommended for detection of HCS in additional species, including equine, canine, bovine, porcine and avian.

Suitable for use as control antibody for HCS siRNA (h): sc-72198, HCS siRNA (m): sc-72199, HCS shRNA Plasmid (h): sc-72198-SH, HCS shRNA Plasmid (m): sc-72199-SH, HCS shRNA (h) Lentiviral Particles: sc-72198-V and HCS shRNA (m) Lentiviral Particles: sc-72199-V.

Molecular Weight of HCS: 86 kDa.

Positive Controls: JAR cell lysate: sc-2276.

## RECOMMENDED SECONDARY REAGENTS

To ensure optimal results, the following support (secondary) reagents are recommended: 1) Western Blotting: use donkey anti-goat IgG-HRP: sc-2020 (dilution range: 1:2000-1:100,000) or Cruz Marker<sup>™</sup> compatible donkey anti-goat IgG-HRP: sc-2033 (dilution range: 1:2000-1:5000), Cruz Marker<sup>™</sup> Molecular Weight Standards: sc-2035, TBS Blotto A Blocking Reagent: sc-2333 and Western Blotting Luminol Reagent: sc-2048. 2) Immunofluorescence: use donkey anti-goat IgG-FITC: sc-2024 (dilution range: 1:100-1:400) or donkey anti-goat IgG-TR: sc-2783 (dilution range: 1:100-1:400) with UltraCruz<sup>™</sup> Mounting Medium: sc-24941.

## STORAGE

Store at 4° C, **\*\*DO NOT FREEZE\*\***. Stable for one year from the date of shipment. Non-hazardous. No MSDS required.

## PROTOCOLS

See our web site at [www.scbt.com](http://www.scbt.com) or our catalog for detailed protocols and support products.