

DAND5 (M-14): sc-240289

BACKGROUND

The DAN (differential screening-selected gene aberrative in neuroblastoma) protein family contains antagonists of bone morphogenetic protein (BMP) signaling that are expressed in the neural crest. All family members are secreted proteins that act as BMP antagonists in embryonic explants and are expressed in the proximal airway epithelium of the lung during embryonic development. DAND5 (DAN domain family, member 5), also known as Cerl-2 (cerberus-like protein 2), CKTSF1B3 (cysteine knot superfamily 1, BMP antagonist 3) or Gremlin-3, is a 189 amino acid member of the DAN protein family. DAND5 is a secreted protein that contains one CTCK (C-terminal cysteine knot-like) domain. DAND5 is thought to play a role in the correct specification of the left-right axis during development. DAND5 may antagonize Nodal and BMP-4 signaling.

REFERENCES

- Ohtori, S., Yamamoto, T., Ino, H., Hanaoka, E., Shinbo, J., Ozaki, T., Takada, N., Nakamura, Y., Chiba, T., Nakagawara, A., Sakiyama, S., Sakashita, Y., Takahashi, K., Tanaka, K., Yamagata, M., Yamazaki, M., Shimizu, S. and Moriya, H. 2002. Differential screening-selected gene aberrative in neuroblastoma protein modulates inflammatory pain in the spinal dorsal horn. *Neuroscience* 110: 579-586.
- Kim, A.S. and Pleasure, S.J. 2003. Expression of the BMP antagonist Dan during murine forebrain development. *Brain Res. Dev. Brain Res.* 145: 159-162.
- Gerlach-Bank, L.M., Cleveland, A.R. and Barald, K.F. 2004. DAN directs endolymphatic sac and duct outgrowth in the avian inner ear. *Dev. Dyn.* 229: 219-230.
- Marques, S., Borges, A.C., Silva, A.C., Freitas, S., Cordenonsi, M. and Belo, J.A. 2004. The activity of the Nodal antagonist Cerl-2 in the mouse node is required for correct L/R body axis. *Genes Dev.* 18: 2342-2347.
- Avsian-Kretchmer, O. and Hsueh, A.J. 2004. Comparative genomic analysis of the eight-membered ring cysteine knot-containing bone morphogenetic protein antagonists. *Mol. Endocrinol.* 18: 1-12.
- Katoh, M. and Katoh, M. 2004. Identification and characterization of human CKTSF1B2 and CKTSF1B3 genes *in silico*. *Oncol. Rep.* 12: 423-427.

CHROMOSOMAL LOCATION

Genetic locus: Dand5 (mouse) mapping to 8 C3.

SOURCE

DAND5 (M-14) is an affinity purified goat polyclonal antibody raised against a peptide mapping near the C-terminus of DAND5 of mouse origin.

PRODUCT

Each vial contains 200 µg IgG in 1.0 ml of PBS with < 0.1% sodium azide and 0.1% gelatin.

Blocking peptide available for competition studies, sc-240289 P, (100 µg peptide in 0.5 ml PBS containing < 0.1% sodium azide and 0.2% BSA).

APPLICATIONS

DAND5 (M-14) is recommended for detection of DAND5 of mouse and rat origin by Western Blotting (starting dilution 1:200, dilution range 1:100-1:1000), immunofluorescence (starting dilution 1:50, dilution range 1:50-1:500) and solid phase ELISA (starting dilution 1:30, dilution range 1:30-1:3000).

Suitable for use as control antibody for DAND5 siRNA (m): sc-142871, DAND5 shRNA Plasmid (m): sc-142871-SH and DAND5 shRNA (m) Lentiviral Particles: sc-142871-V.

Molecular Weight of DAND5: 20 kDa.

RECOMMENDED SECONDARY REAGENTS

To ensure optimal results, the following support (secondary) reagents are recommended: 1) Western Blotting: use donkey anti-goat IgG-HRP: sc-2020 (dilution range: 1:2000-1:100,000) or Cruz Marker™ compatible donkey anti-goat IgG-HRP: sc-2033 (dilution range: 1:2000-1:5000), Cruz Marker™ Molecular Weight Standards: sc-2035, TBS Blotto A Blocking Reagent: sc-2333 and Western Blotting Luminol Reagent: sc-2048. 2) Immunofluorescence: use donkey anti-goat IgG-FITC: sc-2024 (dilution range: 1:100-1:400) or donkey anti-goat IgG-TR: sc-2783 (dilution range: 1:100-1:400) with UltraCruz™ Mounting Medium: sc-24941.

STORAGE

Store at 4° C, ****DO NOT FREEZE****. Stable for one year from the date of shipment. Non-hazardous. No MSDS required.

RESEARCH USE

For research use only, not for use in diagnostic procedures.

PROTOCOLS

See our web site at www.scbt.com or our catalog for detailed protocols and support products.