

SLC7A10 (P-13): sc-240887

BACKGROUND

SLC7A10 (solute carrier family 7, (neutral amino acid transporter, γ^+ system member 10), also known as ASC1 (asc-type amino acid transporter 1), is a 523 amino acid multi-pass membrane protein that belongs to the amino acid-polyamine-organocation (APC) superfamily. SLC7A10 forms a disulfide-linked heterodimer with CD98 to regulate high-affinity transport of several neutral amino acids, including D-serine. Expressed in pancreas, skeletal muscle, placenta, heart, brain, kidney, lung and liver, SLC7A10 is encoded by a gene that maps to human chromosome 19q13.11 and is considered a candidate gene for cystinuria. Chromosome 19 consists of over 63 million bases, houses approximately 1,400 genes and is recognized for having the greatest gene density of the human chromosomes.

REFERENCES

1. Fukasawa, Y., et al. 2000. Identification and characterization of a Na⁺-independent neutral amino acid transporter that associates with the 4F2 heavy chain and exhibits substrate selectivity for small neutral D- and L-amino acids. *J. Biol. Chem.* 275: 9690-9698.
2. Nakauchi, J., et al. 2000. Cloning and characterization of a human brain Na⁺-independent transporter for small neutral amino acids that transports D-serine with high affinity. *Neurosci. Lett.* 287: 231-235.
3. Leclerc, D., et al. 2001. Is the SLC7A10 gene on chromosome 19 a candidate locus for cystinuria? *Mol. Genet. Metab.* 73: 333-339.
4. Dehal, P., et al. 2001. Human chromosome 19 and related regions in mouse: conservative and lineage-specific evolution. *Science* 293: 104-111.
5. Pineda, M., et al. 2004. The amino acid transporter asc-1 is not involved in cystinuria. *Kidney Int.* 66: 1453-1464.
6. Xu, D., et al. 2005. Metabolic activation-related CD147-CD98 complex. *Mol. Cell. Proteomics* 4: 1061-1071.
7. Online Mendelian Inheritance in Man, OMIM[™]. 2007. Johns Hopkins University, Baltimore, MD. MIM Number: 607959. World Wide Web URL: <http://www.ncbi.nlm.nih.gov/omim/>
8. Bröer, S. 2008. Amino acid transport across mammalian intestinal and renal epithelia. *Physiol. Rev.* 88: 249-286.

CHROMOSOMAL LOCATION

Genetic locus: SLC7A10 (human) mapping to 19q13.11; Slc7a10 (mouse) mapping to 7 B1.

SOURCE

SLC7A10 (P-13) is an affinity purified goat polyclonal antibody raised against a peptide mapping within an internal region of SLC7A10 of human origin.

PRODUCT

Each vial contains 200 μ g IgG in 1.0 ml of PBS with < 0.1% sodium azide and 0.1% gelatin.

Blocking peptide available for competition studies, sc-240887 P, (100 μ g peptide in 0.5 ml PBS containing < 0.1% sodium azide and 0.2% BSA).

APPLICATIONS

SLC7A10 (P-13) is recommended for detection of SLC7A10 of mouse, rat and human origin by Western Blotting (starting dilution 1:200, dilution range 1:100-1:1000), immunofluorescence (starting dilution 1:50, dilution range 1:50-1:500) and solid phase ELISA (starting dilution 1:30, dilution range 1:30-1:3000); non cross-reactive with other SLC7A family members.

Suitable for use as control antibody for SLC7A10 siRNA (h): sc-97607, SLC7A10 siRNA (m): sc-153578, SLC7A10 shRNA Plasmid (h): sc-97607-SH, SLC7A10 shRNA Plasmid (m): sc-153578-SH, SLC7A10 shRNA (h) Lentiviral Particles: sc-97607-V and SLC7A10 shRNA (m) Lentiviral Particles: sc-153578-V.

Molecular Weight of SLC7A10: 57 kDa.

RECOMMENDED SECONDARY REAGENTS

To ensure optimal results, the following support (secondary) reagents are recommended: 1) Western Blotting: use donkey anti-goat IgG-HRP: sc-2020 (dilution range: 1:2000-1:100,000) or Cruz Marker[™] compatible donkey anti-goat IgG-HRP: sc-2033 (dilution range: 1:2000-1:5000), Cruz Marker[™] Molecular Weight Standards: sc-2035, TBS Blotto A Blocking Reagent: sc-2333 and Western Blotting Luminol Reagent: sc-2048. 2) Immunofluorescence: use donkey anti-goat IgG-FITC: sc-2024 (dilution range: 1:100-1:400) or donkey anti-goat IgG-TR: sc-2783 (dilution range: 1:100-1:400) with UltraCruz[™] Mounting Medium: sc-24941.

STORAGE

Store at 4° C, ****DO NOT FREEZE****. Stable for one year from the date of shipment. Non-hazardous. No MSDS required.

RESEARCH USE

For research use only, not for use in diagnostic procedures.

PROTOCOLS

See our web site at www.scbt.com or our catalog for detailed protocols and support products.