

ACSF2 (E-20): sc-241776

BACKGROUND

ACSF2 (acyl-CoA synthetase family member 2), also known as PPAR γ binding, long chain fatty acid acyl Co-A ligase like, is a 615 amino acid protein belonging to the ATP-dependent AMP-binding enzyme family. Encoded by a gene that maps to human chromosome 17q21.33, ACSF2 contains a PPAR γ -recognition element in its promoter sequence. Exhibiting mitochondrial subcellular localization, ACSF2 forms a thioester with CoA, thereby catalyzing the initial reaction in fatty acid metabolism. ACSF2 participates in ATP binding, ligase activity and adipocyte differentiation, and displays a preference for medium-chain substrates. ACSF2 may be a potential prognostic molecular marker for breast cancer.

REFERENCES

1. Grützmann, R., et al. 2003. Systematic isolation of genes differentially expressed in normal and cancerous tissue of the pancreas. *Pancreatology* 3: 169-178.
2. Yang, X.J., et al. 2004. Gene expression profiling of renal medullary carcinoma: potential clinical relevance. *Cancer* 100: 976-985.
3. Kim, J.M., et al. 2005. Identification of gastric cancer-related genes using a cDNA microarray containing novel expressed sequence tags expressed in gastric cancer cells. *Clin. Cancer Res.* 11: 473-482.
4. Sørlie, T., et al. 2006. Distinct molecular mechanisms underlying clinically relevant subtypes of breast cancer: gene expression analyses across three different platforms. *BMC Genomics* 7: 127.
5. Yao, J., et al. 2006. Combined cDNA array comparative genomic hybridization and serial analysis of gene expression analysis of breast tumor progression. *Cancer Res.* 66: 4065-4078.
6. Gianazza, E., et al. 2006. Coordinated and reversible reduction of enzymes involved in terminal oxidative metabolism in skeletal muscle mitochondria from a riboflavin-responsive, multiple acyl-CoA dehydrogenase deficiency patient. *Electrophoresis* 27: 1182-1198.
7. Perera, R.J., et al. 2006. Identification of novel PPAR γ target genes in primary human adipocytes. *Gene* 369: 90-99.
8. Watkins, P.A., et al. 2007. Evidence for 26 distinct acyl-coenzyme A synthetase genes in the human genome. *J. Lipid Res.* 48: 2736-2750.
9. Pei, Z., et al. 2009. Acyl-CoA synthetase VL3 knockdown inhibits human glioma cell proliferation and tumorigenicity. *Cancer Res.* 69: 9175-9182.

CHROMOSOMAL LOCATION

Genetic locus: ACSF2 (human) mapping to 17q21.33; Acsf2 (mouse) mapping to 11 D.

SOURCE

ACSF2 (E-20) is an affinity purified goat polyclonal antibody raised against a peptide mapping within an internal region of ACSF2 of human origin.

RESEARCH USE

For research use only, not for use in diagnostic procedures.

PRODUCT

Each vial contains 200 μ g IgG in 1.0 ml of PBS with < 0.1% sodium azide and 0.1% gelatin.

Blocking peptide available for competition studies, sc-241778 P, (100 μ g peptide in 0.5 ml PBS containing < 0.1% sodium azide and 0.2% BSA).

APPLICATIONS

ACSF2 (E-20) is recommended for detection of ACSF2 of mouse, rat and human origin by Western Blotting (starting dilution 1:200, dilution range 1:100-1:1000), immunofluorescence (starting dilution 1:50, dilution range 1:50-1:500) and solid phase ELISA (starting dilution 1:30, dilution range 1:30-1:3000); non cross-reactive with ACSF1, ACSF3 or ACSF4.

ACSF2 (E-20) is also recommended for detection of ACSF2 in additional species, including equine, canine, bovine and porcine.

Suitable for use as control antibody for ACSF2 siRNA (h): sc-93920, ACSF2 siRNA (m): sc-140828, ACSF2 shRNA Plasmid (h): sc-93920-SH, ACSF2 shRNA Plasmid (m): sc-140828-SH, ACSF2 shRNA (h) Lentiviral Particles: sc-93920-V and ACSF2 shRNA (m) Lentiviral Particles: sc-140828-V.

Molecular Weight of ACSF2: 68 kDa.

RECOMMENDED SECONDARY REAGENTS

To ensure optimal results, the following support (secondary) reagents are recommended: 1) Western Blotting: use donkey anti-goat IgG-HRP: sc-2020 (dilution range: 1:2000-1:100,000) or Cruz Marker[™] compatible donkey anti-goat IgG-HRP: sc-2033 (dilution range: 1:2000-1:5000), Cruz Marker[™] Molecular Weight Standards: sc-2035, TBS Blotto A Blocking Reagent: sc-2333 and Western Blotting Luminol Reagent: sc-2048. 2) Immunofluorescence: use donkey anti-goat IgG-FITC: sc-2024 (dilution range: 1:100-1:400) or donkey anti-goat IgG-TR: sc-2783 (dilution range: 1:100-1:400) with UltraCruz[™] Mounting Medium: sc-24941.

STORAGE

Store at 4° C, ****DO NOT FREEZE****. Stable for one year from the date of shipment. Non-hazardous. No MSDS required.

PROTOCOLS

See our web site at www.scbt.com or our catalog for detailed protocols and support products.