

AIM1L (N-16): sc-241795

BACKGROUND

AIM1L (absent in melanoma 1-like), also known as CRYBG2 (β/γ crystallin domain-containing protein 2), is a 616 amino acid protein that participates in sugar binding. AIM1L contains ten β/γ crystallin 'Greek key' domains and one ricin B-type lectin domain. Associated with tumorigenesis, AIM1L is commonly down-regulated in colorectal amilial adenomatous polyposis (FAP) polyps, compared to normal colon tissue. AIM1L is encoded by a gene that maps to human chromosome 1p36.11. As the largest human chromosome, chromosome 1 makes up approximately 8% of the human genome and contains 260 million base pairs encoding 3,000 genes. Numerous diseases are linked to chromosome 1, notably the rare aging disease Hutchinson-Gilford progeria, which is associated with Lamin A. When defective, Lamin A can accumulate in nucleus, causing characteristic nuclear blebs. The MUTYH gene is located on chromosome 1 and is partially responsible for familial adenomatous polyposis. Stickler syndrome, Parkinson's disease, Gaucher disease and Usher syndrome are also associated with chromosome 1. Aberrations in chromosome 1 exist in a variety of cancers, including head and neck cancer, malignant melanoma and multiple myeloma.

REFERENCES

- Teichmann, U., et al. 1998. Cloning and tissue expression of the mouse ortholog of AIM1, a β/γ -crystallin superfamily member. *Mamm. Genome* 9: 715-720.
- Weise, A., et al. 2005. New insights into the evolution of chromosome 1. *Cytogenet. Genome Res.* 108: 217-222.
- Marzin, Y., et al. 2006. Chromosome 1 abnormalities in multiple myeloma. *Anticancer Res.* 26: 953-959.
- Gregory, S.G., et al. 2006. The DNA sequence and biological annotation of human chromosome 1. *Nature* 441: 315-321.
- McClintock, D., et al. 2006. Hutchinson-Gilford progeria mutant lamin A primarily targets human vascular cells as detected by an anti-Lamin A G608G antibody. *Proc. Natl. Acad. Sci. USA* 103: 2154-2159.
- Bowden, N.A., et al. 2007. Gene expression profiling in familial adenomatous polyposis adenomas and desmoid disease. *Hered. Cancer Clin. Pract.* 5: 79-96.
- Jelinsky, S.A., et al. 2008. Molecular analysis of the vaginal response to estrogens in the ovariectomized rat and postmenopausal woman. *BMC Med. Genomics* 1: 27.
- Kirpich, I.A., et al. 2010. Integrated hepatic transcriptome and proteome analysis of mice with high-fat diet-induced nonalcoholic fatty liver disease. *J. Nutr. Biochem.* 22: 38-45.

CHROMOSOMAL LOCATION

Genetic locus: AIM1L (human) mapping to 1p36.11; Aim1l (mouse) mapping to 4 D3.

SOURCE

AIM1L (N-16) is an affinity purified goat polyclonal antibody raised against a peptide mapping near the N-terminus of AIM1L of human origin.

PRODUCT

Each vial contains 200 μ g IgG in 1.0 ml of PBS with < 0.1% sodium azide and 0.1% gelatin.

Blocking peptide available for competition studies, sc-241795 P, (100 μ g peptide in 0.5 ml PBS containing < 0.1% sodium azide and 0.2% BSA).

APPLICATIONS

AIM1L (N-16) is recommended for detection of AIM1L of mouse, rat and human origin by Western Blotting (starting dilution 1:200, dilution range 1:100-1:1000), immunofluorescence (starting dilution 1:50, dilution range 1:50-1:500) and solid phase ELISA (starting dilution 1:30, dilution range 1:30-1:3000); non cross-reactive with AIM1 or AIM2.

AIM1L (N-16) is also recommended for detection of AIM1L in additional species, including bovine.

Suitable for use as control antibody for AIM1L siRNA (h): sc-78904, AIM1L siRNA (m): sc-140967, AIM1L shRNA Plasmid (h): sc-78904-SH, AIM1L shRNA Plasmid (m): sc-140967-SH, AIM1L shRNA (h) Lentiviral Particles: sc-78904-V and AIM1L shRNA (m) Lentiviral Particles: sc-140967-V.

Molecular Weight of AIM1L: 69 kDa.

RECOMMENDED SECONDARY REAGENTS

To ensure optimal results, the following support (secondary) reagents are recommended: 1) Western Blotting: use donkey anti-goat IgG-HRP: sc-2020 (dilution range: 1:2000-1:100,000) or Cruz Marker™ compatible donkey anti-goat IgG-HRP: sc-2033 (dilution range: 1:2000-1:5000), Cruz Marker™ Molecular Weight Standards: sc-2035, TBS Blotto A Blocking Reagent: sc-2333 and Western Blotting Luminol Reagent: sc-2048. 2) Immunofluorescence: use donkey anti-goat IgG-FITC: sc-2024 (dilution range: 1:100-1:400) or donkey anti-goat IgG-TR: sc-2783 (dilution range: 1:100-1:400) with UltraCruz™ Mounting Medium: sc-24941.

STORAGE

Store at 4° C, ****DO NOT FREEZE****. Stable for one year from the date of shipment. Non-hazardous. No MSDS required.

RESEARCH USE

For research use only, not for use in diagnostic procedures.

PROTOCOLS

See our web site at www.scbt.com or our catalog for detailed protocols and support products.