

FADS3 (P-16): sc-242650

BACKGROUND

FADS3 (fatty acid desaturase 3), also known as cytochrome $\beta 5$ -related protein (CYB5RP), is a 445 amino acid protein belonging to the fatty acid desaturase family. Fatty acid desaturase proteins introduce double bonds between defined carbons of the fatty acyl chain to cause unsaturation of fatty acids. Fatty acid desaturase family members are fusion products composed of two conserved histidine motifs: an N-terminal cytochrome $\beta 5$ heme-binding domain and a C-terminal multiple membrane-spanning desaturase portion. The cytochrome $\beta 5$ heme-binding domain may contain an active site and may be responsible for metal ion binding. FADS3 is a multi-pass membrane protein localized to the endoplasmic reticulum and can be found in heart, uterus, lung, brainstem and liver.

REFERENCES

1. Stöhr, H., Marquardt, A., Rivera, A., Cooper, P.R., Nowak, N.J., Shows, T.B., Gerhard, D.S. and Weber, B.H. 1998. A gene map of the Best's vitelliform macular dystrophy region in chromosome 11q12-q13.1. *Genome Res.* 8: 48-56.
2. Marquardt, A., Stöhr, H., White, K. and Weber, B.H. 2000. cDNA cloning, genomic structure, and chromosomal localization of three members of the human fatty acid desaturase family. *Genomics* 66: 175-183.
3. Nakamura, M.T. and Nara, T.Y. 2004. Structure, function, and dietary regulation of $\delta 6$, $\delta 5$, and $\delta 9$ desaturases. *Annu. Rev. Nutr.* 24: 345-376.
4. Online Mendelian Inheritance in Man, OMIM™. 2005. Johns Hopkins University, Baltimore, MD. MIM Number: 606150. World Wide Web URL: <http://www.ncbi.nlm.nih.gov/omim/>
5. Taylor, T.D., Noguchi, H., Totoki, Y., Toyoda, A., Kuroki, Y., Dewar, K., Lloyd, C., Itoh, T., Takeda, T., Kim, D.W., She, X., Barlow, K.F., Bloom, T., Bruford, E., Chang, J.L., Cuomo, C.A., Eichler, E., FitzGerald, M.G., Jaffe, D.B., LaButti, K., Nicol, R., Park, H.S., Seaman, C., Sougnez, C., Yang, X., Zimmer, A.R., Zody, M.C., Birren, B.W., Nusbaum, C., Fujiyama, A., Hattori, M., Rogers, J., Lander, E.S. and Sakaki, Y. 2006. Human chromosome 11 DNA sequence and analysis including novel gene identification. *Nature* 440: 497-500.

CHROMOSOMAL LOCATION

Genetic locus: FADS3 (human) mapping to 11q12.2.

SOURCE

FADS3 (P-16) is an affinity purified goat polyclonal antibody raised against a peptide mapping within a cytoplasmic domain of FADS3 of human origin.

PRODUCT

Each vial contains 200 μ g IgG in 1.0 ml of PBS with < 0.1% sodium azide and 0.1% gelatin.

Blocking peptide available for competition studies, sc-242650 P, (100 μ g peptide in 0.5 ml PBS containing < 0.1% sodium azide and 0.2% BSA).

RESEARCH USE

For research use only, not for use in diagnostic procedures.

APPLICATIONS

FADS3 (P-16) is recommended for detection of FADS3 of human origin by Western Blotting (starting dilution 1:200, dilution range 1:100-1:1000), immunofluorescence (starting dilution 1:50, dilution range 1:50-1:500) and solid phase ELISA (starting dilution 1:30, dilution range 1:30-1:3000); non cross-reactive with other FADS family members.

FADS3 (P-16) is also recommended for detection of FADS3 in additional species, including porcine.

Suitable for use as control antibody for FADS3 siRNA (h): sc-96983, FADS3 shRNA Plasmid (h): sc-96983-SH and FADS3 shRNA (h) Lentiviral Particles: sc-96983-V.

Molecular Weight of FADS3: 51 kDa.

RECOMMENDED SECONDARY REAGENTS

To ensure optimal results, the following support (secondary) reagents are recommended: 1) Western Blotting: use donkey anti-goat IgG-HRP: sc-2020 (dilution range: 1:2000-1:100,000) or Cruz Marker™ compatible donkey anti-goat IgG-HRP: sc-2033 (dilution range: 1:2000-1:5000), Cruz Marker™ Molecular Weight Standards: sc-2035, TBS Blotto A Blocking Reagent: sc-2333 and Western Blotting Luminol Reagent: sc-2048. 2) Immunofluorescence: use donkey anti-goat IgG-FITC: sc-2024 (dilution range: 1:100-1:400) or donkey anti-goat IgG-TR: sc-2783 (dilution range: 1:100-1:400) with UltraCruz™ Mounting Medium: sc-24941.

STORAGE

Store at 4° C, **DO NOT FREEZE**. Stable for one year from the date of shipment. Non-hazardous. No MSDS required.

PROTOCOLS

See our web site at www.scbt.com or our catalog for detailed protocols and support products.