

TTLL13 (L-15): sc-244534

BACKGROUND

A large protein group known as the tubulin tyrosine ligase-like (TTLL) family is implied to catalyze ligations of amino acids to tubulins and other substrates. Each member contains a characteristic TTL domain. TTLL13 (tubulin tyrosine ligase-like family, member 13) is a 815 amino acid protein that belongs to the tubulin—tyrosine ligase family. Containing one TTL domain, TTLL13 may act as a tubulin polyglutamylase that forms polyglutamate side chains on tubulin. TTLL13 exists as 2 alternatively spliced isoforms and is encoded by a gene located on human chromosome 15q26.1. Encoding more than 700 genes, chromosome 15 is made up of approximately 106 million base pairs and is about 3% of the human genome.

REFERENCES

1. Trichet, V., Ruault, M., Roizès and G., De Sario, A. 2000. Characterization of the human tubulin tyrosine ligase-like 1 gene (TTLL1) mapping to 22q13.1. *Gene* 257: 109-117.
2. Janke, C., Rogowski, K., Wloga, D., Regnard, C., Kajava, A.V., Strub, J.M., Temurak, N., van Dijk, J., Boucher, D., van Dorsselaer, A., Suryavanshi, S., Gaertig and J., Edde, B. 2005. Tubulin polyglutamylase enzymes are members of the TTL domain protein family. *Science* 308: 1758-1762.
3. Wloga, D., Rogowski, K., Sharma, N., Van Dijk, J., Janke, C., Edde, B., Bre, M.H., Levilliers, N., Redeker, V., Duan, J., Gorovsky, M.A., Jerka-Dziedzic and M., Gaertig, J. 2008. Glutamylation on α -tubulin is not essential but affects the assembly and functions of a subset of microtubules in *Tetrahymena thermophila*. *Eukaryotic Cell* 7: 1362-1372.
4. Ikegami, K., Horigome, D., Mukai, M., Livnat, I., MacGregor and G.R., Setou, M. 2008. TTLL10 is a protein polyglycylase that can modify nucleosome assembly protein 1. *FEBS Lett.* 582: 1129-1134.
5. Ikegami, K., Setou, M. 2009. TTLL10 can perform tubulin glycylation when co-expressed with TTLL8. *FEBS Lett.* 583: 1957-1963.
6. Kubo, T., Yanagisawa, H.A., Yagi, T., Hirono and M., Kamiya, R. 2010. Tubulin polyglutamylation regulates axonemal motility by modulating activities of inner-arm dyneins. *Curr. Biol.* 20: 441-445.
7. Wasyluk, C., Zambrano, A., Zhao, C., Brants, J., Abecassis, J., Schalken, J.A., Rogatsch, H., Schaefer, G., Pycha, A., Klocker and H., Wasyluk, B. 2010. Tubulin tyrosine ligase like 12, link to prostate cancer through tubulin post-translational modification and chromosome ploidy. *Int. J. Cancer* 127: 2542-53.

CHROMOSOMAL LOCATION

Genetic locus: TTLL13 (human) mapping to 15q26.1; Ttll13 (mouse) mapping to 7 D3.

SOURCE

TTLL13 (L-15) is an affinity purified goat polyclonal antibody raised against a peptide mapping within an internal region of TTLL13 of human origin.

RESEARCH USE

For research use only, not for use in diagnostic procedures.

PRODUCT

Each vial contains 200 μ g IgG in 1.0 ml of PBS with < 0.1% sodium azide and 0.1% gelatin.

Blocking peptide available for competition studies, sc-244534 P, (100 μ g peptide in 0.5 ml PBS containing < 0.1% sodium azide and 0.2% BSA).

APPLICATIONS

TTLL13 (L-15) is recommended for detection of TTLL13 of mouse, rat and human origin by Western Blotting (starting dilution 1:200, dilution range 1:100-1:1000), immunofluorescence (starting dilution 1:50, dilution range 1:50-1:500) and solid phase ELISA (starting dilution 1:30, dilution range 1:30-1:3000); non cross-reactive with other TTLL family members.

Suitable for use as control antibody for TTLL13 siRNA (h): sc-90050, TTLL13 siRNA (m): sc-154790, TTLL13 shRNA Plasmid (h): sc-90050-SH, TTLL13 shRNA Plasmid (m): sc-154790-SH, TTLL13 shRNA (h) Lentiviral Particles: sc-90050-V and TTLL13 shRNA (m) Lentiviral Particles: sc-154790-V.

Molecular Weight of TTLL13 isoforms: 94/53 kDa.

RECOMMENDED SECONDARY REAGENTS

To ensure optimal results, the following support (secondary) reagents are recommended: 1) Western Blotting: use donkey anti-goat IgG-HRP: sc-2020 (dilution range: 1:2000-1:100,000) or Cruz Marker™ compatible donkey anti-goat IgG-HRP: sc-2033 (dilution range: 1:2000-1:5000), Cruz Marker™ Molecular Weight Standards: sc-2035, TBS Blotto A Blocking Reagent: sc-2333 and Western Blotting Luminol Reagent: sc-2048. 2) Immunofluorescence: use donkey anti-goat IgG-FITC: sc-2024 (dilution range: 1:100-1:400) or donkey anti-goat IgG-TR: sc-2783 (dilution range: 1:100-1:400) with UltraCruz™ Mounting Medium: sc-24941.

STORAGE

Store at 4° C, **DO NOT FREEZE**. Stable for one year from the date of shipment. Non-hazardous. No MSDS required.

PROTOCOLS

See our web site at www.scbt.com or our catalog for detailed protocols and support products.