

# VWC1 (N-17): sc-244597

## BACKGROUND

VWC1, also known as VWCE (von Willebrand factor C and EGF domains), HBV X protein up-regulated gene 11 protein or URG11, is a 955 amino acid secreted protein that functions as a regulatory element during the  $\beta$ -catenin signaling pathway. Existing as two alternatively spliced isoforms, VWC1 is expressed in liver and has been shown to be upregulated in uHepG2 cells expressing Hep B xAg (hepatitis B virus X antigen). VWC1 contains four EGF-like domains, six VWFC domains, and is considered a target for prevention of hepatocarcinogenesis. The gene encoding VWC1 maps to human chromosome 11, which houses over 1,400 genes and comprises nearly 4% of the human genome. Jervell and Lange-Nielsen syndrome, Jacobsen syndrome, Niemann-Pick disease, hereditary angioedema and Smith-Lemli-Opitz syndrome are associated with defects in genes that maps to chromosome 11.

## REFERENCES

1. Fabiani, J.E., et al. 2000. Hereditary angioedema. Long-term follow-up of 88 patients. Experience of the Argentine allergy and immunology institute. *Allergol. Immunopathol. (Madr)* 28: 267-271.
2. Lian, Z., et al. 2003. Upregulated expression of a unique gene by hepatitis B x antigen promotes hepatocellular growth and tumorigenesis. *Neoplasia* 5: 229-244.
3. Lian, Z., et al. 2006. Enhanced cell survival of Hep3B cells by the hepatitis B x antigen effector, URG11, is associated with upregulation of  $\beta$ -catenin. *Hepatology* 43: 415-424.
4. Schuchman, E.H. 2007. The pathogenesis and treatment of acid sphingomyelinase-deficient Niemann-Pick disease. *J. Inher. Metab. Dis.* 30: 654-663.
5. Online Mendelian Inheritance in Man, OMIM<sup>™</sup>. 2007. Johns Hopkins University, Baltimore, MD. MIM Number: 611115. World Wide Web URL: <http://www.ncbi.nlm.nih.gov/omim/>
6. Bhuiyan, Z.A., et al. 2008. An intronic mutation leading to incomplete skipping of exon-2 in KCNQ1 rescues hearing in Jervell and Lange-Nielsen syndrome. *Prog. Biophys. Mol. Biol.* 98: 319-327.
7. Coldren, C.D., et al. 2009. Chromosomal microarray mapping suggests a role for BSX and Neurogranin in neurocognitive and behavioral defects in the 11q terminal deletion disorder (Jacobsen syndrome). *Neurogenetics* 10: 89-95.
8. Du, R., et al. 2010. URG11 mediates hypoxia-induced epithelial-to-mesenchymal transition by modulation of E-cadherin and  $\beta$ -catenin. *Biochem. Biophys. Res. Commun.* 391: 135-141.
9. Du, R., et al. 2010. URG11 promotes gastric cancer growth and invasion by activation of  $\beta$ -catenin signalling pathway. *J. Cell. Mol. Med.* 14: 621-635.

## CHROMOSOMAL LOCATION

Genetic locus: VWCE (human) mapping to 11q12.2; Vwce (mouse) mapping to 19 A.

## SOURCE

VWC1 (N-17) is an affinity purified goat polyclonal antibody raised against a peptide mapping near the N-terminus of VWC1 of human origin.

## PRODUCT

Each vial contains 200  $\mu$ g IgG in 1.0 ml of PBS with < 0.1% sodium azide and 0.1% gelatin.

Blocking peptide available for competition studies, sc-244597 P, (100  $\mu$ g peptide in 0.5 ml PBS containing < 0.1% sodium azide and 0.2% BSA).

## APPLICATIONS

VWC1 (N-17) is recommended for detection of VWC1 of mouse, rat and human origin by Western Blotting (starting dilution 1:200, dilution range 1:100-1:1000), immunofluorescence (starting dilution 1:50, dilution range 1:50-1:500) and solid phase ELISA (starting dilution 1:30, dilution range 1:30-1:3000); non cross-reactive with VWC2.

VWC1 (N-17) is also recommended for detection of VWC1 in additional species, including equine, bovine and porcine.

Suitable for use as control antibody for VWC1 siRNA (h): sc-96303, VWCE siRNA (m): sc-155238, VWC1 shRNA Plasmid (h): sc-96303-SH, VWCE shRNA Plasmid (m): sc-155238-SH, VWC1 shRNA (h) Lentiviral Particles: sc-96303-V and VWCE shRNA (m) Lentiviral Particles: sc-155238-V.

Molecular Weight of VWC1 isoforms: 99/23 kDa.

## RECOMMENDED SECONDARY REAGENTS

To ensure optimal results, the following support (secondary) reagents are recommended: 1) Western Blotting: use donkey anti-goat IgG-HRP: sc-2020 (dilution range: 1:2000-1:100,000) or Cruz Marker<sup>™</sup> compatible donkey anti-goat IgG-HRP: sc-2033 (dilution range: 1:2000-1:5000), Cruz Marker<sup>™</sup> Molecular Weight Standards: sc-2035, TBS Blotto A Blocking Reagent: sc-2333 and Western Blotting Luminol Reagent: sc-2048. 2) Immunofluorescence: use donkey anti-goat IgG-FITC: sc-2024 (dilution range: 1:100-1:400) or donkey anti-goat IgG-TR: sc-2783 (dilution range: 1:100-1:400) with UltraCruz<sup>™</sup> Mounting Medium: sc-24941.

## STORAGE

Store at 4° C, **\*\*DO NOT FREEZE\*\***. Stable for one year from the date of shipment. Non-hazardous. No MSDS required.

## RESEARCH USE

For research use only, not for use in diagnostic procedures.

## PROTOCOLS

See our web site at [www.scbt.com](http://www.scbt.com) or our catalog for detailed protocols and support products.