

11 β -HSD1L (S-17): sc-244841

BACKGROUND

11 β -HSD1, a low-affinity NADP-dependent dehydrogenase/oxoreductase, activates cortisol from cortisone and may also play a role in glucose homeostasis and in the pathogenesis of a number of disorders, such as Insulin resistance, obesity and cancer. 11 β -HSD1L (hydroxysteroid (11- β) dehydrogenase 1-like), also known as HSD3, HSD11B1L or SCDR10, is a 315 amino acid secreted protein that belongs to the short-chain dehydrogenases/reductases (SDR) family and is thought to function in a similar manner as 11 β -HSD1. Due to its similarity with 11 β -HSD1, 11 β -HSD1L may be involved in glucose regulation and in the development of glucose-related disease. Eight isoforms of 11 β -HSD1L exist due to alternative splicing events.

REFERENCES

- Tomlinson, J.W., Walker, E.A., Bujalska, I.J., Draper, N., Lavery, G.G., Cooper, M.S., Hewison, M. and Stewart, P.M. 2004. 11 β -hydroxysteroid dehydrogenase type 1: a tissue-specific regulator of glucocorticoid response. *Endocr. Rev.* 25: 831-866.
- Wake, D.J. and Walker, B.R. 2006. Inhibition of 11 β -hydroxysteroid dehydrogenase type 1 in obesity. *Endocrine* 29: 101-108.
- Ayachi, S.E., Paulmyer-Lacroix, O., Verdier, M., Alessi, M.C., Dutour, A. and Grino, M. 2006. 11 β -hydroxysteroid dehydrogenase type 1-driven cortisone reactivation regulates plasminogen activator inhibitor type 1 in adipose tissue of obese women. *J. Thromb. Haemost.* 4: 621-627.
- Komiyama, J., Nishimura, R., Lee, H.Y., Sakumoto, R., Tetsuka, M., Acosta, T.J., Skarzynski, D.J. and Okuda, K. 2008. Cortisol is a suppressor of apoptosis in bovine corpus luteum. *Biol. Reprod.* 78: 888-895.
- Uçkaya, G., Karadurmuş, N., Kutlu, O., Corakçı, A., Kizilda, S., Ural, A.U., Gül, D. and Kutlu, M. 2008. Adipose tissue 11 β -hydroxysteroid dehydrogenase type 1 and hexose-6-phosphate dehydrogenase gene expressions are increased in patients with type 2 diabetes mellitus. *Diabetes Res. Clin. Pract.* 82: S135-S140.
- Swali, A., Walker, E.A., Lavery, G.G., Tomlinson, J.W. and Stewart, P.M. 2008. 11 β -hydroxysteroid dehydrogenase type 1 regulates Insulin and glucagon secretion in pancreatic islets. *Diabetologia* 51: 2003-2011.
- Hughes, K.A., Webster, S.P. and Walker, B.R. 2008. 11 β -hydroxysteroid dehydrogenase type 1 (11 β -HSD1) inhibitors in type 2 diabetes mellitus and obesity. *Expert Opin. Investig. Drugs* 17: 481-496.
- Su, X., Vicker, N., Trusselle, M., Halem, H., Culler, M.D. and Potter, B.V. 2009. Discovery of novel inhibitors of human 11 β -hydroxysteroid dehydrogenase type 1. *Mol. Cell. Endocrinol.* 301: 169-173.

CHROMOSOMAL LOCATION

Genetic locus: HSD11B1L (human) mapping to 19p13.3.

SOURCE

11 β -HSD1L (S-17) is an affinity purified goat polyclonal antibody raised against a peptide mapping within an internal region of 11 β -HSD1L of human origin.

PRODUCT

Each vial contains 200 μ g IgG in 1.0 ml of PBS with < 0.1% sodium azide and 0.1% gelatin.

Blocking peptide available for competition studies, sc-244841 P, (100 μ g peptide in 0.5 ml PBS containing < 0.1% sodium azide and 0.2% BSA).

APPLICATIONS

11 β -HSD1L (S-17) is recommended for detection of 11 β -HSD1L of human origin by Western Blotting (starting dilution 1:200, dilution range 1:100-1:1000), immunofluorescence (starting dilution 1:50, dilution range 1:50-1:500) and solid phase ELISA (starting dilution 1:30, dilution range 1:30-1:3000); non cross-reactive with 11 β -HSD1.

11 β -HSD1L (S-17) is also recommended for detection of 11 β -HSD1L in additional species, including bovine and porcine.

Suitable for use as control antibody for 11 β -HSD1L siRNA (h): sc-97222, 11 β -HSD1L shRNA Plasmid (h): sc-97222-SH and 11 β -HSD1L shRNA (h) Lentiviral Particles: sc-97222-V.

Molecular Weight of 11 β -HSD1L: 34 kDa.

RECOMMENDED SECONDARY REAGENTS

To ensure optimal results, the following support (secondary) reagents are recommended: 1) Western Blotting: use donkey anti-goat IgG-HRP: sc-2020 (dilution range: 1:2000-1:100,000) or Cruz Marker™ compatible donkey anti-goat IgG-HRP: sc-2033 (dilution range: 1:2000-1:5000), Cruz Marker™ Molecular Weight Standards: sc-2035, TBS Blotto A Blocking Reagent: sc-2333 and Western Blotting Luminol Reagent: sc-2048. 2) Immunofluorescence: use donkey anti-goat IgG-FITC: sc-2024 (dilution range: 1:100-1:400) or donkey anti-goat IgG-TR: sc-2783 (dilution range: 1:100-1:400) with UltraCruz™ Mounting Medium: sc-24941.

STORAGE

Store at 4° C, ****DO NOT FREEZE****. Stable for one year from the date of shipment. Non-hazardous. No MSDS required.

RESEARCH USE

For research use only, not for use in diagnostic procedures.

PROTOCOLS

See our web site at www.scbt.com or our catalog for detailed protocols and support products.