

# ANKIB1 (N-18): sc-244913

## BACKGROUND

ANKIB1 (ankyrin repeat and IBR domain containing 1) is a 1,089 amino acid, ubiquitin-interacting motif (UIM) protein belonging to the RBR family. Encoded by a gene that maps to human chromosome 7q21.2, ANKIB1 is conserved in chimpanzee, canine, bovine, mouse, rat, chicken and zebrafish. ANKIB1 contains two ANK repeats, two RING-type zinc fingers, one UIM repeat and one IBR-type zinc finger. However, ANKIB1 lacks one Cys residue in its IBR-type zinc finger domain. ANKIB1 participates in protein and zinc ion binding, and may act as an E3 ubiquitin-protein ligase, or as part of the E3 complex, which receives ubiquitin from E2 ubiquitin-conjugating enzymes and moves it to substrates. Deletion of human chromosome 7q21.2 (KRIT1 locus), which includes partial loss of the ANKIB1 gene, is linked to cerebral cavernous malformations, a common autosomal dominant disorder that predisposes patients to haemorrhagic strokes and focal neurological signs.

## REFERENCES

- Jasnow, A.M., Mong, J.A., Romeo, R.D. and Pfaff, D.W. 2007. Estrogenic regulation of gene and protein expression within the amygdala of female mice. *Endocrine* 32: 271-279.
- Chadwick, R.B., Bu, L., Yu, H., Hu, Y., Wergedal, J.E., Mohan, S. and Baylink, D.J. 2007. Digit tip regrowth and differential gene expression in MRL/Mpj, DBA/2, and C57BL/6 mice. *Wound repair regen.* 15: 275-284.
- Cappelletti, V., Gariboldi, M., De Cecco, L., Toffanin, S., Reid, J.F., Lusa, L., Bajetta, E., Celio, L., Greco, M., Fabbri, A., Pierotti, M.A. and Daidone, M.G. 2008. Patterns and changes in gene expression following neo-adjuvant anti-estrogen treatment in estrogen receptor-positive breast cancer. *Endocr. Relat. Cancer* 15: 439-449.
- Thomas, J.W., Cáceres, M., Lowman, J.J., Morehouse, C.B., Short, M.E., Baldwin, E.L., Maney, D.L. and Martin, C.L. 2008. The chromosomal polymorphism linked to variation in social behavior in the white-throated sparrow (*Zonotrichia albicollis*) is a complex rearrangement and suppressor of recombination. *Genetics* 179: 1455-1468.
- Marín, I. 2009. RBR ubiquitin ligases: Diversification and streamlining in animal lineages. *J. Mol. Evol.* 69: 54-64.
- Muscarella, L.A., Guarnieri, V., Coco, M., Belli, S., Parrella, P., Pulcrano, G., Catapano, D., D'Angelo, V.A., Zelante, L. and D'Agruma, L. 2010. Small deletion at the 7q21.2 locus in a CCM family detected by real-time quantitative PCR. *J. Biomed. Biotechnol.*
- Wang, X., Gao, P., Long, M., Lin, F., Wei, J.X., Ren, J.H., Yan, L., He, T., Han, Y. and Zhang, H.Z. 2010. Essential role of cell cycle regulatory genes p21 and p27 expression in inhibition of breast cancer cells by arsenic trioxide. *Med Oncol.*

## CHROMOSOMAL LOCATION

Genetic locus: ANKIB1 (human) mapping to 7q21.2; Ankib1 (mouse) mapping to 5 A1.

## SOURCE

ANKIB1 (N-18) is an affinity purified goat polyclonal antibody raised against a peptide mapping near the N-terminus of ANKIB1 of human origin.

## PRODUCT

Each vial contains 200 µg IgG in 1.0 ml of PBS with < 0.1% sodium azide and 0.1% gelatin.

Blocking peptide available for competition studies, sc-244913 P, (100 µg peptide in 0.5 ml PBS containing < 0.1% sodium azide and 0.2% BSA).

## APPLICATIONS

ANKIB1 (N-18) is recommended for detection of ANKIB1 of mouse, rat and human origin by Western Blotting (starting dilution 1:200, dilution range 1:100-1:1000), immunofluorescence (starting dilution 1:50, dilution range 1:50-1:500) and solid phase ELISA (starting dilution 1:30, dilution range 1:30-1:3000).

ANKIB1 (N-18) is also recommended for detection of ANKIB1 in additional species, including equine, canine, bovine, porcine and avian.

Suitable for use as control antibody for ANKIB1 siRNA (h): sc-89667, ANKIB1 siRNA (m): sc-141068, ANKIB1 shRNA Plasmid (h): sc-89667-SH, ANKIB1 shRNA Plasmid (m): sc-141068-SH, ANKIB1 shRNA (h) Lentiviral Particles: sc-89667-V and ANKIB1 shRNA (m) Lentiviral Particles: sc-141068-V.

Molecular Weight of ANKIB1: 122 kDa.

## RECOMMENDED SECONDARY REAGENTS

To ensure optimal results, the following support (secondary) reagents are recommended: 1) Western Blotting: use donkey anti-goat IgG-HRP: sc-2020 (dilution range: 1:2000-1:100,000) or Cruz Marker™ compatible donkey anti-goat IgG-HRP: sc-2033 (dilution range: 1:2000-1:5000), Cruz Marker™ Molecular Weight Standards: sc-2035, TBS Blotto A Blocking Reagent: sc-2333 and Western Blotting Luminol Reagent: sc-2048. 2) Immunofluorescence: use donkey anti-goat IgG-FITC: sc-2024 (dilution range: 1:100-1:400) or donkey anti-goat IgG-TR: sc-2783 (dilution range: 1:100-1:400) with UltraCruz™ Mounting Medium: sc-24941.

## STORAGE

Store at 4° C, \*\*DO NOT FREEZE\*\*. Stable for one year from the date of shipment. Non-hazardous. No MSDS required.

## RESEARCH USE

For research use only, not for use in diagnostic procedures.

## PROTOCOLS

See our web site at [www.scbt.com](http://www.scbt.com) or our catalog for detailed protocols and support products.