

ANKRD56 (C-17): sc-244967

BACKGROUND

Ankyrins are membrane adaptor molecules that play important roles in coupling integral membrane proteins to the spectrin-based cytoskeleton network. Mutations of ankyrin genes lead to severe genetic diseases, such as fatal cardiac arrhythmias and hereditary spherocytosis. ANKRD56 (ankyrin repeat domain 56) is a 793 amino acid protein that contains 2 ANK repeats. Conserved in bovine, mouse, rat and chicken, ANKRD56 is encoded by a gene that maps to human chromosome 4, which represents approximately 6% of the human genome and contains nearly 900 genes. Notably, the Huntingtin gene, which encodes an expanded glutamine tract in cases of Huntington's disease, is located on chromosome 4. FGFR-3 is also encoded by a gene that maps to human chromosome 4 and is associated with thanatophoric dwarfism, achondroplasia, Muenke syndrome and bladder cancer. Chromosome 4 is also linked to Ellis-van Creveld syndrome, methylmalonic acidemia and polycystic kidney disease.

REFERENCES

1. Cowan, C.M. and Raymond, L.A. 2006. Selective neuronal degeneration in Huntington's disease. *Curr. Top. Dev. Biol.* 75: 25-71.
2. Chandler, R.J., Sloan, J., Fu, H., Tsai, M., Stabler, S., Allen, R., Kaestner, K.H., Kazazian, H.H. and Venditti, C.P. 2007. Metabolic phenotype of methylmalonic acidemia in mice and humans: the role of skeletal muscle. *BMC Med. Genet.* 8: 64.
3. de Frutos, C.A., Vega, S., Manzanares, M., Flores, J.M., Huertas, H., Martínez-Frías, M.L. and Nieto, M.A. 2007. Snail1 is a transcriptional effector of FGFR3 signaling during chondrogenesis and achondroplasias. *Dev. Cell* 13: 872-883.
4. Ruiz-Perez, V.L., Blair, H.J., Rodriguez-Andres, M.E., Blanco, M.J., Wilson, A., Liu, Y.N., Miles, C., Peters, H. and Goodship, J.A. 2007. Evc is a positive mediator of Ihh-regulated bone growth that localises at the base of chondrocyte cilia. *Development* 134: 2903-2912.
5. Shin, J.W., Jin, P. and Stroncek, D. 2008. Effect of leukapheresis on gene expression profiles of donor's peripheral blood mononuclear cells. *Korean J. Lab. Med.* 28: 130-135.
6. Silva, G.L., Junta, C.M., Sakamoto-Hojo, E.T., Donadi, E.A., Louzada-Junior, P. and Passos, G.A. 2009. Genetic susceptibility loci in rheumatoid arthritis establish transcriptional regulatory networks with other genes. *Ann. N.Y. Acad. Sci.* 1173: 521-537.
7. SWISS-PROT/TrEMBL (A6NEL2). World Wide Web URL: <http://www.uniprot.org/uniprot/A6NEL2>

CHROMOSOMAL LOCATION

Genetic locus: SOWAHB (human) mapping to 4q21.1.

SOURCE

ANKRD56 (C-17) is an affinity purified goat polyclonal antibody raised against a peptide mapping at the C-terminus of ANKRD56 of human origin.

PRODUCT

Each vial contains 200 µg IgG in 1.0 ml of PBS with < 0.1% sodium azide and 0.1% gelatin.

Blocking peptide available for competition studies, sc-244967 P, (100 µg peptide in 0.5 ml PBS containing < 0.1% sodium azide and 0.2% BSA).

APPLICATIONS

ANKRD56 (C-17) is recommended for detection of ANKRD56 of human origin by Western Blotting (starting dilution 1:200, dilution range 1:100-1:1000), immunofluorescence (starting dilution 1:50, dilution range 1:50-1:500) and solid phase ELISA (starting dilution 1:30, dilution range 1:30-1:3000); non cross-reactive with other ANKRD family members.

ANKRD56 (C-17) is also recommended for detection of ANKRD56 in additional species, including canine and porcine.

Suitable for use as control antibody for ANKRD56 siRNA (h): sc-89009, ANKRD56 shRNA Plasmid (h): sc-89009-SH and ANKRD56 shRNA (h) Lentiviral Particles: sc-89009-V.

Molecular Weight of ANKRD56: 86 kDa.

RECOMMENDED SECONDARY REAGENTS

To ensure optimal results, the following support (secondary) reagents are recommended: 1) Western Blotting: use donkey anti-goat IgG-HRP: sc-2020 (dilution range: 1:2000-1:100,000) or Cruz Marker™ compatible donkey anti-goat IgG-HRP: sc-2033 (dilution range: 1:2000-1:5000), Cruz Marker™ Molecular Weight Standards: sc-2035, TBS Blotto A Blocking Reagent: sc-2333 and Western Blotting Luminol Reagent: sc-2048. 2) Immunofluorescence: use donkey anti-goat IgG-FITC: sc-2024 (dilution range: 1:100-1:400) or donkey anti-goat IgG-TR: sc-2783 (dilution range: 1:100-1:400) with UltraCruz™ Mounting Medium: sc-24941.

STORAGE

Store at 4° C, ****DO NOT FREEZE****. Stable for one year from the date of shipment. Non-hazardous. No MSDS required.

RESEARCH USE

For research use only, not for use in diagnostic procedures.

PROTOCOLS

See our web site at www.scbt.com or our catalog for detailed protocols and support products.