

TSPAN11 (N-13): sc-249105

BACKGROUND

Tetraspanins are a group of hydrophobic membrane proteins that interact with a wide variety of proteins including intracellular signaling molecules, integrins and membrane receptors. Members of the tetraspanin family are characterized by the presence of four hydrophobic domains and play a role in cell development, activation, growth and motility. TSPAN11 (tetraspanin-11) is a 253 amino acid multi-pass membrane protein belonging to the tetraspanin (TM4SF) family. The gene encoding TSPAN11 maps to human chromosome 12, which encodes over 1,100 genes and comprises approximately 4.5% of the human genome. Chromosome 12 is associated with a variety of diseases and afflictions, including Hypochondrogenesis, achondrogenesis, Kniest dysplasia, Noonan syndrome and Trisomy 12p, which causes facial developmental defects and seizure disorders.

REFERENCES

1. Todd, S.C., Doctor, V.S. and Levy, S. 1998. Sequences and expression of six new members of the tetraspanin/TM4SF family. *Biochim. Biophys. Acta* 1399: 101-104.
2. Domínguez-Jiménez, C., Yáñez-Mó, M., Carreira, A., Tejedor, R., González-Amaro, R., Alvarez, V. and Sánchez-Madrid, F. 2001. Involvement of α 3 integrin/tetraspanin complexes in the angiogenic response induced by angiotensin II. *FASEB J.* 15: 1457-1459.
3. Berditchevski, F. 2001. Complexes of tetraspanins with integrins: more than meets the eye. *J. Cell Sci.* 114: 4143-4151.
4. Yokoyama, T., Nakatani, S. and Murakami, A. 2003. A case of Kniest dysplasia with retinal detachment and the mutation analysis. *Am. J. Ophthalmol.* 136: 1186-1188.
5. Chen, L., Li, X., Wang, G.L., Wang, Y., Zhu, Y.Y. and Zhu, J. 2008. Clinico-pathological significance of overexpression of TSPAN1, Ki67 and CD34 in gastric carcinoma. *Tumori* 94: 531-538.
6. Benussi, D.G., Costa, P., Zollino, M., Murdolo, M., Petix, V., Carrozzi, M. and Pecile, V. 2009. Trisomy 12p and monosomy 4p: phenotype-genotype correlation. *Genet. Test. Mol. Biomarkers* 13: 199-204.
7. Scholz, C.J., Kurzeder, C., Koretz, K., Windisch, J., Kreienberg, R., Sauer, G. and Deissler, H. 2009. Tspan-1 is a tetraspanin preferentially expressed by mucinous and endometrioid subtypes of human ovarian carcinomas. *Cancer Lett.* 275: 198-203.
8. Chen, L., Zhu, Y.Y., Zhang, X.J., Wang, G.L., Li, X.Y., He, S., Zhang, J.B. and Zhu, J.W. 2009. TSPAN1 protein expression: a significant prognostic indicator for patients with colorectal adenocarcinoma. *World J. Gastroenterol.* 15: 2270-2276.
9. Bennett, G., Sadlier, D., Doran, P.P., Macmathuna, P. and Murray, D.W. 2011. A functional and transcriptomic analysis of NET1 bioactivity in gastric cancer. *BMC Cancer* 11: 50.

STORAGE

Store at 4° C, ****DO NOT FREEZE****. Stable for one year from the date of shipment. Non-hazardous. No MSDS required.

CHROMOSOMAL LOCATION

Genetic locus: TSPAN11 (human) mapping to 12p11.21.

SOURCE

TSPAN11 (N-13) is an affinity purified goat polyclonal antibody raised against a peptide mapping at the N-terminus of TSPAN11 of human origin.

PRODUCT

Each vial contains 200 μ g IgG in 1.0 ml of PBS with < 0.1% sodium azide and 0.1% gelatin.

Blocking peptide available for competition studies, sc-249105 P, (100 μ g peptide in 0.5 ml PBS containing < 0.1% sodium azide and 0.2% BSA).

APPLICATIONS

TSPAN11 (N-13) is recommended for detection of TSPAN11 of human origin by Western Blotting (starting dilution 1:200, dilution range 1:100-1:1000), immunofluorescence (starting dilution 1:50, dilution range 1:50-1:500) and solid phase ELISA (starting dilution 1:30, dilution range 1:30-1:3000); non cross-reactive with other TSPAN family members.

Molecular Weight of TSPAN11: 28 kDa.

RECOMMENDED SECONDARY REAGENTS

To ensure optimal results, the following support (secondary) reagents are recommended: 1) Western Blotting: use donkey anti-goat IgG-HRP: sc-2020 (dilution range: 1:2000-1:100,000) or Cruz Marker™ compatible donkey anti-goat IgG-HRP: sc-2033 (dilution range: 1:2000-1:5000), Cruz Marker™ Molecular Weight Standards: sc-2035, TBS Blotto A Blocking Reagent: sc-2333 and Western Blotting Luminol Reagent: sc-2048. 2) Immunofluorescence: use donkey anti-goat IgG-FITC: sc-2024 (dilution range: 1:100-1:400) or donkey anti-goat IgG-TR: sc-2783 (dilution range: 1:100-1:400) with UltraCruz™ Mounting Medium: sc-24941.

RESEARCH USE

For research use only, not for use in diagnostic procedures.

PROTOCOLS

See our web site at www.scbt.com or our catalog for detailed protocols and support products.