

TSPYL6 (F-16): sc-249114

BACKGROUND

TSPYL6 (Testis-specific Y-encoded-like protein 6) is a 410 amino acid member of the nucleosome assembly protein (NAP) family. TSPYL6 is believed to be similar to Testis-specific Y-encoded protein 1 in form and function. The gene that encodes TSPYL6 is found on chromosome 2 which consists of 237 million bases encoding over 1,400 genes and making up approximately 8% of the human genome. A number of genetic diseases are linked to genes on chromosome 2. The lipid metabolic disorder sitosterolemia is associated with ABCG5 and ABCG8. An extremely rare recessive genetic disorder, Alström syndrome is due to mutations in the ALMS1 gene. Interestingly, chromosome 2 contains what appears to be a vestigial second centromere and vestigial telomeres which gives credence to the hypothesis that human chromosome 2 is the result of an ancient fusion of two ancestral chromosomes seen in modern form today in apes.

REFERENCES

1. Kido, T. and Lau, Y.F. 2006. The rat Tspy is preferentially expressed in elongated spermatids and interacts with the core histones. *Biochem. Biophys. Res. Commun.* 350: 56-67.
2. Kido, T. and Lau, Y.F. 2008. The human Y-encoded testis-specific protein interacts functionally with eukaryotic translation elongation factor eEF1A, a putative oncoprotein. *Int. J. Cancer* 123: 1573-1585.
3. Li, Y. and Lau, Y.F. 2008. TSPYL6 and its X-encoded homologue interact with cyclin B but exert contrasting functions on cyclin-dependent kinase 1 activities. *Oncogene* 27: 6141-6150.
4. Lau, Y.F., Li, Y. and Kido, T. 2009. Gonadoblastoma locus and the TSPYL6 gene on the human Y chromosome. *Birth Defects Res. C Embryo Today* 87: 114-122.
5. Vallangeon, B.D., Eble, J.N. and Ulbright, T.M. 2010. Macroscopic sertoli cell nodule: a study of 6 cases that presented as testicular masses. *Am. J. Surg. Pathol.* 34: 1874-1880.
6. Nickkholgh, B., Noordam, M.J., Hovingh, S.E., van Pelt, A.M., van der Veen, F. and Repping, S. 2010. Y chromosome TSPYL6 copy numbers and semen quality. *Fertil. Steril.* 94: 1744-1747.
7. Kumar, S., Gupta, R., Kumar, S. and Ali, S. 2011. Molecular mining of alleles in water buffalo *Bubalus bubalis* and characterization of the TSPYL6 and COL6A1 genes. *PLoS ONE* 6: e24958.

CHROMOSOMAL LOCATION

Genetic locus: TSPYL6 (human) mapping to 2p16.2.

SOURCE

TSPYL6 (F-16) is an affinity purified goat polyclonal antibody raised against a peptide mapping within an internal region of TSPYL6 of human origin.

STORAGE

Store at 4° C, ****DO NOT FREEZE****. Stable for one year from the date of shipment. Non-hazardous. No MSDS required.

PRODUCT

Each vial contains 200 µg IgG in 1.0 ml of PBS with < 0.1% sodium azide and 0.1% gelatin.

Blocking peptide available for competition studies, sc-249114 P, (100 µg peptide in 0.5 ml PBS containing < 0.1% sodium azide and 0.2% BSA).

APPLICATIONS

TSPYL6 (F-16) is recommended for detection of TSPYL6 of human origin by Western Blotting (starting dilution 1:200, dilution range 1:100-1:1000), immunofluorescence (starting dilution 1:50, dilution range 1:50-1:500) and solid phase ELISA (starting dilution 1:30, dilution range 1:30-1:3000); non cross-reactive with other TSPYL family members.

Suitable for use as control antibody for TSPYL6 siRNA (h): sc-94760, TSPYL6 shRNA Plasmid (h): sc-94760-SH and TSPYL6 shRNA (h) Lentiviral Particles: sc-94760-V.

Molecular Weight of TSPYL6: 35 kDa.

RECOMMENDED SECONDARY REAGENTS

To ensure optimal results, the following support (secondary) reagents are recommended: 1) Western Blotting: use donkey anti-goat IgG-HRP: sc-2020 (dilution range: 1:2000-1:100,000) or Cruz Marker™ compatible donkey anti-goat IgG-HRP: sc-2033 (dilution range: 1:2000-1:5000), Cruz Marker™ Molecular Weight Standards: sc-2035, TBS Blotto A Blocking Reagent: sc-2333 and Western Blotting Luminol Reagent: sc-2048. 2) Immunofluorescence: use donkey anti-goat IgG-FITC: sc-2024 (dilution range: 1:100-1:400) or donkey anti-goat IgG-TR: sc-2783 (dilution range: 1:100-1:400) with UltraCruz™ Mounting Medium: sc-24941.

RESEARCH USE

For research use only, not for use in diagnostic procedures.

PROTOCOLS

See our web site at www.scbt.com or our catalog for detailed protocols and support products.