

Rpb1 (γ-80): sc-25758

BACKGROUND

Commitment to cell division occurs at a point late in the G₁ phase of the cell cycle, termed Start. Passage through start requires the activation of the Cdc28 protein kinase by the cell cycle-regulated G₁ cyclins. Maximal expression of these G₁ cyclins is induced by the heterodimeric transcription factor complex composed of Swi4 (also designated Art1) and Swi6. Swi4 is the DNA-binding subunit of this complex. In addition to binding Swi4, Swi6 forms a complex with Mbp1. This complex activates S-phase cyclins and genes involved in DNA synthesis. Rpb1 is the largest subunit of the yeast RNA polymerase II. Srb4 is a basal transcription factor that is essential for the establishment of the transcription initiation apparatus. Stress factors induce transcription through the induction of various transcription factors. Yap1 activates expression in response to oxidative stress, while Msn2 and Msn4 mediate transcription via the stress response element (STRE).

REFERENCES

- Koch, T., et al. 1993. A role for the transcription factors Mbp1 and Swi4 in progression from G₁ to S phase. *Science* 261: 1551-1557.
- Koch, C., et al. 1996. Switching transcription on and off during the yeast cell cycle: Cln/Cdc28 kinases activate bound transcription factor SBF (Swi4/Swi6) at start, whereas Clb/Cdc28 kinases displace it from the promoter in G₂. *Genes Dev.* 10: 129-141.
- Siegmund, R.F. and Nasmyth, K.A. 1996. The *Saccharomyces cerevisiae* Start-specific transcription factor Swi4 interacts through the ankyrin repeats with the mitotic Clb2/Cdc28 kinase and through its conserved carboxy terminus with Swi6. *Mol. Cell. Biol.* 16: 2647-2655.
- Harrington, L.A. and Andrews, B.J. 1996. Binding to the yeast Swi4,6-dependent cell cycle box, CACGAAA, is cell cycle regulated *in vivo*. *Nucleic Acids Res.* 24: 558-565.

SOURCE

Rpb1 (γ-80) is a rabbit polyclonal antibody raised against amino acids 1-80 of Rpb1 of *Saccharomyces cerevisiae* origin.

PRODUCT

Each vial contains 200 μg IgG in 1.0 ml of PBS with < 0.1% sodium azide and 0.1% gelatin.

STORAGE

Store at 4° C, ****DO NOT FREEZE****. Stable for one year from the date of shipment. Non-hazardous. No MSDS required.

RESEARCH USE

For research use only, not for use in diagnostic procedures.

PROTOCOLS

See our web site at www.scbt.com or our catalog for detailed protocols and support products.

APPLICATIONS

Rpb1 (γ-80) is recommended for detection of Rpb1 of *Saccharomyces cerevisiae* origin by Western Blotting (starting dilution 1:200, dilution range 1:100-1:1000), immunoprecipitation [1-2 μg per 100-500 μg of total protein (1 ml of cell lysate)] and solid phase ELISA (starting dilution 1:30, dilution range 1:30-1:3000).

RECOMMENDED SECONDARY REAGENTS

To ensure optimal results, the following support (secondary) reagents are recommended: 1) Western Blotting: use goat anti-rabbit IgG-HRP: sc-2004 (dilution range: 1:2000-1:100,000) or Cruz Marker™ compatible goat anti-rabbit IgG-HRP: sc-2030 (dilution range: 1:2000-1:5000), Cruz Marker™ Molecular Weight Standards: sc-2035, TBS Blotto A Blocking Reagent: sc-2333 and Western Blotting Luminol Reagent: sc-2048. 2) Immunoprecipitation: use Protein A/G PLUS-Agarose: sc-2003 (0.5 ml agarose/2.0 ml).

SELECT PRODUCT CITATIONS

- Kizer, K.O., et al. 2005. A novel domain in Set2 mediates RNA polymerase II interaction and couples Histone H3 K36 methylation with transcript elongation. *Mol. Cell. Biol.* 25: 3305-3316.
- Radonjic, M., et al. 2005. Genome-wide analyses reveal RNA polymerase II located upstream of genes poised for rapid response upon *S. cerevisiae* stationary phase exit. *Mol. Cell* 18: 171-183.
- Larabee, R.N., et al. 2005. BUR kinase selectively regulates H3 K4 trimethylation and H2B ubiquitylation through recruitment of the PAF elongation complex. *Curr. Biol.* 15: 1487-1493.
- Xiao, T., et al. 2007. The RNA polymerase II kinase Ctk1 regulates positioning of a 5' histone methylation boundary along genes. *Mol. Cell. Biol.* 27: 721-731.
- Larabee, R.N., et al. 2007. CCR4/NOT complex associates with the proteasome and regulates histone methylation. *Proc. Natl. Acad. Sci. USA* 104: 5836-5841.
- Jensen, M.M., et al. 2008. Requirements for chromatin reassembly during transcriptional downregulation of a heat shock gene in *Saccharomyces cerevisiae*. *FEBS J.* 275: 2956-2964.
- Berretta, J., et al. 2008. A cryptic unstable transcript mediates transcriptional *trans*-silencing of the Ty1 retrotransposon in *S. cerevisiae*. *Genes Dev.* 22: 615-626.
- Rougemaille, M., et al. 2008. THO/Sub2p functions to coordinate 3'-end processing with gene-nuclear pore association. *Cell* 135: 308-321.
- Tan-Wong, S.M., et al. 2009. Gene loops function to maintain transcriptional memory through interaction with the nuclear pore complex. *Genes Dev.* 23: 2610-2624.
- Assenholt, J., et al. 2011. Implication of Ccr4-Not complex function in mRNA quality control in *Saccharomyces cerevisiae*. *RNA* 17: 1788-1794.
- Garrido-Godino, A.I., et al. 2013. Correct assembly of RNA polymerase II depends on the foot domain and is required for multiple steps of transcription in *Saccharomyces cerevisiae*. *Mol. Cell. Biol.* 33: 3611-3626.