


CETP (H-300): sc-25833



BACKGROUND

Cholesteryl ester transfer protein (CETP) is a circulating metabolic enzyme that transports cholesterol from the arteries to the liver. CETP converts cholesteryl esters from antiatherogenic high density lipoproteins (HDLs) to proatherogenic apolipoprotein B (apoB)-containing lipoproteins, including very low-, intermediate- and low-DLs. Efficient transfer and exchange of cholesteryl esters and triglycerides between the lipoprotein classes of human plasma is necessary to maintain healthy arteries.

REFERENCES

- MacLean, P.S. and Barakat, H.A. 2000. Insulin does not regulate the promoter of cholesteryl ester transfer protein (CETP) in HIRc/pCETP-CAT cells. *Mol. Cell Biochem.* 211: 1-7.
- Barter, P. 2000. CETP and atherosclerosis. *Arterioscler. Thromb. Vasc. Biol.* 20: 2029-2031.
- Wang, J., Qiang, H., Chen, D., Zhang, C. and Zhuang, Y. 2002. CETP gene mutation (D442G) increases low-density lipoprotein particle size in patients with coronary heart disease. *Clin. Chim. Acta* 322: 85-90.
- Cuchel, M., Wolfe, M.L., deLemos, A.S. and Rader, D.J. 2002. The frequency of the cholesteryl ester transfer protein-Taql B2 allele is lower in African Americans than in Caucasians. *Atherosclerosis* 163: 169-174.
- Barter, P.J., Brewer, H.B., Jr., Chapman, M.J., Hennekens, C.H., Rader, D.J. and Tall, A.R. 2003. Cholesteryl ester transfer protein: a novel target for raising HDL and inhibiting atherosclerosis. *Arterioscler. Thromb. Vasc. Biol.* 23: 160-167.
- Lira, M.E., Lloyd, D.B., Hallowell, S., Milos, P.M. and Thompson, J.F. 2004. Highly polymorphic repeat region in the CETP promoter induces unusual DNA structure. *Biochim. Biophys. Acta* 1684: 38-45.
- Zak, Z., Gautier, T., Dumont, L., Masson, D., Deckert, V., Duverneuil, L., Pais De Barros, J.P., Le Guern, N., Schneider, M., Moulin, P., Bataillard, A. and Lagrost L. 2005. Effect of cholesteryl ester transfer protein (CETP) expression on diet-induced hyperlipidemias in transgenic rats. *Atherosclerosis* 178: 279-286.

CHROMOSOMAL LOCATION

Genetic locus: CETP (human) mapping to 16q13.

SOURCE

CETP (H-300) is a rabbit polyclonal antibody raised against amino acids 194-493 mapping at the C-terminus of CETP of human origin.

PRODUCT

Each vial contains 200 µg IgG in 1.0 ml of PBS with < 0.1% sodium azide and 0.1% gelatin.

STORAGE

Store at 4° C, **DO NOT FREEZE**. Stable for one year from the date of shipment. Non-hazardous. No MSDS required.

APPLICATIONS

CETP (H-300) is recommended for detection of CETP of human origin by Western Blotting (starting dilution 1:200, dilution range 1:100-1:1000), immunoprecipitation [1–2 µg per 100–500 µg of total protein (1 ml of cell lysate)], immunofluorescence (starting dilution 1:50, dilution range 1:50-1:500) and solid phase ELISA (starting dilution 1:30, dilution range 1:30-1:3000).

Suitable for use as control antibody for CETP siRNA (h): sc-45913.


Molecular Weight of CETP: 129-154 kDa.

Positive Controls: Jurkat whole cell lysate: sc-2204.


RECOMMENDED SECONDARY REAGENTS

To ensure optimal results, the following support (secondary) reagents are recommended: 1) Western Blotting: use goat anti-rabbit IgG-HRP: sc-2004 (dilution range: 1:2000-1:100,000) or Cruz Marker™ compatible goat anti-rabbit IgG-HRP: sc-2030 (dilution range: 1:2000-1:5000), Cruz Marker™ Molecular Weight Standards: sc-2035, TBS Blotto A Blocking Reagent: sc-2333 and Western Blotting Luminol Reagent: sc-2048. 2) Immunoprecipitation: use Protein A/G PLUS-Agarose: sc-2003 (0.5 ml agarose/2.0 ml). 3) Immunofluorescence: use goat anti-rabbit IgG-FITC: sc-2012 (dilution range: 1:100-1:400) or goat anti-rabbit IgG-TR: sc-2780 (dilution range: 1:100-1:400) with UltraCruz™ Mounting Medium: sc-24941.

DATA



CETP (H-300): sc-25833. Western blot analysis of CETP expression in Jurkat whole cell lysate.



CETP (H-300): sc-25833. Western blot analysis of CETP expression in Jurkat whole cell lysate.

RESEARCH USE

For research use only, not for use in diagnostic procedures.

PROTOCOLS

See our web site at www.scbt.com or our catalog for detailed protocols and support products.