

KIN17 (C-16): sc-26955

BACKGROUND

The KIN17 protein binds to bent or curved double-stranded DNA fragments found at illegitimate recombination sites. KIN17 is ubiquitously expressed, with the highest levels of expression in muscle, heart and testis. Low doses of ionizing radiation increase KIN17 expression in mammalian cells. In keratinocytes, KIN17 expression increases during periods of hyperproliferation. UVC irradiation also increases KIN17 expression when functional XPA and XPC proteins are present. Antisense studies indicate that a decrease in KIN17 correlates with a decrease in cell proliferation and an accumulation of cells in early and mid-S phase. SV40-transformed fibroblasts overexpress KIN17, which interacts with large T antigen and reduces T antigen-dependent DNA replication. The gene encoding human KIN17 maps to chromosome 10p14.

REFERENCES

- Mazin, A., et al. 1994. KIN17, a mouse nuclear zinc finger protein that binds preferentially to curved DNA. *Nucleic Acids Res.* 22: 4335-4341.
- Mazin, A., et al. 1994. KIN17, a mouse nuclear protein, binds to bent DNA fragments that are found at illegitimate recombination junctions in mammalian cells. *Mol. Gen. Genet.* 244: 435-438.
- Biard, D.S., et al. 1997. Enhanced expression of the KIN17 protein immediately after low doses of ionizing radiation. *Radiat. Res.* 147: 442-450.
- Biard, D.S., et al. 1997. Differential expression of the HsKIN17 protein during differentiation of *in vitro* reconstructed human skin. *Arch. Dermatol. Res.* 289: 448-456.
- Kannouche, P., et al. 2000. Molecular cloning and characterization of the human KIN17 cDNA encoding a component of the UVC response that is conserved among metazoans. *Carcinogenesis* 21:1701-1710.
- Biard, D.S., et al. 2002. Ionizing radiation triggers chromatin-bound KIN17 complex formation in human cells. *J. Biol. Chem.* 277: 19156-19165.
- Miccoli, L., et al. 2002. Human KIN17 protein directly interacts with the simian virus 40 large T-antigen and inhibits DNA replication. *Cancer Res.* 62: 5425-5435.

CHROMOSOMAL LOCATION

Genetic locus: KIN (human) mapping to 10p14; Kin (mouse) mapping to 2 A1.

SOURCE

KIN17 (C-16) is an affinity purified goat polyclonal antibody raised against a peptide mapping at the C-terminus of KIN17 of human origin.

PRODUCT

Each vial contains 200 µg IgG in 1.0 ml of PBS with < 0.1% sodium azide and 0.1% gelatin.

Blocking peptide available for competition studies, sc-26955 P, (100 µg peptide in 0.5 ml PBS containing < 0.1% sodium azide and 0.2% BSA).

Available as TransCruz reagent for Gel Supershift and ChIP applications, sc-26955 X, 200 µg/0.1 ml.

APPLICATIONS

KIN17 (C-16) is recommended for detection of KIN17 of mouse, rat and human origin by Western Blotting (starting dilution 1:200, dilution range 1:100-1:1000), immunoprecipitation [1-2 µg per 100-500 µg of total protein (1 ml of cell lysate)], immunofluorescence (starting dilution 1:50, dilution range 1:50-1:500) and solid phase ELISA (starting dilution 1:30, dilution range 1:30-1:3000).

KIN17 (C-16) is also recommended for detection of KIN17 in additional species, including equine, canine and avian.

Suitable for use as control antibody for KIN17 siRNA (h): sc-45958, KIN17 siRNA (m): sc-45959, KIN17 shRNA Plasmid (h): sc-45958-SH, KIN17 shRNA Plasmid (m): sc-45959-SH, KIN17 shRNA (h) Lentiviral Particles: sc-45958-V and KIN17 shRNA (m) Lentiviral Particles: sc-45959-V.

KIN17 (C-16) X TransCruz antibody is recommended for Gel Supershift and ChIP applications.


Molecular Weight of KIN17: 45 kDa.

Positive Controls: RKO whole cell lysate: sc-36479, HCT-116 whole cell lysate: sc-364175 or Jurkat whole cell lysate: sc-2204.

RECOMMENDED SECONDARY REAGENTS

To ensure optimal results, the following support (secondary) reagents are recommended: 1) Western Blotting: use donkey anti-goat IgG-HRP: sc-2020 (dilution range: 1:2000-1:100,000) or Cruz Marker™ compatible donkey anti-goat IgG-HRP: sc-2033 (dilution range: 1:2000-1:5000), Cruz Marker™ Molecular Weight Standards: sc-2035, TBS Blotto A Blocking Reagent: sc-2333 and Western Blotting Luminol Reagent: sc-2048. 2) Immunoprecipitation: use Protein A/G PLUS-Agarose: sc-2003 (0.5 ml agarose/2.0 ml). 3) Immunofluorescence: use donkey anti-goat IgG-FITC: sc-2024 (dilution range: 1:100-1:400) or donkey anti-goat IgG-TR: sc-2783 (dilution range: 1:100-1:400) with UltraCruz™ Mounting Medium: sc-24941.

DATA



KIN17 (C-16): sc-26955. Western blot analysis of KIN17 expression in Jurkat whole cell lysate.

STORAGE

Store at 4° C, **DO NOT FREEZE**. Stable for one year from the date of shipment. Non-hazardous. No MSDS required.

RESEARCH USE

For research use only, not for use in diagnostic procedures.