

# Arnt 1 (H-6): sc-271801

## BACKGROUND

AhR, Arnt 1, Arnt 2 and BMAL1 are members of a family of transcription factors that contain a basic helix-loop-helix motif and a common "PAS" motif. The aromatic (aryl) hydrocarbon receptor, AhR, is a ligand dependent transcription factor that interacts with specific DNA sequences termed xenobiotic responsive elements (XREs) to activate several genes including CYP1A1, glutathione S-transferase Ya subunit and DT-diaphorase. The Ah receptor nuclear translocator proteins (Arnt 1 or Arnt 2) are required for ligand-dependent nuclear translocation of the Ah receptor and are also necessary for Ah Receptor binding to the XRE element. Arnt 1 (aryl hydrocarbon receptor nuclear translocator), also known as HIF1B, TANGO, bHLHe2, HIF1BETA, HIF-1 $\beta$  or ARNT, is a 789 amino acid nuclear protein that contains a basic helix-loop-helix (bHLH) domain, a PAC (PAS-associated C-terminal) domain and two PAS (PER-ARNT-SIM) domains.

## REFERENCES

1. Reyes, H., et al. 1992. Identification of the Ah Receptor nuclear translocator protein (Arnt) as a component of the DNA binding form of the Ah Receptor. *Science* 256: 1193-1195.
2. Sogawa, K., et al. 1995. Transcriptional activation domains of the Ah Receptor and Ah Receptor nuclear translocator. *J. Cancer Res. Clin. Oncol.* 121: 612-620.
3. Drutel, G., et al. 1996. Cloning and selective expression in brain and kidney of Arnt 2 homologous to the Ah Receptor nuclear translocator (Arnt). *Biochem. Biophys. Res. Commun.* 225: 333-339.
4. Hirose, K., et al. 1996. cDNA cloning and tissue-specific expression of a novel basic helix-loop-helix/PAS factor (Arnt 2) with close sequence similarity to the aryl hydrocarbon receptor nuclear translocator (Arnt). *Mol. Cell. Biol.* 16: 1706-1713.
5. Ikeda, M., et al. 1997. cDNA cloning and tissue-specific expression of a novel basic helix-loop-helix/PAS protein (BMAL1) and identification of alternatively spliced variants with alternative translation initiation site usage. *Biochem. Biophys. Res. Commun.* 233: 258-264.

## CHROMOSOMAL LOCATION

Genetic locus: ARNT (human) mapping to 1q21.3; Arnt (mouse) mapping to 3 F2.1.

## SOURCE

Arnt 1 (H-6) is a mouse monoclonal antibody specific for an epitope mapping between amino acids 759-789 at the C-terminus of Arnt 1 of human origin.

## PRODUCT

Each vial contains 200  $\mu$ g IgG<sub>2a</sub> kappa light chain in 1.0 ml of PBS with < 0.1% sodium azide and 0.1% gelatin. Also available as TransCruz reagent for Gel Supershift and ChIP applications, sc-271801 X, 200  $\mu$ g/0.1 ml.

Blocking peptide available for competition studies, sc-271801 P, (100  $\mu$ g peptide in 0.5 ml PBS containing < 0.1% sodium azide and 0.2% stabilizer protein).

## APPLICATIONS

Arnt 1 (H-6) is recommended for detection of Arnt 1 of mouse, rat and human origin by Western Blotting (starting dilution 1:100, dilution range 1:100-1:1000), immunoprecipitation [1-2  $\mu$ g per 100-500  $\mu$ g of total protein (1 ml of cell lysate)], immunofluorescence (starting dilution 1:50, dilution range 1:50-1:500) and solid phase ELISA (starting dilution 1:30, dilution range 1:30-1:3000).

Suitable for use as control antibody for Arnt 1 siRNA (h): sc-29733, Arnt 1 siRNA (m): sc-29734, Arnt 1 siRNA (r): sc-156041, Arnt 1 shRNA Plasmid (h): sc-29733-SH, Arnt 1 shRNA Plasmid (m): sc-29734-SH, Arnt 1 shRNA Plasmid (r): sc-156041-SH, Arnt 1 shRNA (h) Lentiviral Particles: sc-29733-V, Arnt 1 shRNA (m) Lentiviral Particles: sc-29734-V and Arnt 1 shRNA (r) Lentiviral Particles: sc-156041-V.

Arnt 1 (H-6) X TransCruz antibody is recommended for Gel Supershift and ChIP applications.


Molecular Weight of Arnt 1: 95 kDa.

Positive Controls: K-562 whole cell lysate: sc-2203, HeLa whole cell lysate: sc-2200 or NIH/3T3 whole cell lysate: sc-2210.


## RECOMMENDED SUPPORT REAGENTS

To ensure optimal results, the following support reagents are recommended:  
 1) Western Blotting: use m-IgG $\kappa$  BP-HRP: sc-516102 or m-IgG $\kappa$  BP-HRP (Cruz Marker): sc-516102-CM (dilution range: 1:1000-1:10000), Cruz Marker™ Molecular Weight Standards: sc-2035, UltraCruz® Blocking Reagent: sc-516214 and Western Blotting Luminol Reagent: sc-2048. 2) Immunoprecipitation: use Protein A/G PLUS-Agarose: sc-2003 (0.5 ml agarose/2.0 ml). 3) Immunofluorescence: use m-IgG $\kappa$  BP-FITC: sc-516140 or m-IgG $\kappa$  BP-PE: sc-516141 (dilution range: 1:50-1:200) with UltraCruz® Mounting Medium: sc-24941 or UltraCruz® Hard-set Mounting Medium: sc-359850.

## DATA



Arnt 1 (H-6): sc-271801. Western blot analysis of Arnt 1 expression in HeLa (A), RAW 264.7 (B) and RPE-J (C) whole cell lysates.



Arnt 1 (H-6): sc-271801. Western blot analysis of Arnt 1 expression in K-562 (A) and NIH/3T3 (B) whole cell lysates.

## STORAGE

Store at 4° C, \*\*DO NOT FREEZE\*\*. Stable for one year from the date of shipment. Non-hazardous. No MSDS required.

## RESEARCH USE

For research use only, not for use in diagnostic procedures.