



## Rad55 (yC-16): sc-27182

### BACKGROUND

The process of homologous recombination is a major DNA repair pathway that operates on DNA double-strand breaks to promote error-free repair. Central to the process of homologous recombination are the RAD52 group genes (RAD50, RAD51, RAD52, RAD54, RDH54/TID1, RAD55, RAD57, RAD59, MRE11, and XRS2), most of which were identified by their requirement for the repair of ionizing-radiation-induced DNA damage in *Saccharomyces cerevisiae*. The *Saccharomyces cerevisiae*.RAD51, RAD55, and RAD57 genes, required for genetic recombination and DNA double-strand-break repair, encode proteins homologous to one another and to the *Escherichia coli* RecA protein (sung97). The RAD55 and RAD57 encoded products exist as a stable heterodimer.

### REFERENCES

1. Lovett, S.T., et al. 1994. Sequence of the RAD55 gene of *Saccharomyces cerevisiae*. similarity of RAD55 to prokaryotic RecA and other RecA-like proteins. *Gene* 142: 103-106.
2. Hays, S.L., et al. 1995. Complex formation in yeast double-strand break repair: participation of Rad51, Rad52, Rad55, and Rad57 proteins. *Proc. Natl. Acad. Sci. U.S.A.* 92: 6925-6929.
3. Johnson, R.D., et al. 1995. Functional differences and interactions among the putative RecA homologs Rad51, Rad55, and Rad57. *Mol. Cell. Biol.* 15: 4843-4850.
4. Sung, P. 1997. Yeast Rad55 and Rad57 proteins form a heterodimer that functions with replication protein A to promote DNA strand exchange by Rad51 recombinase. *Genes Dev.* 11: 1111-1121.
5. Bashkirov, V.I., et al. 2000. DNA repair protein Rad55 is a terminal substrate of the DNA damage checkpoints. *Mol. Cell. Biol.* 20: 4393-4404.
6. Symington, L.S. 2002. Role of RAD52 epistasis group genes in homologous recombination and double-strand break repair. *Microbiol. Mol. Biol. Rev.* 66: 630-670.
7. Sugawara, N., et al. 2003. *In vivo* roles of Rad52, Rad54, and Rad55 proteins in Rad51-mediated recombination. *Mol. Cell.* 12: 209-219.
8. SWISS-PROT/TrEMBL (P38953). World Wide Web URL: <http://www.expasy.ch/sprot/sprot-top.html>

### SOURCE

Rad55 (yC-16) is an affinity purified goat polyclonal antibody raised against a peptide mapping at the C-terminus of Rad55 of *Saccharomyces cerevisiae* origin.

### STORAGE

Store at 4° C, **\*\*DO NOT FREEZE\*\***. Stable for one year from the date of shipment. Non-hazardous. No MSDS required.

### PROTOCOLS

See our web site at [www.scbt.com](http://www.scbt.com) or our catalog for detailed protocols and support products.

### PRODUCT

Each vial contains 200 µg IgG in 1.0 ml of PBS with < 0.1% sodium azide and 0.1% gelatin.

Blocking peptide available for competition studies, sc-27182 P, (100 µg peptide in 0.5 ml PBS containing < 0.1% sodium azide and 0.2% BSA).

### APPLICATIONS

Rad55 (yC-16) is recommended for detection of Rad55 of *Saccharomyces cerevisiae* origin by Western Blotting (starting dilution 1:200, dilution range 1:100-1:1000) and solid phase ELISA (starting dilution 1:30, dilution range 1:30-1:3000).

### RECOMMENDED SECONDARY REAGENTS

To ensure optimal results, the following support (secondary) reagents are recommended: 1) Western Blotting: use donkey anti-goat IgG-HRP: sc-2020 (dilution range: 1:2000-1:100,000) or Cruz Marker™ compatible donkey anti-goat IgG-HRP: sc-2033 (dilution range: 1:2000-1:5000), Cruz Marker™ Molecular Weight Standards: sc-2035, TBS Blotto A Blocking Reagent: sc-2333 and Western Blotting Luminol Reagent: sc-2048.

### RESEARCH USE

For research use only, not for use in diagnostic procedures.