



GAGA Factor (dY-20): sc-27381

BACKGROUND

Chromatin structure plays a critical role in the regulation of transcription. *Drosophila* GAGA factor (Trithorax-like protein, GAF) directs chromatin remodeling to its binding sites and is known to activate the expression of many genes in *Drosophila*.

REFERENCES

- Pagans, S., Ortiz-Lombardia, M., Espinas, M.L., Bernues, J. and Azorin, F. 2002. The *Drosophila* transcription factor tramtrack (TTK) interacts with Trithorax-like (GAGA) and represses GAGA-mediated activation. *Nucleic Acids Res.* 30: 4406-4413.
- Shimajima, T., Okada, M., Nakayama, T., Ueda, H., Okawa, K., Iwamatsu, A., Handa, H. and Hirose, S. 2003. *Drosophila* FACT contributes to Hox gene expression through physical and functional interactions with GAGA factor. *Genes Dev.* 17: 1605-1616.
- van Steensel, B., Delrow, J. and Bussemaker, H.J. 2003. Genomewide analysis of *Drosophila* GAGA factor target genes reveals context-dependent DNA binding. *Proc. Natl. Acad. Sci USA.* 100: 2580-2585.
- Schweinsberg, S., Hagstrom, K., Gohl, D., Schedl, P., Kumar, R.P., Mishra, R. and Karch, F. 2004. The enhancer-blocking activity of the Fab-7 boundary from the *Drosophila* bithorax complex requires GAGA-factor-binding sites. *Genetics.* 168:1371-1384.
- Melnikova, L., Juge, F., Gruzdeva, N., Mazur, A., Cavalli, G. and Georgiev, P. 2004. Interaction between the GAGA factor and Mod(mdg4) proteins promotes insulator bypass in *Drosophila*. *Proc. Natl. Acad. Sci. USA.* 101: 14806-14811.
- Bonet, C., Fernandez, I., Aran, X., Bernues, J., Giral, E. and Azorin, F. 2005. The GAGA protein of *Drosophila* is phosphorylated by CK2. *J. Mol. Biol.* 351: 562-572.
- Petrasccheck, M., Escher, D., Mahmoudi, T., Verrijzer, C.P., Schaffner, W. and Barberis, A. 2005. DNA looping induced by a transcriptional enhancer *in vivo*. *Nucleic Acids Res.* 33: 3743-3750.
- Wang, Y.V., Tang, H. and Gilmour, D.S. 2005. Identification *in vivo* of different rate-limiting steps associated with transcriptional activators in the presence and absence of a GAGA element. *Mol. Cell. Biol.* 25:3543-3552.

SOURCE

GAGA Factor (dY-20) is an affinity purified goat polyclonal antibody raised against a peptide mapping within an internal region of GAGA Factor of *Drosophila melanogaster* origin.

PRODUCT

Each vial contains 200 µg IgG in 1.0 ml of PBS with < 0.1% sodium azide and 0.1% gelatin.

Blocking peptide available for competition studies, sc-27381 P, (100 µg peptide in 0.5 ml PBS containing < 0.1% sodium azide and 0.2% BSA).

RESEARCH USE

For research use only, not for use in diagnostic procedures.

APPLICATIONS

GAGA Factor (dY-20) is recommended for detection of GAGA Factor of *Drosophila melanogaster* origin by Western Blotting (starting dilution 1:200, dilution range 1:100-1:1000), immunofluorescence (starting dilution 1:50, dilution range 1:50-1:500) and solid phase ELISA (starting dilution 1:30, dilution range 1:30-1:3000).

RECOMMENDED SECONDARY REAGENTS

To ensure optimal results, the following support (secondary) reagents are recommended: 1) Western Blotting: use donkey anti-goat IgG-HRP: sc-2020 (dilution range: 1:2000-1:100,000) or Cruz Marker™ compatible donkey anti-goat IgG-HRP: sc-2033 (dilution range: 1:2000-1:5000), Cruz Marker™ Molecular Weight Standards: sc-2035, TBS Blotto A Blocking Reagent: sc-2333 and Western Blotting Luminol Reagent: sc-2048. 2) Immunofluorescence: use donkey anti-goat IgG-FITC: sc-2024 (dilution range: 1:100-1:400) or donkey anti-goat IgG-TR: sc-2783 (dilution range: 1:100-1:400) with UltraCruz™ Mounting Medium: sc-24941.

STORAGE

Store at 4° C, **DO NOT FREEZE**. Stable for one year from the date of shipment. Non-hazardous. No MSDS required.

PROTOCOLS

See our web site at www.scbt.com or our catalog for detailed protocols and support products.