

# Prostasin (L-19): sc-27592

## BACKGROUND

Prostasin, a serine protease first identified in prostate tissue, activates epithelial sodium channels in a variety of tissues. Though typically a membrane-anchored protein, free Prostasin is also found in physiologic fluids and tissue culture media, indicating a mechanism for secretion from the cells as well. Aprotinin and other protease inhibitors suppress the channel-activating capacity of Prostasin, while aldosterone increases Prostasin expression and stimulates sodium uptake. In addition, DNA methylation negatively correlates with Prostasin expression in cancer cells, while enforced reexpression decreases invasiveness as well as metastatic potential, implying that Prostasin activity reflects epithelial cell physiology.

## REFERENCES

1. Liu, L., Hering-Smith, K.S., Schiro, F.R. and Hamm, L.L. 2002. Serine protease activity in M-1 cortical collecting duct cells. *Hypertension* 39: 860-864.
2. Narikiyo, T., Kitamura, K., Adachi, M., Miyoshi, T., Iwashita, K., Shiraishi, N., Nonoguchi, H., Chen, L.M., Chai, K.X., Chao, J. and Tomita, K. 2002. Regulation of Prostasin by aldosterone in the kidney. *J. Clin. Invest.* 109: 401-408.
3. Chen, L.M. and Chai, K.X. 2002. Prostasin serine protease inhibits breast cancer invasiveness and is transcriptionally regulated by promoter DNA methylation. *Int. J. Cancer* 97: 323-329.
4. Wang, C., Chao, J. and Chao, L. 2003. Adenovirus-mediated human Prostasin gene delivery is linked to increased aldosterone production and hypertension in rats. *Am. J. Physiol. Regul. Integr. Comp. Physiol.* 284: 1031-1036.
5. Iwashita, K., Kitamura, K., Narikiyo, T., Adachi, M., Shiraishi, N., Miyoshi, T., Nagano, J., Tuyen, D.G., Nonoguchi, H. and Tomita, K. 2003. Inhibition of Prostasin secretion by serine protease inhibitors in the kidney. *J. Am. Soc. Nephrol.* 14: 11-16.

## CHROMOSOMAL LOCATION

Genetic locus: Prss8 (mouse) mapping to 7 F3.

## SOURCE

Prostasin (L-19) is an affinity purified goat polyclonal antibody raised against a peptide mapping within an internal region of Prostasin of mouse origin.

## PRODUCT

Each vial contains 200 µg IgG in 1.0 ml of PBS with < 0.1% sodium azide and 0.1% gelatin.

Blocking peptide available for competition studies, sc-27592 P, (100 µg peptide in 0.5 ml PBS containing < 0.1% sodium azide and 0.2% BSA).

## STORAGE

Store at 4° C, **\*\*DO NOT FREEZE\*\***. Stable for one year from the date of shipment. Non-hazardous. No MSDS required.

## APPLICATIONS

Prostasin (L-19) is recommended for detection of precursor and mature heavy chain Prostasin of mouse, rat and, to a lesser extent, human origin by Western Blotting (starting dilution 1:200, dilution range 1:100-1:1000), immunoprecipitation [1-2 µg per 100-500 µg of total protein (1 ml of cell lysate)], immunofluorescence (starting dilution 1:50, dilution range 1:50-1:500) and solid phase ELISA (starting dilution 1:30, dilution range 1:30-1:3000).

Prostasin (L-19) is also recommended for detection of precursor and mature heavy chain Prostasin in additional species, including canine.

Suitable for use as control antibody for Prostasin siRNA (m): sc-63315, Prostasin shRNA Plasmid (m): sc-63315-SH and Prostasin shRNA (m) Lentiviral Particles: sc-63315-V.


Molecular Weight of Prostasin: 40 kDa.

Positive Controls: Prostasin (h): 293T Lysate: sc-170169.

## RECOMMENDED SECONDARY REAGENTS

To ensure optimal results, the following support (secondary) reagents are recommended: 1) Western Blotting: use donkey anti-goat IgG-HRP: sc-2020 (dilution range: 1:2000-1:100,000) or Cruz Marker™ compatible donkey anti-goat IgG-HRP: sc-2033 (dilution range: 1:2000-1:5000), Cruz Marker™ Molecular Weight Standards: sc-2035, TBS Blotto A Blocking Reagent: sc-2333 and Western Blotting Luminol Reagent: sc-2048. 2) Immunoprecipitation: use Protein A/G PLUS-Agarose: sc-2003 (0.5 ml agarose/2.0 ml). 3) Immunofluorescence: use donkey anti-goat IgG-FITC: sc-2024 (dilution range: 1:100-1:400) or donkey anti-goat IgG-TR: sc-2783 (dilution range: 1:100-1:400) with UltraCruz™ Mounting Medium: sc-24941.

## DATA



Prostasin (L-19): sc-27592. Western blot analysis of Prostasin expression in non-transfected: sc-117752 (A) and human Prostasin transfected: sc-170169 (B) 293T whole cell lysates.

## RESEARCH USE

For research use only, not for use in diagnostic procedures.

**MONOS**  
Satisfaction  
Guaranteed

Try **Prostasin (C-10): sc-514983** or **Prostasin (2): sc-136272**, our highly recommended monoclonal alternatives to Prostasin (L-19).