

rabbit anti-goat IgG-HRP: sc-2768

BACKGROUND

Santa Cruz Biotechnology's secondary antibodies are available conjugated to either an enzyme, biotin or fluorophore for use in a variety of antibody-based applications including Western Blot, immunostaining, flow cytometry and ELISA. Secondary antibodies are commonly affinity purified against immobilized whole IgG isotypes, including IgG₁, IgG_{2a}, IgG_{2b}, IgG₃ and IgG₄. Santa Cruz Biotechnology offers a wide selection of secondary antibodies, which are used in conjunction with our Cruz Marker™ molecular weight standards. We also provide specialized secondaries, such as pre-adsorbed secondary antibodies, which are pre-adsorbed with human IgG and mouse IgG for immunoglobulin-rich tissues and cells, F(ab')₂ fragment secondary antibodies that reduce non-specific secondary antibody binding to Fc receptors on the cell surface, and isotype-specific secondary antibodies against IgM, IgA and IgY.

SOURCE

rabbit anti-goat IgG-HRP is an affinity purified secondary antibody raised in rabbit against goat IgG and conjugated to HRP (horseradish peroxidase).

PRODUCT

Each vial contains 200 µg IgG in 0.5 ml of 1X PBS containing 40% glycerol.

APPLICATIONS

rabbit anti-goat IgG-HRP is recommended for detection of goat IgG by Western Blotting (starting dilution: 1:5000, dilution range 1:5000-1:10000; optimal dilution to be determined by titration).

RECOMMENDED SUPPORT PRODUCTS

- Western Blotting Luminol Reagent, for 2,000 cm² membrane area: sc-2048
- RIPA Lysis Buffer, 50 ml, cell lysis buffer with protease inhibitors: sc-24948
- Electrophoresis Sample Buffer, 2X, 25 ml, reducing buffer: sc-24945
- Complete™ Protease Inhibitor Cocktail Tablet, 20 tablets: sc-29130
- Running Buffer, 10X, 1 L, TRIS-Glycine WB running buffer, pH 8.3: sc-24949
- Towbin, with SDS, 10X, 1 L, WB transfer buffer pH 8.3: sc-24954
- Bovine Serum Albumin (BSA), 100 g, blocking/incubation agent: sc-2323
- TBS Blotto A, lyophilized powder in single-use bottle: sc-2333
- UltraCruz™ PVDF Transfer Membrane, 0.45 µm, 30 cm x 3 m roll: sc-3723
- UltraCruz™ Nitrocellulose Pure Transfer Membrane, 0.22 µm, 30 cm x 3 m roll: sc-3718
- UltraCruz™ Tissue Culture Dish, 100 mm polystyrene dish: sc-200286
- UltraCruz™ Cell Scrapers, 25 cm, sterile, 100 per case: sc-213229
- UltraCruz™ Electrophoresis Cell: sc-201625 : runs up to 10 or 15 sample by SDS – PAGE protein electrophoresis
- UltraCruz™ Autoradiography Film, Blue, 8 x 1, 100 sheets: sc-201697
- UltraCruz™ Gel Incubation Trays, 100 per pack: sc-201755 (blue), sc-201756 (green), sc-201757 (pink), sc-201758 (yellow), sc-201759 (orange)


RESEARCH USE

For research use only, not for use in diagnostic procedures.

STORAGE

Store at 4° C, **DO NOT FREEZE**. Stable for one year from the date of shipment. Non-hazardous. No MSDS required.

DATA



rabbit anti-goat IgG-HRP: sc-2768. Western blot analysis of p53 expression in Lac-Z whole cell lysate. Antibody tested: p53 (FL-393)-G: sc-6243-G.

SELECT PRODUCT CITATIONS

1. Chang, D.K., et al. 2000. Steady-state regulation of the human DNA mismatch repair system. *J. Biol. Chem.* 275: 18424-18431.
2. Sinha, S., et al. 2008. Alterations in candidate genes PHF2, FANCC, PTCH1 and XPA at chromosomal 9q22.3 region: pathological significance in early- and late-onset breast carcinoma. *Mol. Cancer* 7: 84.
3. Priyadarsini, R.V., et al. 2010. The neem limonoids azadirachtin and nimbolide induce cellcycle arrest and mitochondria-mediated apoptosis in human cervical cancer (HeLa) cells. *Free Radic. Res.* 44: 624-634.
4. Brito, M.A., et al. 2010. N-Methyl-aspartate receptor and neuronal nitric oxide synthase activation mediate bilirubin-induced neurotoxicity. *Mol. Med.* 16: 372-380.
5. Stenger, C., et al. 2010. Up-regulation of hepatic lipolysis stimulated lipoprotein receptor by leptin: a potential lever for controlling lipid clearance during the postprandial phase. *FASEB J.* 24: 4218-4228.
6. Frinchi, M., et al. 2010. FGF-2/FGFR1 neurotrophic system expression level and its basal activation do not account for the age-dependent decline of precursor cell proliferation in the subventricular zone of rat brain. *Brain Res.* 1358: 39-45.
7. Kourtidis, A., et al. 2010. An RNA interference screen identifies metabolic regulators NR1D1 and PBP as novel survival factors for breast cancer cells with the ERBB2 signature. *Cancer Res.* 70: 1783-1792.
8. Dutta, A., et al. 2010. Culture of K562 human myeloid leukemia cells in presence of fibronectin expresses and secretes MMP-9 in serum-free culture medium. *Int. J. Clin. Exp. Pathol.* 3: 288-302.
9. Mitra, S., et al. 2010. Amplification of CyclinL1 in uterine cervical carcinoma has prognostic implications. *Mol. Carcinog.* 49: 935-943.
10. Teupser, D., et al. 2011. No reduction of atherosclerosis in C-reactive protein (CRP)-deficient mice. *J. Biol. Chem.* 286: 6272-6279.
11. Wang, K., et al. 2012. Cardiac hypertrophy is positively regulated by MicroRNA miR-23a. *J. Biol. Chem.* 287: 589-599.