

mouse anti-Armenian hamster IgG-HRP: sc-2789

BACKGROUND

Santa Cruz Biotechnology's high quality, well characterized monoclonal secondary antibodies are available conjugated to either an enzyme, biotin or fluorophore for use in a variety of antibody-based applications, including Western blotting, immunostaining and flow cytometry. Santa Cruz secondary antibodies are commonly affinity purified against immobilized whole IgG isotypes, including IgG₁, IgG_{2a}, IgG_{2b}, IgG₃ and IgG₄. Monoclonal secondary antibodies are available conjugated to HRP for Western blotting (WB) and immunohistochemistry (IHC); (CM) or Cruz Marker form of HRP conjugated secondary antibodies are suitable for use with our Cruz Marker™ molecular weight standards; FITC (fluorescein isothiocyanate), PE (phycoerythrin), R (TRITC: tetramethyl rhodamine isothiocyanate), TR (Texas Red[®]), PerCP (peridinin chlorophyll protein complex), PerCP-Cy5.5 (peridinin chlorophyll protein complex with cyanin-5.5), and CruzFluor™ (488, 555 and 594) for immunofluorescence (IF), immunohistochemistry (IHC) and flow cytometry (FCM); B (biotin) for immunohistochemistry (IHC); AP (alkaline phosphatase) for Western blotting (WB); and CruzFluor[®] 680 and 790 for near-infrared (NIR) Western blotting (WB), immunofluorescence (IF), immunohistochemistry (IHC) and flow cytometry (FCM).

SOURCE

mouse anti-Armenian hamster IgG-HRP is an affinity purified secondary antibody raised in mouse against Armenian hamster IgG and conjugated to HRP (horseradish peroxidase).

PRODUCT

Each vial contains 200 µg mouse IgG in 0.5 ml of PBS containing 40% glycerol, 1% stabilizer protein and < 0.01% thimerosal.

APPLICATIONS

mouse anti-Armenian hamster IgG-HRP is recommended for detection of Armenian hamster IgG by Western Blotting (starting dilution: 1:5000, dilution range: 1:5000-1:10000) and immunohistochemical staining (starting dilution: 1:100, dilution range: 1:100-1:400). Optimal dilution to be determined by titration.

STORAGE

Store at 4° C, ****DO NOT FREEZE****. Stable for one year from the date of shipment. Non-hazardous. No MSDS required.


PROTOCOLS

See our web site at www.scbt.com for detailed protocols and support products.


RECOMMENDED SUPPORT PRODUCTS

- Western Blotting Luminol Reagent, for 2,000 cm² membrane area: sc-2048
- RIPA Lysis Buffer, 50 ml, cell lysis buffer with protease inhibitors: sc-24948
- Electrophoresis Sample Buffer, 2X, 25 ml, reducing buffer: sc-24945
- Complete™ Protease Inhibitor Cocktail Tablet, 20 tablets: sc-29130
- Running Buffer, 10X, 1 L, TRIS-Glycine WB running buffer, pH 8.3: sc-24949
- Towbin, with SDS, 10X, 1 L, WB transfer buffer pH 8.3: sc-24954
- TBS Blotto A, lyophilized powder in single-use bottle: sc-2333
- UltraCruz[®] PVDF Transfer Membrane, 0.45 µm, 30 cm x 3 m roll: sc-3723
- UltraCruz[®] Nitrocellulose Pure Transfer Membrane, 0.22 µm, 30 cm x 3 m roll: sc-3718
- UltraCruz[®] Tissue Culture Dish, 100 mm polystyrene dish: sc-200286
- UltraCruz[®] Electrophoresis Cell: sc-201625: runs up to 10 or 15 sample by SDS – PAGE protein electrophoresis
- UltraCruz[®] Autoradiography Film, Blue, 8 x 1, 100 sheets: sc-201697
- UltraCruz[®] Gel Incubation Trays, 100 per pack: sc-201755 (blue), sc-201756 (green), sc-201757 (pink), sc-201758 (yellow), sc-201759 (orange)

DATA



IL-10R α (BD-9) sc-74144. Western blot analysis of IL-10R α expression in non-transfected 293: sc-110760 (A), human IL-10R α transfected 293: sc-114689 (B) and Ramos (C) whole cell lysates. Detection reagent used: mouse anti-Armenian hamster IgG-HRP: sc-2789.



liver FBPase (N-13): sc-32430. Western blot analysis of liver FBPase expression in non-transfected: sc-117752 (A) and human liver FBPase transfected: sc-113796 (B) 293T whole cell lysates.

SELECT PRODUCT CITATIONS

- Wei, W., et al. 2014. YCF1-mediated cadmium resistance in yeast is dependent on copper metabolism and antioxidant enzymes. *Antiox. Redox Signal.* 21: 1475-89.

RESEARCH USE

For research use only, not for use in diagnostic procedures.