

# Pericentrin 1 (N-17): sc-28140

## BACKGROUND

Pericentrin 1 (nuclear pore complex protein Nup85) is a 656 amino acid member of the nucleoporin Nup85 family. Pericentrin 1 is an essential component of the nuclear pore complex (NPC) and is likely required for NPC assembly and maintenance. Bidirectional transport of macromolecules between the cytoplasm and nucleus occurs through NPCs embedded in the nuclear envelope. The Pericentrin 1 containing NPC subcomplex, Nup107-160, works to organize RNA export and in tethering NUP98/Nup98 and NUP153 to the nucleus. The Nup107-160 complex is likely required for spindle assembly during mitosis and for membrane clustering of CCL2-activated CCR2. Pericentrin 1 also acts as a regulator of microtubule nucleation by associating with  $\gamma$ -tubulin and anchors the PKA holoenzyme at the association site. The gene encoding human pericentrin 1 is found at chromosome 17q25.1.

## REFERENCES

1. Cronshaw, J.M., Krutchinsky, A.N., Zhang, W., Chait, B.T. and Matunis, M.J. 2002. Proteomic analysis of the mammalian nuclear pore complex. *J. Cell Biol.* 158: 915-927.
2. Harel, A., Orjalo, A.V., Vincent, T., Lachish-Zalait, A., Vasu, S., Shah, S., Zimmerman, E., Elbaum, M. and Forbes, D.J. 2003. Removal of a single pore subcomplex results in vertebrate nuclei devoid of nuclear pores. *Mol. Cell* 11: 853-864.
3. Loidice, I., Alves, A., Rabut, G., Van Overbeek, M., Ellenberg, J., Sibarita, J.B. and Doye, V. 2004. The entire Nup107-160 complex, including three new members, is targeted as one entity to kinetochores in mitosis. *Mol. Biol. Cell* 15: 3333-3344.
4. Terashima, Y., Onai, N., Murai, M., Enomoto, M., Poonpiriya, V., Hamada, T., Motomura, K., Suwa, M., Ezaki, T., Haga, T., Kanegasaki, S. and Matsushima, K. 2005. Pivotal function for cytoplasmic protein FROUNT in CCR2-mediated monocyte chemotaxis. *Nat. Immunol.* 6: 827-835.
5. Orjalo, A.V., Arnautov, A., Shen, Z., Boyarchuk, Y., Zeitlin, S.G., Fontoura, B., Briggs, S., Dasso, M. and Forbes, D.J. 2006. The Nup107-160 nucleoporin complex is required for correct bipolar spindle assembly. *Mol. Biol. Cell* 17: 3806-3818.
6. Daub, H., Olsen, J.V., Bairlein, M., Gnad, F., Oppermann, F.S., Körner, R., Greff, Z., Kéri, G., Stemmann, O. and Mann, M. 2008. Kinase-selective enrichment enables quantitative phosphoproteomics of the kinome across the cell cycle. *Mol. Cell* 31: 438-448.

## CHROMOSOMAL LOCATION

Genetic locus: NUP85 (human) mapping to 17q25.1; Nup85 (mouse) mapping to 11 E2.

## SOURCE

Pericentrin 1 (N-17) is an affinity purified goat polyclonal antibody raised against a peptide mapping at the N-terminus of Pericentrin 1 of human origin.

## STORAGE

Store at 4° C, **\*\*DO NOT FREEZE\*\***. Stable for one year from the date of shipment. Non-hazardous. No MSDS required.

## PRODUCT

Each vial contains 200  $\mu$ g IgG in 1.0 ml of PBS with < 0.1% sodium azide and 0.1% gelatin.

Blocking peptide available for competition studies, sc-28140 P, (100  $\mu$ g peptide in 0.5 ml PBS containing < 0.1% sodium azide and 0.2% BSA).

## APPLICATIONS

Pericentrin 1 (N-17) is recommended for detection of Pericentrin 1 of mouse, rat and human origin by Western Blotting (starting dilution 1:200, dilution range 1:100-1:1000), immunofluorescence (starting dilution 1:50, dilution range 1:50-1:500) and solid phase ELISA (starting dilution 1:30, dilution range 1:30-1:3000).

Pericentrin 1 (N-17) is also recommended for detection of Pericentrin 1 in additional species, including equine, canine, bovine and porcine.

Suitable for use as control antibody for Pericentrin 1 siRNA (h): sc-106396, Pericentrin 1 siRNA (m): sc-152168, Pericentrin 1 shRNA Plasmid (h): sc-106396-SH, Pericentrin 1 shRNA Plasmid (m): sc-152168-SH, Pericentrin 1 shRNA (h) Lentiviral Particles: sc-106396-V and Pericentrin 1 shRNA (m) Lentiviral Particles: sc-152168-V.

Molecular Weight of Pericentrin 1: 75 kDa.

## RECOMMENDED SECONDARY REAGENTS

To ensure optimal results, the following support (secondary) reagents are recommended: 1) Western Blotting: use donkey anti-goat IgG-HRP: sc-2020 (dilution range: 1:2000-1:100,000) or Cruz Marker™ compatible donkey anti-goat IgG-HRP: sc-2033 (dilution range: 1:2000-1:5000), Cruz Marker™ Molecular Weight Standards: sc-2035, TBS Blotto A Blocking Reagent: sc-2333 and Western Blotting Luminol Reagent: sc-2048. 2) Immunofluorescence: use donkey anti-goat IgG-FITC: sc-2024 (dilution range: 1:100-1:400) or donkey anti-goat IgG-TR: sc-2783 (dilution range: 1:100-1:400) with UltraCruz™ Mounting Medium: sc-24941.

## RESEARCH USE

For research use only, not for use in diagnostic procedures.

## PROTOCOLS

See our web site at [www.scbt.com](http://www.scbt.com) or our catalog for detailed protocols and support products.



Try **Pericentrin 1 (D-4): sc-376111**, our highly recommended monoclonal alternative to Pericentrin 1 (N-17).