

# CD33 (M-90): sc-28810

## BACKGROUND

CD33 is a type I transmembrane glycoprotein that is found on granulocyte and macrophage precursors in the bone marrow, but is absent from pluripotent stem cells. CD33 is also expressed on monocytes in peripheral blood. It is used as a marker to distinguish myelogenous leukemia cells from lymphoid or erythroid leukemias. CD33 may function as a sialic acid-dependent cell adhesion molecule.

## REFERENCES

1. Griffin, J.D., et al. 1984. A monoclonal antibody reactive with normal and leukemic human myeloid progenitor cells. *Leuk. Res.* 8: 521-534.
2. Favaloro, E.J., et al. 1987. Characterization of monoclonal antibodies to the human myeloid-differentiation antigen, "gp67" (CD-33). *Dis. Markers* 5: 215-225.
3. Andrews, R.G., et al. 1989. Precursors of colony-forming cells in humans can be distinguished from colony-forming cells by expression of the CD33 and CD34 antigens and light scatter properties. *J. Exp. Med.* 169: 1721-1731.
4. Handgretinger, R., et al. 1993. Expression of an early myelopoietic antigen (CD33) of a subset of human umbilical cord blood-derived natural killer cells. *Immunol. Lett.* 37: 223-228.
5. Pierelli, L., et al. 1993. Further investigations on the expression of HLA-DR, CD33 and CD13 surface antigens in purified bone marrow and peripheral blood CD34<sup>+</sup> haematopoietic progenitor cells. *Br. J. Haematol.* 84: 24-30.
6. Freeman, S.D., et al. 1995. Characterization of CD33 as a new member of the sialoadhesin family of cellular interaction molecules. *Blood* 85: 2005-2012.
7. Kelm, S., et al. 1996. The Sialoadhesins: a family of sialic-acid-dependent cellular recognition molecules within the immunoglobulin superfamily. *Glycoconj. J.* 13: 913-926.
8. Sgroi, D., et al. 1996. A single N-linked glycosylation site is implicated in the regulation of ligand recognition by the I-type lectins CD22 and CD33. *J. Biol. Chem.* 271: 18803-18809.

## CHROMOSOMAL LOCATION

Genetic locus: CD33 (human) mapping to 19q13.41; Cd33 (mouse) mapping to 7 B4.

## SOURCE

CD33 (M-90) is a rabbit polyclonal antibody raised against amino acids 18-107 mapping within an extracellular domain of CD33 of mouse origin.

## PRODUCT

Each vial contains 200 µg IgG in 1.0 ml of PBS with < 0.1% sodium azide and 0.1% gelatin.

## STORAGE

Store at 4° C, \*\*DO NOT FREEZE\*\*. Stable for one year from the date of shipment. Non-hazardous. No MSDS required.

## APPLICATIONS

CD33 (M-90) is recommended for detection of precursor and mature forms of CD33 of mouse, rat and, to a lesser extent, human origin by Western Blotting (starting dilution 1:200, dilution range 1:100-1:1000), immunoprecipitation [1-2 µg per 100-500 µg of total protein (1 ml of cell lysate)], immunofluorescence (starting dilution 1:50, dilution range 1:50-1:500) and solid phase ELISA (starting dilution 1:30, dilution range 1:30-1:3000); partially cross reactive with other Siglec family members.

Suitable for use as control antibody for CD33 siRNA (h): sc-42782, CD33 siRNA (m): sc-42783, CD33 shRNA Plasmid (h): sc-42782-SH, CD33 shRNA Plasmid (m): sc-42783-SH, CD33 shRNA (h) Lentiviral Particles: sc-42782-V and CD33 shRNA (m) Lentiviral Particles: sc-42783-V.


Molecular Weight of CD33: 67 kDa.

Positive Controls: AML-193 whole cell lysate: sc-364182.


## RECOMMENDED SECONDARY REAGENTS

To ensure optimal results, the following support (secondary) reagents are recommended: 1) Western Blotting: use goat anti-rabbit IgG-HRP: sc-2004 (dilution range: 1:2000-1:100,000) or Cruz Marker™ compatible goat anti-rabbit IgG-HRP: sc-2030 (dilution range: 1:2000-1:5000), Cruz Marker™ Molecular Weight Standards: sc-2035, TBS Blotto A Blocking Reagent: sc-2333 and Western Blotting Luminol Reagent: sc-2048. 2) Immunoprecipitation: use Protein A/G PLUS-Agarose: sc-2003 (0.5 ml agarose/2.0 ml). 3) Immunofluorescence: use goat anti-rabbit IgG-FITC: sc-2012 (dilution range: 1:100-1:400) or goat anti-rabbit IgG-TR: sc-2780 (dilution range: 1:100-1:400) with UltraCruz™ Mounting Medium: sc-24941.

## DATA



CD33 (M-90): sc-28810. Western blot analysis of CD33 expression in rat PBL whole cell lysate.



CD33 (M-90): sc-28810. Immunofluorescence staining of normal mouse intestine frozen section showing membrane staining.

## RESEARCH USE

For research use only, not for use in diagnostic procedures.

**MONOS**  
Satisfaction  
Guaranteed

Try **CD33 (6C5/2): sc-53199** or **CD33 (B-3): sc-514119**, our highly recommended monoclonal alternatives to CD33 (M-90). Also, for AC, HRP, FITC, PE, Alexa Fluor<sup>®</sup> 488 and Alexa Fluor<sup>®</sup> 647 conjugates, see **CD33 (6C5/2): sc-53199**.