

PKA II α reg (H-80): sc-28894

BACKGROUND

The second messenger cyclic AMP (cAMP) mediates diverse cellular responses to external signals such as proliferation, ion transport, regulation of metabolism and gene transcription by activation of the cAMP-dependent protein kinase (cAPK or PKA). Activation of PKA occurs when cAMP binds to the two regulatory subunits of the tetrameric PKA holoenzyme resulting in release of active catalytic subunits. Three catalytic (C) subunits have been identified, designated C α , C β and C γ , that each represent specific gene products. C α and C β are closely related (93% amino acid sequence similarity), whereas C γ displays 83% and 79% similarity to C α and C β , respectively. Activation of transcription upon elevation of cAMP levels results from translocation of PKA to the nucleus where it phosphorylates the transcription factor cAMP response element binding protein (CREB) on serine 133 which in turn leads to TFIIIB binding to TATA-box-binding protein TBP1, thus linking phospho-CREB to the pol II transcription initiation complex.

REFERENCES

1. Beavo, J.A., et al. 1974. Activation of protein kinase by physiological concentrations of cyclic AMP. Proc. Natl. Acad. Sci. USA 71: 3580-3583.
2. Krebs, E.G., et al. 1980. Phosphorylation and dephosphorylation of enzymes. Annu. Rev. Biochem. 48: 923-959.
3. Maldonado, F., et al. 1988. cAMP-dependent protein kinase, α -catalytic subunit. Nucleic Acids Res. 16: 8189-8190.
4. Gonzalez, G.A., et al. 1989. Cyclic AMP stimulates somatostatin gene transcription by phosphorylation of CREB at serine 133. Cell 59: 675-680.
5. Beebe, S.J., et al. 1990. cAMP-dependent protein kinase, β -catalytic subunit. Mol. Endocrinol. 4: 465-475.
6. Meinkoth, J.L., et al. 1993. Signal transduction through the cAMP-dependent protein kinase. Mol. Cell. Biochem. 127/128: 179-186.
7. Nordheim, A. 1994. CREB takes CBP to tango. Nature 370: 177-178.

CHROMOSOMAL LOCATION

Genetic locus: PRKAR2A (human) mapping to 3p21.31.

SOURCE

PKA II α reg (H-80) is a rabbit polyclonal antibody raised against amino acids 21-100 (deletion 39-54) mapping near the N-terminus of PKA II α reg of human origin.

PRODUCT

Each vial contains 200 μ g IgG in 1.0 ml of PBS with < 0.1% sodium azide and 0.1% gelatin.

STORAGE

Store at 4° C, **DO NOT FREEZE**. Stable for one year from the date of shipment. Non-hazardous. No MSDS required.

RESEARCH USE

For research use only, not for use in diagnostic procedures.

APPLICATIONS

PKA II α reg (H-80) is recommended for detection of PKA II α reg of human origin by Western Blotting (starting dilution 1:200, dilution range 1:100-1:1000), immunoprecipitation [1-2 μ g per 100-500 μ g of total protein (1 ml of cell lysate)], immunofluorescence (starting dilution 1:50, dilution range 1:50-1:500) and solid phase ELISA (starting dilution 1:30, dilution range 1:30-1:3000).

Suitable for use as control antibody for PKA II α reg siRNA (h): sc-39164, PKA II α reg shRNA Plasmid (h): sc-39164-SH and PKA II α reg shRNA (h) Lentiviral Particles: sc-39164-V.


Molecular Weight of PKA II α reg: 50 kDa.

Positive Controls: NTERA-2 cl.D1 whole cell lysate: sc-364181, HeLa whole cell lysate: sc-2200 or PC-3 cell lysate: sc-2220.


RECOMMENDED SECONDARY REAGENTS

To ensure optimal results, the following support (secondary) reagents are recommended: 1) Western Blotting: use goat anti-rabbit IgG-HRP: sc-2004 (dilution range: 1:2000-1:100,000) or Cruz Marker™ compatible goat anti-rabbit IgG-HRP: sc-2030 (dilution range: 1:2000-1:5000), Cruz Marker™ Molecular Weight Standards: sc-2035, TBS Blotto A Blocking Reagent: sc-2333 and Western Blotting Luminol Reagent: sc-2048. 2) Immunoprecipitation: use Protein A/G PLUS-Agarose: sc-2003 (0.5 ml agarose/2.0 ml). 3) Immunofluorescence: use goat anti-rabbit IgG-FITC: sc-2012 (dilution range: 1:100-1:400) or goat anti-rabbit IgG-TR: sc-2780 (dilution range: 1:100-1:400) with UltraCruz™ Mounting Medium: sc-24941.

DATA



PKA II α reg (H-80): sc-28894. Western blot analysis of PKA II α reg expression in NTERA-2 cl.D1 whole cell lysate.



PKA II α reg (H-80): sc-28894. Western blot analysis of PKA II α reg expression in PC-3 whole cell lysate.

PROTOCOLS

See our web site at www.scbt.com or our catalog for detailed protocols and support products.

MONOS
Satisfaction
Guaranteed

Try **PKA II α reg (H-12): sc-137220** or **PKA II α reg (40): sc-136262**, our highly recommended monoclonal alternatives to PKA II α reg (H-80).