

# NOXO1 (H-300): sc-292093

## BACKGROUND

NADPH oxidase (NOX) proteins include a membrane-bound flavocytochrome containing two subunits (gp91<sup>phox</sup> and p22<sup>phox</sup>) and the cytosolic proteins p47<sup>phox</sup> and p67<sup>phox</sup>. NOX activation leads to the formation of a complex that catalyzes the transfer of electrons from NADPH to molecular oxygen, therefore generating reactive oxygen species (ROS). NOXO1 (NADPH oxidase organizer 1), also designated SH3 and PX domain-containing protein 5 and Nox-organizing protein 1, is a 376 amino acid protein that targets NOX to different subcellular compartments and also targets NOX activators to NOX. Interestingly, NOXO1 is required for the synthesis of otoliths, crystalline structures of the inner ear that are involved in the perception of gravity. There are four isoforms of NOXO1 that are produced as a result of alternative splicing events.

## REFERENCES

- Bánfi, B., Clark, R.A., Steger, K. and Krause, K.H. 2003. Two novel proteins activate superoxide generation by the NADPH oxidase NOX1. *J. Biol. Chem.* 278: 3510-3513.
- Geiszt, M., Lekstrom, K., Witta, J. and Leto, T.L. 2003. Proteins homologous to p47<sup>phox</sup> and p67<sup>phox</sup> support superoxide production by NAD(P)H oxidase 1 in colon epithelial cells. *J. Biol. Chem.* 278: 20006-20012.
- Takeya, R., Ueno, N., Kami, K., Taura, M., Kohjima, M., Izaki, T., Nunoi, H. and Sumimoto, H. 2003. Novel human homologues of p47<sup>phox</sup> and p67<sup>phox</sup> participate in activation of superoxide-producing NADPH oxidases. *J. Biol. Chem.* 278: 25234-25246.
- Cheng, G. and Lambeth, J.D. 2004. NOXO1, regulation of lipid binding, localization, and activation of Nox1 by the Phox homology (PX) domain. *J. Biol. Chem.* 279: 4737-4742.
- Bánfi, B., Malgrange, B., Knisz, J., Steger, K., Dubois-Dauphin, M. and Krause, K.H. 2004. NOX3, a superoxide-generating NADPH oxidase of the inner ear. *J. Biol. Chem.* 279: 46065-46072.
- Cheng, G. and Lambeth, J.D. 2005. Alternative mRNA splice forms of NOXO1: differential tissue expression and regulation of Nox1 and Nox3. *Gene* 356: 118-126.
- Ueno, N., Takeya, R., Miyano, K., Kikuchi, H. and Sumimoto, H. 2005. The NADPH oxidase Nox3 constitutively produces superoxide in a p22<sup>phox</sup>-dependent manner: its regulation by oxidase organizers and activators. *J. Biol. Chem.* 280: 23328-23339.
- Park, H.S., Park, D. and Bae, Y.S. 2006. Molecular interaction of NADPH oxidase 1 with βPix and Nox Organizer 1. *Biochem. Biophys. Res. Commun.* 339: 985-990.
- Yamamoto, A., Kami, K., Takeya, R. and Sumimoto, H. 2007. Interaction between the SH3 domains and C-terminal proline-rich region in NADPH oxidase organizer 1 (Noxo1). *Biochem. Biophys. Res. Commun.* 352: 560-565.

## CHROMOSOMAL LOCATION

Genetic locus: NOXO1 (human) mapping to 16p13.3.

## SOURCE

NOXO1 (H-300) is a rabbit polyclonal antibody raised against amino acids 77-376 mapping at the C-terminus of NOXO1 of human origin.

## PRODUCT

Each vial contains 200 µg IgG in 1.0 ml of PBS with < 0.1% sodium azide and 0.1% gelatin.

## APPLICATIONS

NOXO1 (H-300) is recommended for detection of NOXO1 of human origin by Western Blotting (starting dilution 1:200, dilution range 1:100-1:1000), immunoprecipitation [1-2 µg per 100-500 µg of total protein (1 ml of cell lysate)], immunofluorescence (starting dilution 1:50, dilution range 1:50-1:500) and solid phase ELISA (starting dilution 1:30, dilution range 1:30-1:3000).

Suitable for use as control antibody for NOXO1 siRNA (h): sc-75949, NOXO1 shRNA Plasmid (h): sc-75949-SH and NOXO1 shRNA (h) Lentiviral Particles: sc-75949-V.

Molecular Weight of human NOXO1 isoforms γ/δ/β/α: 41 kDa.

Molecular Weight of mouse NOXO1 isoforms 1/2: 39/27 kDa.

## RECOMMENDED SECONDARY REAGENTS

To ensure optimal results, the following support (secondary) reagents are recommended: 1) Western Blotting: use goat anti-rabbit IgG-HRP: sc-2004 (dilution range: 1:2000-1:100,000) or Cruz Marker™ compatible goat anti-rabbit IgG-HRP: sc-2030 (dilution range: 1:2000-1:5000), Cruz Marker™ Molecular Weight Standards: sc-2035, TBS Blotto A Blocking Reagent: sc-2333 and Western Blotting Luminol Reagent: sc-2048. 2) Immunoprecipitation: use Protein A/G PLUS-Agarose: sc-2003 (0.5 ml agarose/2.0 ml). 3) Immunofluorescence: use goat anti-rabbit IgG-FITC: sc-2012 (dilution range: 1:100-1:400) or goat anti-rabbit IgG-TR: sc-2780 (dilution range: 1:100-1:400) with UltraCruz™ Mounting Medium: sc-24941.

## STORAGE

Store at 4° C, \*\*DO NOT FREEZE\*\*. Stable for one year from the date of shipment. Non-hazardous. No MSDS required.

## RESEARCH USE

For research use only, not for use in diagnostic procedures.

## PROTOCOLS

See our web site at [www.scbt.com](http://www.scbt.com) or our catalog for detailed protocols and support products.