

JMJD1A (H-300): sc-292213

BACKGROUND

JMJD1A (jumoni domain containing 1A), also known as TSGA (testis-specific protein A), JMJD1, KDM3A, JHDM2A (JMJC domain-containing histone demethylation protein 2A) or JHDM2A, is a member of the JHDM2 histone demethylase family of proteins that is predominantly expressed in testis. Containing one JMJC domain and a C-terminal C2HC4 zinc finger, JMJD1A functions as a mono- and dimethylation-specific demethylase, binding iron and α -ketoglutarate as cofactors and demethylating Lysine 9 of Histone H3. This suggests that JMJD1A plays a central role in the histone code and participates in nuclear hormone receptor-based transcriptional regulation. In addition, JMJD1A plays an important role in the regulation of cell growth during development and in chromatin regulation. JMJD1A directly regulates the expression of TNP1 and Protamine 1 (proteins required for the proper packaging and condensation of sperm chromatin) and, therefore, plays an essential role in spermatogenesis.

REFERENCES

- Nagase, T., et al. 1998. Prediction of the coding sequences of unidentified human genes. XI. The complete sequences of 100 new cDNA clones from brain which code for large proteins *in vitro*. DNA Res. 5: 277-286.
- Online Mendelian Inheritance in Man, OMIM™. 2002. Johns Hopkins University, Baltimore, MD. MIM Number: 611512. World Wide Web URL: <http://www.ncbi.nlm.nih.gov/omim/>
- Knebel, J., et al. 2006. Repression of transcription by TSGA/JMJD1A, a novel interaction partner of the ETS protein ER71. J. Cell. Biochem. 99: 319-329.
- Yamane, K., et al. 2006. JHDM2A, a JMJC-containing H3K9 demethylase, facilitates transcription activation by androgen receptor. Cell 125: 483-495.
- Ko, S.Y., et al. 2006. Identification of JMJD1A as a Stat3 downstream gene in mES cells. Cell Struct. Funct. 31: 53-62.
- Lockman, K., et al. 2007. The histone demethylase, JMJD1A, interacts with the Myocardin factors to regulate SMC differentiation marker gene expression. Circ. Res. 101: e115-e123.
- Okada, Y., et al. 2007. Histone demethylase JHDM2A is critical for Tnp1 and Prm1 transcription and spermatogenesis. Nature 450: 119-123.
- Loh, Y.H., et al. 2007. JMJD1A and JMJD2C Histone H3 Lys 9 demethylases regulate self-renewal in embryonic stem cells. Genes Dev. 21: 2545-2557.

CHROMOSOMAL LOCATION

Genetic locus: KDM3A (human) mapping to 2p11.2; Kdm3a (mouse) mapping to 6 C1.

SOURCE

JMJD1A (H-300) is a rabbit polyclonal antibody raised against amino acids 121-420 mapping near the N-terminus of JMJD1A of human origin.

PRODUCT

Each vial contains 200 μ g IgG1 in 1.0 ml of PBS with < 0.1% sodium azide and 0.1% gelatin.

APPLICATIONS

JMJD1A (H-300) is recommended for detection of JMJD1A of mouse, rat and human origin by Western Blotting (starting dilution 1:200, dilution range 1:100-1:1000), immunoprecipitation [1-2 μ g per 100-500 μ g of total protein (1 ml of cell lysate)], immunofluorescence (starting dilution 1:50, dilution range 1:50-1:500) and solid phase ELISA (starting dilution 1:30, dilution range 1:30-1:3000).

JMJD1A (H-300) is also recommended for detection of JMJD1A in additional species, including equine, canine, bovine and porcine.

Suitable for use as control antibody for JMJD1A siRNA (h): sc-94627, JMJD1A siRNA (m): sc-146322, JMJD1A shRNA Plasmid (h): sc-94627-SH, JMJD1A shRNA Plasmid (m): sc-146322-SH, JMJD1A shRNA (h) Lentiviral Particles: sc-94627-V and JMJD1A shRNA (m) Lentiviral Particles: sc-146322-V.

JMJD1A (H-300) X TransCruz antibody is recommended for Gel Supershift and ChIP applications.

Molecular Weight of JMJD1A monomer: 150 kDa.

Molecular Weight of JMJD1A homodimer: 300 kDa.

RECOMMENDED SECONDARY REAGENTS

To ensure optimal results, the following support (secondary) reagents are recommended: 1) Western Blotting: use goat anti-rabbit IgG-HRP: sc-2004 (dilution range: 1:2000-1:100,000) or Cruz Marker™ compatible goat anti-rabbit IgG-HRP: sc-2030 (dilution range: 1:2000-1:5000), Cruz Marker™ Molecular Weight Standards: sc-2035, TBS Blotto A Blocking Reagent: sc-2333 and Western Blotting Luminol Reagent: sc-2048. 2) Immunoprecipitation: use Protein A/G PLUS-Agarose: sc-2003 (0.5 ml agarose/2.0 ml). 3) Immunofluorescence: use goat anti-rabbit IgG-FITC: sc-2012 (dilution range: 1:100-1:400) or goat anti-rabbit IgG-TR: sc-2780 (dilution range: 1:100-1:400) with UltraCruz™ Mounting Medium: sc-24941.

STORAGE

Store at 4° C, **DO NOT FREEZE**. Stable for one year from the date of shipment. Non-hazardous. No MSDS required.

RESEARCH USE

For research use only, not for use in diagnostic procedures.



Try **JMJD1A (H-1): sc-398946** or **JMJD1A (C-6): sc-376608**, our highly recommended monoclonal alternatives to JMJD1A (H-300).