

FNDC3B (H-42): sc-292241

BACKGROUND

Adipogenesis, the process of transforming pre-adipocytes into mature fat cells, is of particular interest due to the role adipocytes play in obesity and type II diabetes. Adipocytes have been shown to affect a variety of functions, including hemostasis, angiogenesis and energy balance, by secreting hormones and bioactive peptides. The FNDC3B protein, also designated FAD104 (factor for adipocyte differentiation 104) or HCV NS5A-binding protein 37, is expressed during early adipogenesis. Belonging to the FNDC3 family of proteins, FNDC3B is a 1,204 amino acid protein that contains 9 fibronectin type-III domains. FNDC3B-deficient mice die within one day of birth, suggesting that FNDC3B is crucial for postpartum survival. Mouse embryonic fibroblasts (MEFs) with loss of FNDC3B function displayed a reduction in stress fiber formation, indicating a role for FNDC3B in cell proliferation, adhesion, spreading and migration.

REFERENCES

1. Gregoire, F.M., Smas, C.M. and Sul, H.S. 1998. Understanding adipocyte differentiation. *Physiol. Rev.* 78: 783-809.
2. Rosen, E.D. 2002. The molecular control of adipogenesis, with special reference to lymphatic pathology. *Ann. N.Y. Acad. Sci.* 979: 143-158.
3. Tominaga, K., Kondo, C., Johmura, Y., Nishizuka, M. and Imagawa, M. 2004. The novel gene fad104, containing a fibronectin type III domain, has a significant role in adipogenesis. *FEBS Lett.* 577: 49-54.
4. Tominaga, K., Johmura, Y., Nishizuka, M. and Imagawa, M. 2004. Fad24, a mammalian homolog of Noc3p, is a positive regulator in adipocyte differentiation. *J. Cell Sci.* 117: 6217-6226.
5. Ailhaud, G. 2006. Adipose tissue as a secretory organ: from adipogenesis to the metabolic syndrome. *C. R. Biol.* 329: 570-577.
6. Johmura, Y. 2007. Characterization of novel genes regulating adipocyte differentiation. *Yakugaku Zasshi* 127: 135-142.
7. Nishizuka, M., Kishimoto, K., Kato, A., Ikawa, M., Okabe, M., Sato, R., Niida, H., Nakanishi, M., Osada, S. and Imagawa, M. 2008. Disruption of the novel gene fad104 causes rapid postnatal death and attenuation of cell proliferation, adhesion, spreading and migration. *Exp. Cell Res.* 315: 809-819.
8. Hishida, T., Eguchi, T., Osada, S., Nishizuka, M. and Imagawa, M. 2008. A novel gene, fad49, plays a crucial role in the immediate early stage of adipocyte differentiation via involvement in mitotic clonal expansion. *FEBS J.* 275: 5576-5588.
9. Szeliga, M., Obara-Michlewska, M., Matyja, E., Lazarczyk, M., Lobo, C., Hilgier, W., Alonso, F.J., Márquez, J. and Albrecht, J. 2008. Transfection with liver-type glutaminase cDNA alters gene expression and reduces survival, migration and proliferation of T98G glioma cells. *Glia* 57: 1014-1023

CHROMOSOMAL LOCATION

Genetic locus: FNDC3B (human) mapping to 3q26.31; Fndc3b (mouse) mapping to 3 A3.

RESEARCH USE

For research use only, not for use in diagnostic procedures.

SOURCE

FNDC3B (H-42) is a rabbit polyclonal antibody raised against amino acids 1-42 mapping at the N-terminus of FNDC3B of human origin.

PRODUCT

Each vial contains 200 µg IgG in 1.0 ml of PBS with < 0.1% sodium azide and 0.1% gelatin.

APPLICATIONS

FNDC3B (H-42) is recommended for detection of FNDC3B of mouse, rat and human origin by Western Blotting (starting dilution 1:200, dilution range 1:100-1:1000), immunoprecipitation [1-2 µg per 100-500 µg of total protein (1 ml of cell lysate)], immunofluorescence (starting dilution 1:50, dilution range 1:50-1:500) and solid phase ELISA (starting dilution 1:30, dilution range 1:30-1:3000); non cross-reactive with other FNDC family members.

FNDC3B (H-42) is also recommended for detection of FNDC3B in additional species, including equine, canine, bovine and avian.

Suitable for use as control antibody for FNDC3B siRNA (h): sc-78339, FNDC3B siRNA (m): sc-145212, FNDC3B shRNA Plasmid (h): sc-78339-SH, FNDC3B shRNA Plasmid (m): sc-145212-SH, FNDC3B shRNA (h) Lentiviral Particles: sc-78339-V and FNDC3B shRNA (m) Lentiviral Particles: sc-145212-V.

Molecular Weight of FNDC3B isoforms: 133/70/8 kDa.

Positive Controls: 3T3-L1 cell lysate: sc-2243.

RECOMMENDED SECONDARY REAGENTS

To ensure optimal results, the following support (secondary) reagents are recommended: 1) Western Blotting: use goat anti-rabbit IgG-HRP: sc-2004 (dilution range: 1:2000-1:100,000) or Cruz Marker™ compatible goat anti-rabbit IgG-HRP: sc-2030 (dilution range: 1:2000-1:5000), Cruz Marker™ Molecular Weight Standards: sc-2035, TBS Blotto A Blocking Reagent: sc-2333 and Western Blotting Luminol Reagent: sc-2048. 2) Immunoprecipitation: use Protein A/G PLUS-Agarose: sc-2003 (0.5 ml agarose/2.0 ml). 3) Immunofluorescence: use goat anti-rabbit IgG-FITC: sc-2012 (dilution range: 1:100-1:400) or goat anti-rabbit IgG-TR: sc-2780 (dilution range: 1:100-1:400) with UltraCruz™ Mounting Medium: sc-24941.

STORAGE

Store at 4° C, **DO NOT FREEZE**. Stable for one year from the date of shipment. Non-hazardous. No MSDS required.

PROTOCOLS

See our web site at www.scbt.com or our catalog for detailed protocols and support products.