

MAN2A1 (M-224): sc-292352

BACKGROUND

The α -mannosidases (designated MAN1A1, MAN1A2, MAN2A1 and MAN2A2) comprise a group of soluble proteins that localize to the endoplasmic reticulum, the Golgi apparatus or the cytoplasm. Depending on their cellular location, these proteins are involved in either the processing or the degradation of newly synthesized N-glycans. MAN2A1 (mannosidase α class 2A member 1) is a single-pass type II membrane protein that localizes to the cisternae of the Golgi and is involved in protein modification pathways. More specifically, MAN2A1 uses zinc as a cofactor to catalyze the first committed step in the formation of N-glycans, namely the hydrolysis of the terminal α -D-mannose residues in the oligosaccharide Man5(GlcNAc)₃.

REFERENCES

1. Fukuda, M.N., Masri, K.A., Dell, A., Luzzatto, L. and Moremen, K.W. 1990. Incomplete synthesis of N-glycans in congenital dyserythropoietic anemia type II caused by a defect in the gene encoding α -mannosidase II. Proc. Natl. Acad. Sci. USA 87: 7443-7447.
2. Yoshihisa, T. and Anraku, Y. 1990. A novel pathway of import of α -mannosidase, a marker enzyme of vacuolar membrane, in *S. cerevisiae*. J. Biol. Chem. 265: 22418-22425.
3. Moremen, K.W., Touster, O. and Robbins, P.W. 1991. Novel purification of the catalytic domain of Golgi α -mannosidase II. Characterization and comparison with the intact enzyme. J. Biol. Chem. 266: 16876-16885.
4. Rabouille, C., Kuntz, D.A., Lockyer, A., Watson, R., Signorelli, T., Rose, D.R., van den Heuvel, M. and Roberts, D.B. 1999. The *Drosophila* GMII gene encodes a Golgi α -mannosidase II. J. Cell Sci. 112: 3319-3330.
5. Hutchins, M.U. and Klionsky, D.J. 2001. Vacuolar localization of oligomeric α -mannosidase requires the cytoplasm to vacuole targeting and autophagy pathway components in *Saccharomyces cerevisiae*. J. Biol. Chem. 276: 20491-20498.
6. van den Elsen, J.M., Kuntz, D.A. and Rose, D.R. 2001. Structure of Golgi α -mannosidase II: a target for inhibition of growth and metastasis of cancer cells. EMBO J. 20: 3008-3017.
7. Huang, W.P. and Klionsky, D.J. 2002. Autophagy in yeast: a review of the molecular machinery. Cell Struct. Funct. 27: 409-420.
8. Thumm, M. 2002. Hitchhikers guide to the vacuole-mechanisms of cargo sequestration in the Cvt and autophagic pathways. Mol. Cell 10: 1257-1258.
9. Shintani, T., Huang, W.P., Stromhaug, P.E. and Klionsky, D.J. 2002. Mechanism of cargo selection in the cytoplasm to vacuole targeting pathway. Dev. Cell 3: 825-837.

CHROMOSOMAL LOCATION

Genetic locus: Man2a1 (mouse) mapping to 17 E1.1.

STORAGE

Store at 4° C, **DO NOT FREEZE**. Stable for one year from the date of shipment. Non-hazardous. No MSDS required.

SOURCE

MAN2A1 (M-224) is a rabbit polyclonal antibody raised against amino acids 583-806 mapping within an internal region of MAN2A1 of mouse origin.

PRODUCT

Each vial contains 200 μ g IgG in 1.0 ml of PBS with < 0.1% sodium azide and 0.1% gelatin.

APPLICATIONS

MAN2A1 (M-224) is recommended for detection of MAN2A1 of mouse and rat origin by Western Blotting (starting dilution 1:200, dilution range 1:100-1:1000), immunoprecipitation [1-2 μ g per 100-500 μ g of total protein (1 ml of cell lysate)], immunofluorescence (starting dilution 1:50, dilution range 1:50-1:500) and solid phase ELISA (starting dilution 1:30, dilution range 1:30-1:3000); non cross-reactive with other MAN family members.

Suitable for use as control antibody for MAN2A1 siRNA (m): sc-61927, MAN2A1 shRNA Plasmid (m): sc-61927-SH and MAN2A1 shRNA (m) Lentiviral Particles: sc-61927-V.

Molecular Weight of MAN2A1: 131 kDa.

RECOMMENDED SECONDARY REAGENTS

To ensure optimal results, the following support (secondary) reagents are recommended: 1) Western Blotting: use goat anti-rabbit IgG-HRP: sc-2004 (dilution range: 1:2000-1:100,000) or Cruz Marker™ compatible goat anti-rabbit IgG-HRP: sc-2030 (dilution range: 1:2000-1:5000), Cruz Marker™ Molecular Weight Standards: sc-2035, TBS Blotto A Blocking Reagent: sc-2333 and Western Blotting Luminol Reagent: sc-2048. 2) Immunoprecipitation: use Protein A/G PLUS-Agarose: sc-2003 (0.5 ml agarose/2.0 ml). 3) Immunofluorescence: use goat anti-rabbit IgG-FITC: sc-2012 (dilution range: 1:100-1:400) or goat anti-rabbit IgG-TR: sc-2780 (dilution range: 1:100-1:400) with UltraCruz™ Mounting Medium: sc-24941.

RESEARCH USE

For research use only, not for use in diagnostic procedures.

PROTOCOLS

See our web site at www.scbt.com or our catalog for detailed protocols and support products.