

APOBEC3A/B (H-89): sc-292434

BACKGROUND

APOBEC3A (apolipoprotein B mRNA editing enzyme, catalytic polypeptide-like 3A), also known as ARP3 or PHRBN is a 199 amino acid protein that belongs to the cytidine deaminase family. Expressed in keratinocytes, as well as in peripheral leukocytes, APOBEC3A may be implicated in cell growth and cell cycle control. The gene encoding APOBEC3A is one of seven related genes which map to a cluster on chromosome 22 and are thought to play a role in tumor transformation and metastasis, as well as in host-associated antiviral activity. APOBEC3B (apolipoprotein B mRNA editing enzyme, catalytic polypeptide-like 3B), a 382 amino acid protein that is expressed in spleen, testis, ovary and heart, is also encoded by a gene that maps to the oncogenic cluster on chromosome 22. APOBEC3B uses zinc to catalyze reactions that may play important roles in cell growth and development.

REFERENCES

- Jarmuz, A., Chester, A., Bayliss, J., Gisbourne, J., Dunham, I., Scott, J. and Navaratnam, N. 2002. An anthropoid-specific locus of orphan C to U RNA-editing enzymes on chromosome 22. *Genomics* 79: 285-296.
- Mariani, R., Chen, D., Schrüfelbauer, B., Navarro, F., König, R., Bollman, B., Münk, C., Nymark-McMahon, H. and Landau, N.R. 2003. Species-specific exclusion of APOBEC3G from HIV-1 virions by Vif. *Cell* 114: 21-31.
- Wedekind, J.E., Dance, G.S., Sowden, M.P. and Smith, H.C. 2003. Messenger RNA editing in mammals: new members of the APOBEC family seeking roles in the family business. *Trends Genet.* 19: 207-216.
- Bogerd, H.P., Wiegand, H.L., Doehle, B.P., Lueders, K.K. and Cullen, B.R. 2006. APOBEC3A and APOBEC3B are potent inhibitors of LTR-retrotransposon function in human cells. *Nucleic Acids Res.* 34: 89-95.
- Goila-Gaur, R., Khan, M.A., Miyagi, E., Kao, S. and Strebel, K. 2007. Targeting APOBEC3A to the viral nucleoprotein complex confers antiviral activity. *Retrovirology* 4: 61.
- Marin, M., Golem, S., Rose, K.M., Kozak, S.L. and Kabat, D. 2008. Human immunodeficiency virus type 1 Vif functionally interacts with diverse APOBEC3 cytidine deaminases and moves with them between cytoplasmic sites of mRNA metabolism. *J. Virol.* 82: 987-998.
- Vartanian, J.P., Guétard, D., Henry, M. and Wain-Hobson, S. 2008. Evidence for editing of human papillomavirus DNA by APOBEC3 in benign and pre-cancerous lesions. *Science* 320: 230-233.
- Online Mendelian Inheritance in Man, OMIM[™]. 2008. Johns Hopkins University, Baltimore, MD. MIM Number: 607109. World Wide Web URL: <http://www.ncbi.nlm.nih.gov/omim/>

CHROMOSOMAL LOCATION

Genetic locus: APOBEC3A (human) mapping to 22q13.1.

STORAGE

Store at 4° C, ****DO NOT FREEZE****. Stable for one year from the date of shipment. Non-hazardous. No MSDS required.

SOURCE

APOBEC3A/B (H-89) is a rabbit polyclonal antibody raised against amino acids 193-281 mapping within an internal region of APOBEC3B of human origin.

PRODUCT

Each vial contains 200 µg IgG in 1.0 ml of PBS with < 0.1% sodium azide and 0.1% gelatin.

APPLICATIONS

APOBEC3A/B (H-89) is recommended for detection of APOBEC3A, APOBEC3B of human origin by Western Blotting (starting dilution 1:200, dilution range 1:100-1:1000), immunoprecipitation [1-2 µg per 100-500 µg of total protein (1 ml of cell lysate)], immunofluorescence (starting dilution 1:50, dilution range 1:50-1:500) and solid phase ELISA (starting dilution 1:30, dilution range 1:30-1:3000).

Molecular Weight of APOBEC3A: 23 kDa.

Molecular Weight of APOBEC3B: 46 kDa.

RECOMMENDED SECONDARY REAGENTS

To ensure optimal results, the following support (secondary) reagents are recommended: 1) Western Blotting: use goat anti-rabbit IgG-HRP: sc-2004 (dilution range: 1:2000-1:100,000) or Cruz Marker[™] compatible goat anti-rabbit IgG-HRP: sc-2030 (dilution range: 1:2000-1:5000), Cruz Marker[™] Molecular Weight Standards: sc-2035, TBS Blotto A Blocking Reagent: sc-2333 and Western Blotting Luminol Reagent: sc-2048. 2) Immunofluorescence: use goat anti-rabbit IgG-FITC: sc-2012 (dilution range: 1:100-1:400) or goat anti-rabbit IgG-TR: sc-2780 (dilution range: 1:100-1:400) with UltraCruz[™] Mounting Medium: sc-24941.

RESEARCH USE

For research use only, not for use in diagnostic procedures.

PROTOCOLS

See our web site at www.scbt.com or our catalog for detailed protocols and support products.