

RNase Z1 (H-110): sc-292763

BACKGROUND

RNase Z1 (ribonuclease Z 1), also known as ELAC1 or D29, is a 363 amino acid nuclear protein that belongs to the RNase Z family of RNA-processing enzymes. Expressed throughout the body, with highest expression in lung, brain, placenta, heart, liver, kidney and pancreas, RNase Z1 is a zinc phosphodiesterase that plays a role in tRNA maturation by removing 3'-trailer segments from precursor tRNA. In addition, RNase Z has a variety of other zinc-dependent RNA processing activities, including generation of 3' termini on tRNA, endonucleolytic cleavage of RNA and removal of extra 3' nucleotides from unprocessed tRNA. Defects in the gene encoding RNase Z1 may be involved in the pathogenesis of prostate cancer, suggesting a possible role for RNase Z1 in carcinogenesis.

REFERENCES

- Yanaihara, N., Kohno, T., Takakura, S., Takei, K., Otsuka, A., Sunaga, N., Takahashi, M., Yamazaki, M., Tashiro, H., Fukuzumi, Y., Fujimori, Y., Hagiwara, K., Tanaka, T. and Yokota, J. 2001. Physical and transcriptional map of a 311-kb segment of chromosome 18q21, a candidate lung tumor suppressor locus. *Genomics* 72: 169-179.
- Tavtigian, S.V., Simard, J., Teng, D.H., Abtin, V., Baumgard, M., Beck, A., Camp, N.J., Carillo, A.R., Chen, Y., Dayananth, P., Desrochers, M., Dumont, M., Farnham, J.M., Frank, D., Frye, C., Ghaffari, S., Gupte, J.S., Hu, R., et al. 2001. A candidate prostate cancer susceptibility gene at chromosome 17p. *Nat. Genet.* 27: 172-180.
- Schiffer, S., Rösch, S. and Marchfelder, A. 2002. Assigning a function to a conserved group of proteins: the tRNA 3'-processing enzymes. *EMBO J.* 21: 2769-2777.
- Online Mendelian Inheritance in Man, OMIM™. 2002. Johns Hopkins University, Baltimore, MD. MIM Number: 608079. World Wide Web URL: <http://www.ncbi.nlm.nih.gov/omim/>
- Takaku, H., Minagawa, A., Takagi, M. and Nashimoto, M. 2003. A candidate prostate cancer susceptibility gene encodes tRNA 3' processing endoribonuclease. *Nucleic Acids Res.* 31: 2272-2278.
- Schilling, O., Rüggeberg, S., Vogel, A., Rittner, N., Weichert, S., Schmidt, S., Doig, S., Franz, T., Benes, V., Andrews, S.C., Baum, M. and Meyer-Klaucke, W. 2004. Characterization of an *Escherichia coli* *elaC* deletion mutant. *Biochem. Biophys. Res. Commun.* 320: 1365-1373.

CHROMOSOMAL LOCATION

Genetic locus: ELAC1 (human) mapping to 18q21.2; Elac1 (mouse) mapping to 18 E2.

SOURCE

RNase Z1 (H-110) is a rabbit polyclonal antibody raised against amino acids 161-270 mapping within an internal region of RNase Z1 of human origin.

STORAGE

Store at 4° C, **DO NOT FREEZE**. Stable for one year from the date of shipment. Non-hazardous. No MSDS required.

PRODUCT

Each vial contains 200 µg IgG in 1.0 ml of PBS with < 0.1% sodium azide and 0.1% gelatin.

Available as TransCruz reagent for Gel Supershift and ChIP applications, sc-292763 X, 200 µg/0.1 ml.

APPLICATIONS

RNase Z1 (H-110) is recommended for detection of RNase Z1 of mouse, rat and human origin by Western Blotting (starting dilution 1:200, dilution range 1:100-1:1000), immunoprecipitation [1-2 µg per 100-500 µg of total protein (1 ml of cell lysate)], immunofluorescence (starting dilution 1:50, dilution range 1:50-1:500) and solid phase ELISA (starting dilution 1:30, dilution range 1:30-1:3000).

RNase Z1 (H-110) is also recommended for detection of RNase Z1 in additional species, including equine, canine, bovine and porcine.

Suitable for use as control antibody for RNase Z1 siRNA (h): sc-76411, RNase Z1 siRNA (m): sc-144626, RNase Z1 shRNA Plasmid (h): sc-76411-SH, RNase Z1 shRNA Plasmid (m): sc-144626-SH, RNase Z1 shRNA (h) Lentiviral Particles: sc-76411-V and RNase Z1 shRNA (m) Lentiviral Particles: sc-144626-V.

RNase Z1 (H-110) X TransCruz antibody is recommended for Gel Supershift and ChIP applications.

Molecular Weight of RNase Z1: 40 kDa.

Positive Controls: PC-12 cell lysate: sc-2250.

RECOMMENDED SECONDARY REAGENTS

To ensure optimal results, the following support (secondary) reagents are recommended: 1) Western Blotting: use goat anti-rabbit IgG-HRP: sc-2004 (dilution range: 1:2000-1:100,000) or Cruz Marker™ compatible goat anti-rabbit IgG-HRP: sc-2030 (dilution range: 1:2000-1:5000), Cruz Marker™ Molecular Weight Standards: sc-2035, TBS Blotto A Blocking Reagent: sc-2333 and Western Blotting Luminol Reagent: sc-2048. 2) Immunoprecipitation: use Protein A/G PLUS-Agarose: sc-2003 (0.5 ml agarose/2.0 ml). 3) Immunofluorescence: use goat anti-rabbit IgG-FITC: sc-2012 (dilution range: 1:100-1:400) or goat anti-rabbit IgG-TR: sc-2780 (dilution range: 1:100-1:400) with UltraCruz™ Mounting Medium: sc-24941.

RESEARCH USE

For research use only, not for use in diagnostic procedures.

PROTOCOLS

See our web site at www.scbt.com or our catalog for detailed protocols and support products.