

PRRG2 (H-60): sc-292956

BACKGROUND

PRRG2 (proline-rich γ -carboxyglutamic acid protein 2), also known as PRGP2 or TMG2 (transmembrane γ -carboxyglutamic acid protein 2), is a 202 amino acid single-pass type I membrane protein that is widely expressed, with highest levels found in thyroid. PRRG2 contains an N-terminal propeptide sequence that allows for vitamin K-dependent γ -carboxylation of various glutamic acid residues, followed by an N-terminal Gla (γ -carboxy-glutamate) domain and transmembrane region. PRRG2 has a proline-rich C-terminus with PPXY and PXXP motifs. The gene encoding PRRG2 maps to human chromosome 19, which consists of over 63 million bases, houses approximately 1,400 genes and is recognized for having the greatest gene density of the human chromosomes. It is the genetic home for a number of immunoglobulin (Ig) superfamily members, including the killer cell and leukocyte Ig-like receptors, a number of ICAMs, the CEACAM and PSG family and Fc receptors (FcRs).

REFERENCES

- Olsen, A., Teglund, S., Nelson, D., Gordon, L., Copeland, A., Georgescu, A., Carrano, A. and Hammarström, S. 1994. Gene organization of the pregnancy-specific glycoprotein region on human chromosome 19: assembly and analysis of a 700-kb cosmid contig spanning the region. *Genomics* 23: 659-668.
- Kulman, J.D., Harris, J.E., Haldeman, B.A. and Davie, E.W. 1997. Primary structure and tissue distribution of two novel proline-rich γ -carboxyglutamic acid proteins. *Proc. Natl. Acad. Sci. USA* 94: 9058-9062.
- Wang, L., Lin, S.H., Wu, W.G., Kemp, B.L., Walsh, G.L., Hong, W.K. and Mao, L. 2000. C-CAM1, a candidate tumor suppressor gene, is abnormally expressed in primary lung cancers. *Clin. Cancer Res.* 6: 2988-2993.
- Trowsdale, J., Barten, R., Haude, A., Stewart, C.A., Beck, S. and Wilson, M.J. 2001. The genomic context of natural killer receptor extended gene families. *Immunol. Rev.* 181: 20-38.
- Leeb, T. and Müller, M. 2004. Comparative human-mouse-rat sequence analysis of the ICAM gene cluster on HSA 19p13.2 and a 185-kb porcine region from SSC 2q. *Gene* 343: 239-244.
- Kulman, J.D., Harris, J.E., Xie, L. and Davie, E.W. 2007. Proline-rich Gla protein 2 is a cell-surface vitamin K-dependent protein that binds to the transcriptional coactivator Yes-associated protein. *Proc. Natl. Acad. Sci. USA* 104: 8767-8772.
- Barrow, A.D. and Trowsdale, J. 2008. The extended human leukocyte receptor complex: diverse ways of modulating immune responses. *Immunol. Rev.* 224: 98-123.
- Online Mendelian Inheritance in Man, OMIM™. 2008. Johns Hopkins University, Baltimore, MD. MIM Number: 604429. World Wide Web URL: <http://www.ncbi.nlm.nih.gov/omim/>

CHROMOSOMAL LOCATION

Genetic locus: PRRG2 (human) mapping to 19q13.33; Prrg2 (mouse) mapping to 7 B4.

SOURCE

PRRG2 (H-60) is a rabbit polyclonal antibody raised against amino acids 81-140 mapping within an internal region of PRRG2 of human origin.

PRODUCT

Each vial contains 200 μ g IgG in 1.0 ml of PBS with < 0.1% sodium azide and 0.1% gelatin.

APPLICATIONS

PRRG2 (H-60) is recommended for detection of PRRG2 of mouse, rat and human origin by Western Blotting (starting dilution 1:200, dilution range 1:100-1:1000), immunoprecipitation [1-2 μ g per 100-500 μ g of total protein (1 ml of cell lysate)], immunofluorescence (starting dilution 1:50, dilution range 1:50-1:500) and solid phase ELISA (starting dilution 1:30, dilution range 1:30-1:3000).

PRRG2 (H-60) is also recommended for detection of PRRG2 in additional species, including equine, canine and porcine.

Suitable for use as control antibody for PRRG2 siRNA (h): sc-97344, PRRG2 siRNA (m): sc-152517, PRRG2 shRNA Plasmid (h): sc-97344-SH, PRRG2 shRNA Plasmid (m): sc-152517-SH, PRRG2 shRNA (h) Lentiviral Particles: sc-97344-V and PRRG2 shRNA (m) Lentiviral Particles: sc-152517-V.

Molecular Weight of PRRG2: 22 kDa.

RECOMMENDED SECONDARY REAGENTS

To ensure optimal results, the following support (secondary) reagents are recommended: 1) Western Blotting: use goat anti-rabbit IgG-HRP: sc-2004 (dilution range: 1:2000-1:100,000) or Cruz Marker™ compatible goat anti-rabbit IgG-HRP: sc-2030 (dilution range: 1:2000-1:5000), Cruz Marker™ Molecular Weight Standards: sc-2035, TBS Blotto A Blocking Reagent: sc-2333 and Western Blotting Luminol Reagent: sc-2048. 2) Immunoprecipitation: use Protein A/G PLUS-Agarose: sc-2003 (0.5 ml agarose/2.0 ml). 3) Immunofluorescence: use goat anti-rabbit IgG-FITC: sc-2012 (dilution range: 1:100-1:400) or goat anti-rabbit IgG-TR: sc-2780 (dilution range: 1:100-1:400) with UltraCruz™ Mounting Medium: sc-24941.

STORAGE

Store at 4° C, **DO NOT FREEZE**. Stable for one year from the date of shipment. Non-hazardous. No MSDS required.

RESEARCH USE

For research use only, not for use in diagnostic procedures.

PROTOCOLS

See our web site at www.scbt.com or our catalog for detailed protocols and support products.