

PKC ζ siRNA (h): sc-29451

BACKGROUND

Members of the protein kinase C (PKC) family play a key regulatory role in a variety of cellular functions, including cell growth and differentiation, gene expression, hormone secretion and membrane function. PKCs were originally identified as serine/threonine protein kinases whose activity was dependent on calcium and phospholipids. Diacylglycerols (DAG) and tumor promoting phorbol esters bind to and activate PKC. PKCs can be subdivided into at least two major classes, including conventional (c) PKC isoforms (α , β I, β II and γ) and novel (n) PKC isoforms (δ , ϵ , ζ , η , θ , λ / ι , μ and ν). Patterns of expression for each PKC isoform differ among tissues and PKC family members exhibit clear differences in their cofactor dependencies. For instance, the kinase activities of PKC δ and ϵ are independent of Ca^{2+} . On the other hand, most of the other PKC members possess phorbol ester-binding activities and kinase activities.

CHROMOSOMAL LOCATION

Genetic locus: PRKCZ (human) mapping to 1p36.33.

PRODUCT

PKC ζ siRNA (h) is a target-specific 19-25 nt siRNA designed to knock down gene expression. Each vial contains 3.3 nmol of lyophilized siRNA, sufficient for a 10 μ M solution once resuspended using protocol below. Suitable for 50-100 transfections. Also see PKC ζ shRNA Plasmid (h): sc-29451-SH and PKC ζ shRNA (h) Lentiviral Particles: sc-29451-V as alternate gene silencing products.

STORAGE AND RESUSPENSION

Store lyophilized siRNA duplex at -20° C with desiccant. Stable for at least one year from the date of shipment. Once resuspended, store at -20° C, avoid contact with RNAses and repeated freeze thaw cycles.

Resuspend lyophilized siRNA duplex in 330 μ l of the RNase-free water provided. Resuspension of the siRNA duplex in 330 μ l of RNase-free water makes a 10 μ M solution in a 10 μ M Tris-HCl, pH 8.0, 20 mM NaCl, 1 mM EDTA buffered solution.

APPLICATIONS

PKC ζ siRNA (h) is recommended for the inhibition of PKC ζ expression in human cells.

SUPPORT REAGENTS

For optimal siRNA transfection efficiency, Santa Cruz Biotechnology's siRNA Transfection Reagent: sc-29528 (0.3 ml), siRNA Transfection Medium: sc-36868 (20 ml) and siRNA Dilution Buffer: sc-29527 (1.5 ml) are recommended. Control siRNAs or Fluorescein Conjugated Control siRNAs are available as 10 μ M in 66 μ l. Each contain a scrambled sequence that will not lead to the specific degradation of any known cellular mRNA. Fluorescein Conjugated Control siRNAs include: sc-36869, sc-44239, sc-44240 and sc-44241. Control siRNAs include: sc-37007, sc-44230, sc-44231, sc-44232, sc-44233, sc-44234, sc-44235, sc-44236, sc-44237 and sc-44238.

GENE EXPRESSION MONITORING

PKC ζ (H-1): sc-17781 is recommended as a control antibody for monitoring of PKC ζ gene expression knockdown by Western Blotting (starting dilution 1:200, dilution range 1:100-1:1000) or immunofluorescence (starting dilution 1:50, dilution range 1:50-1:500).

RT-PCR REAGENTS

Semi-quantitative RT-PCR may be performed to monitor PKC ζ gene expression knockdown using RT-PCR Primer: PKC ζ (h)-PR: sc-29451-PR (20 μ l, 449 bp). Annealing temperature for the primers should be 55-60 $^{\circ}$ C and the extension temperature should be 68-72 $^{\circ}$ C.

SELECT PRODUCT CITATIONS

- Xin, M., et al. 2007. Protein kinase C ζ abrogates the proapoptotic function of Bax through phosphorylation. *J. Biol. Chem.* 282: 21268-21277.
- Song, P., et al. 2008. Protein kinase C ζ -dependent LKB1 serine 428 phosphorylation increases LKB1 nucleus export and apoptosis in endothelial cells. *J. Biol. Chem.* 283: 12446-12455.
- Choi, H.C., et al. 2008. Reactive nitrogen species is required for the activation of the AMP-activated protein kinase by statin *in vivo*. *J. Biol. Chem.* 283: 20186-20197.
- Minshall, R.D., et al. 2010. Role of protein kinase C ζ in thrombin-induced RhoA activation and inter-endothelial gap formation of human dermal microvessel endothelial cell monolayers. *Microvasc. Res.* 80: 240-249.
- Géczy, T., et al. 2012. Protein kinase C isoforms have differential roles in the regulation of human sebocyte biology. *J. Invest. Dermatol.* 132: 1988-1997.
- Romano, S., et al. 2013. Synergy between enzastaurin doxorubicin in inducing melanoma apoptosis. *Pigment Cell Melanoma Res.* 26: 900-911.
- Lin, C.C., et al. 2013. Lysophosphatidic acid induces reactive oxygen species generation by activating protein kinase C in PC-3 human prostate cancer cells. *Biochem. Biophys. Res. Commun.* 440: 564-569.
- Wu, J.F., et al. 2014. Growth differentiation factor-15 induces expression of ATP-binding cassette transporter A1 through PI3-K/PKC ζ /SP1 pathway in THP-1 macrophages. *Biochem. Biophys. Res. Commun.* 444: 325-331.
- Shen-Tu, G., et al. 2014. Protein kinase C mediates enterohemorrhagic *Escherichia coli* O157:H7-induced attaching and effacing lesions. *Infect. Immun.* 82: 1648-1656.
- Miklos, W., et al. 2015. Triapine-mediated ABCB1 induction via PKC induces widespread therapy unresponsiveness but is not underlying acquired triapine resistance. *Cancer Lett.* 361: 112-120.
- Kuo, W.T., et al. 2015. LPS receptor subunits have antagonistic roles in epithelial apoptosis and colonic carcinogenesis. *Cell Death Differ.* 22: 1590-1604.

RESEARCH USE

For research use only, not for use in diagnostic procedures.