

SCD (N-19): sc-30436

BACKGROUND

Stearoyl-CoA desaturase (SCD) is a microsomal enzyme required for the synthesis of oleate and palmitoleate, which are the major monounsaturated fatty acids of membrane phospholipids, triglycerides and cholesterol esters. SCD plays a major role in the triacylglycerol and phospholipid secretion process and in mechanisms of cellular cholesterol homeostasis. It is subject to rapid turnover in the cell and, as such, represents a model for studying selective degradation of short-lived proteins of the ER. SCD is also an important regulator of membrane fluidity. An increase in expression levels of SCD is observed in cells which are induced to differentiate into adipocytes and in certain tumor cell lines. Due to gene duplication events, the number of genes in the SCD family differs between species. Their expression patterns are affected by the level of unsaturated fatty acids in the diet of the animal.

REFERENCES

1. Ntambi, J.M., Buhrow, S.A., Kaestner, K.H., Christy, R.J., Sibley, E., Kelly, T.J., Jr. and Lane, M.D. 1988. Differentiation-induced gene expression in 3T3-L1 preadipocytes. Characterization of a differentially expressed gene encoding stearoyl-CoA desaturase. *J. Biol. Chem.* 263: 17291-17300.
2. Kaestner, K.H., Ntambi, J.M., Kelly, T.J., Jr. and Lane, M.D. 1989. Differentiation-induced gene expression in 3T3-L1 preadipocytes. A second differentially expressed gene encoding stearoyl-CoA desaturase. *J. Biol. Chem.* 264: 14755-14761.
3. Li, J., Ding, S.F., Habib, N.A., Fermor, B.F., Wood, C.B. and Gilmour, R.S. 1994. Partial characterization of a cDNA for human stearoyl-CoA desaturase and changes in its mRNA expression in some normal and malignant tissues. *Int. J. Cancer* 57: 348-352.
4. Diot, C., Lefevre, P., Herve, C., Belloir, B., Narce, M., Damon, M., Poisson, J.P., Mallard, J. and Douaire, M. 2000. Stearoyl-CoA desaturase 1 coding sequences and antisense RNA affect lipid secretion in transfected chicken LMH hepatoma cells. *Arch. Biochem. Biophys.* 380: 243-250.
5. Kim, Y.C., Gomez, F.E., Fox, B.G. and Ntambi, J.M. 2000. Differential regulation of the stearoyl-CoA desaturase genes by thiazolidinediones in 3T3-L1 adipocytes. *J. Lipid Res.* 41: 1310-1316.
6. Miyazaki, M., Kim, Y.C., Gray-Keller, M.P., Attie, A.D. and Ntambi, J.M. 2000. The biosynthesis of hepatic cholesterol esters and triglycerides is impaired in mice with a disruption of the gene for stearoyl-CoA desaturase 1. *J. Biol. Chem.* 275: 30132-30138.

CHROMOSOMAL LOCATION

Genetic locus: SCD (human) mapping to 10q24.31; Scd1/Scd2/Scd3/Scd4 (mouse) mapping to 19 C3.

SOURCE

SCD (N-19) is an affinity purified goat polyclonal antibody raised against a peptide mapping near the C-terminus of SCD1 of mouse origin.

STORAGE

Store at 4° C, **DO NOT FREEZE**. Stable for one year from the date of shipment. Non-hazardous. No MSDS required.

PRODUCT

Each vial contains 200 µg IgG in 1.0 ml of PBS with < 0.1% sodium azide and 0.1% gelatin.

Blocking peptide available for competition studies, sc-30435 P, (100 µg peptide in 0.5 ml PBS containing < 0.1% sodium azide and 0.2% BSA).

APPLICATIONS

SCD (N-19) is recommended for detection of SCD of human origin, and SCD1, SCD2, SCD3 and SCD4 of mouse and rat origin by Western Blotting (starting dilution 1:200, dilution range 1:100-1:1000), immunofluorescence (starting dilution 1:50, dilution range 1:50-1:500) and solid phase ELISA (starting dilution 1:30, dilution range 1:30-1:3000); may cross-react with SCD5.

SCD (N-19) is also recommended for detection of SCD in additional species, including equine, canine, bovine and porcine.

Suitable for use as control antibody for SCD siRNA (h): sc-36464, SCD1/2/3/4 siRNA (m): sc-63288, SCD shRNA Plasmid (h): sc-36464-SH, SCD1/2/3/4 shRNA Plasmid (m): sc-63288-SH, SCD shRNA (h) Lentiviral Particles: sc-36464-V and SCD1/2/3/4 shRNA (m) Lentiviral Particles: sc-63288-V.

Molecular Weight of SCD: 40 kDa.

RECOMMENDED SECONDARY REAGENTS

To ensure optimal results, the following support (secondary) reagents are recommended: 1) Western Blotting: use donkey anti-goat IgG-HRP: sc-2020 (dilution range: 1:2000-1:100,000) or Cruz Marker™ compatible donkey anti-goat IgG-HRP: sc-2033 (dilution range: 1:2000-1:5000), Cruz Marker™ Molecular Weight Standards: sc-2035, TBS Blotto A Blocking Reagent: sc-2333 and Western Blotting Luminol Reagent: sc-2048. 2) Immunofluorescence: use donkey anti-goat IgG-FITC: sc-2024 (dilution range: 1:100-1:400) or donkey anti-goat IgG-TR: sc-2783 (dilution range: 1:100-1:400) with UltraCruz™ Mounting Medium: sc-24941.

RESEARCH USE

For research use only, not for use in diagnostic procedures.

PROTOCOLS

See our web site at www.scbt.com or our catalog for detailed protocols and support products.



Try **SCD (A00093.01): sc-81776**, our highly recommended monoclonal alternative to SCD (N-19).