

HCAM (T-20): sc-31043

BACKGROUND

Cell adhesion molecules (CAMs) are a family of closely related, cell surface glycoproteins that are involved in cell-cell interactions and are thought to play an important role in embryogenesis and development. HCAM, also known as CD44, LHR, MDU2, MDU3, MIC4, Pgp1, HCELL, MUTCH-I or ECMR-III, is a 742 amino acid single-pass type I membrane protein that is involved in hematopoiesis, lymphocyte activation and tumor metastasis. Functioning as a receptor for hyaluronic acid (HA) and interacting with ligands such as osteopontin (OPN), HCAM mediates both cell-cell and cell-matrix interactions, thereby playing an essential role in cell adhesion and cell migration. HCAM contains one Link domain and, due to alternative splicing events, is expressed as multiple isoforms, some of which are designated CD44R, CDw44, CD44S, CD44H (hematopoietic) and CD44E (epithelial). While most of the HCAM splice variants are expressed in tissues throughout the body, one specific isoform, namely CD44H, is expressed at high levels in cancer tissue, suggesting an important role for the CD44H splice variant in tumor progression.

REFERENCES

1. McVoy, L.A. and Kew, R.R. 2005. CD44 and annexin A2 mediate the C5a chemotactic cofactor function of the vitamin D binding protein. *J. Immunol.* 175: 4754-4760.
2. Hanley, W.D., Napier, S.L., Burdick, M.M., Schnaar, R.L., Sackstein, R. and Konstantopoulos, K. 2006. Variant isoforms of CD44 are P- and L-selectin ligands on colon carcinoma cells. *FASEB J.* 20: 337-339.
4. Sugahara, K.N., Hirata, T., Hayasaka, H., Stern, R., Murai, T. and Miyasaka, M. 2006. Tumor cells enhance their own CD44 cleavage and motility by generating hyaluronan fragments. *J. Biol. Chem.* 281: 5861-5868.
3. Zhuo, L., Kanamori, A., Kannagi, R., Itano, N., Wu, J., Hamaguchi, M., Ishiguro, N. and Kimata, K. 2006. SHAP potentiates the CD44-mediated leukocyte adhesion to the hyaluronan substratum. *J. Biol. Chem.* 281: 20303-20314.
5. Mielgo, A., Brondani, V., Landmann, L., Glaser-Ruhm, A., Erb, P., Stupack, D. and Günther, U. 2007. The CD44 standard/ezrin complex regulates Fas-mediated apoptosis in Jurkat cells. *Apoptosis* 12: 2051-2061.
6. Singleton, P.A., Salgia, R., Moreno-Vinasco, L., Moitra, J., Sammani, S., Mirzapooiazova, T. and Garcia, J.G. 2007. CD44 regulates hepatocyte growth factor-mediated vascular integrity. Role of c-Met, Tiam1/Rac1, dynamin 2, and cortactin. *J. Biol. Chem.* 282: 30643-30657.
7. Desai, B., Rogers, M.J. and Chellaiah, M.A. 2007. Mechanisms of osteopontin and CD44 as metastatic principles in prostate cancer cells. *Mol. Cancer* 6: 18.

CHROMOSOMAL LOCATION

Genetic locus: CD44 (human) mapping to 11p13; Cd44 (mouse) mapping to 2 E2.

SOURCE

HCAM (T-20) is an affinity purified goat polyclonal antibody raised against a peptide mapping within a C-terminal cytoplasmic domain of HCAM of human origin.

PRODUCT

Each vial contains 200 µg IgG in 1.0 ml of PBS with < 0.1% sodium azide and 0.1% gelatin.

Blocking peptide available for competition studies, sc-31043 P, (100 µg peptide in 0.5 ml PBS containing < 0.1% sodium azide and 0.2% BSA).

APPLICATIONS

HCAM (T-20) is recommended for detection of HCAM of mouse, rat and human origin by Western Blotting (starting dilution 1:200, dilution range 1:100-1:1000), immunofluorescence (starting dilution 1:50, dilution range 1:50-1:500) and solid phase ELISA (starting dilution 1:30, dilution range 1:30-1:3000).

HCAM (T-20) is also recommended for detection of HCAM in additional species, including equine, canine, bovine and porcine.

Suitable for use as control antibody for HCAM siRNA (h): sc-29342, HCAM siRNA (m): sc-35534, HCAM shRNA Plasmid (h): sc-29342-SH, HCAM shRNA Plasmid (m): sc-35534-SH, HCAM shRNA (h) Lentiviral Particles: sc-29342-V and HCAM shRNA (m) Lentiviral Particles: sc-35534-V.

Molecular Weight of HCAM: 90-95 kDa.

Positive Controls: HISM cell lysate: sc-2229, HL-60 whole cell lysate: sc-2209 or HeLa whole cell lysate: sc-2200.

RECOMMENDED SECONDARY REAGENTS

To ensure optimal results, the following support (secondary) reagents are recommended: 1) Western Blotting: use donkey anti-goat IgG-HRP: sc-2020 (dilution range: 1:2000-1:100,000) or Cruz Marker™ compatible donkey anti-goat IgG-HRP: sc-2033 (dilution range: 1:2000-1:5000), Cruz Marker™ Molecular Weight Standards: sc-2035, TBS Blotto A Blocking Reagent: sc-2333 and Western Blotting Luminol Reagent: sc-2048. 2) Immunofluorescence: use donkey anti-goat IgG-FITC: sc-2024 (dilution range: 1:100-1:400) or donkey anti-goat IgG-TR: sc-2783 (dilution range: 1:100-1:400) with UltraCruz™ Mounting Medium: sc-24941.

STORAGE

Store at 4° C, ****DO NOT FREEZE****. Stable for one year from the date of shipment. Non-hazardous. No MSDS required.

RESEARCH USE

For research use only, not for use in diagnostic procedures.

PROTOCOLS

See our web site at www.scbt.com or our catalog for detailed protocols and support products.