

NHE-3 (D-18): sc-31384

BACKGROUND

Na⁺/H⁺ exchangers-1-8 (also designated Na⁺/H⁺ antiporters or NHE-1-8) are integral membrane proteins that are expressed in most mammalian tissues, where they regulate intracellular pH and cell volume. NHEs mediate the transport of hydrogen (H⁺) ions out of cells in exchange for extracellular sodium (Na⁺) ions. While NHE-1 is ubiquitously expressed, the NHE isoforms 2-8 have distinct tissue- and cell type-dependent expression and inhibitory characteristics. NHE-3 localizes to the apical membrane of renal proximal tubules where it is responsible for most of the sodium transport and fluid reabsorption. NHE-3 translocates to internal pools where it mediates natriuresis when blood pressure increases abruptly. NHE-3 is also expressed in the stomach and functions to protect the mucosa by secreting protons that diffuse into the mucous cells.

REFERENCES

- Orlowski, J., Kandasamy, R.A. and Shull, G.E. 1992. Molecular cloning of putative members of the Na/H exchanger gene family. cDNA cloning, deduced amino acid sequence, and mRNA tissue expression of the rat Na/H exchanger NHE-1 and two structurally related proteins. *J. Biol. Chem.* 267: 9331-9339.
- Harris, S.P., Strong, T.V., Wys, N., Richards, N.W., Pouyssegur, J., Ernst, S.A. and Dawson, D.C. 1997. Epithelial localization of a reptilian Na⁺/H⁺ exchanger homologous to NHE-1. *Am. J. Physiol.* 272: 1594-1606.
- Kulaksiz, H., Bektas, H. and Cetin, Y. 2001. Expression and cell-specific and membrane-specific localization of NHE-3 in the human and guinea pig upper gastrointestinal tract. *Cell Tissue Res.* 303: 337-343.
- LaPointe, M.S., Sodhi, C., Sahai, A. and Battle, D. 2002. Na⁺/H⁺ exchange activity and NHE-3 expression in renal tubules from the spontaneously hypertensive rat. *Kidney Int.* 62: 157-165.
- Sangan, P., Rajendran, V.M., Geibel, J.P. and Binder, H.J. 2002. Cloning and expression of a chloride-dependent Na⁺/H⁺ exchanger. *J. Biol. Chem.* 277: 9668-9675.
- Goyal, S., Vanden Heuvel, G. and Aronson, P.S. 2003. Renal expression of novel Na⁺/H⁺ exchanger isoform NHE-8. *Am. J. Physiol. Renal Physiol.* 284: 467-473.
- Goyal, S., Mentone, S. and Aronson, P.S. 2005. Immunolocalization of NHE-8 in rat kidney. *Am. J. Physiol. Renal Physiol.* 288: 530-538.
- Xu, H., Chen, R. and Ghishan, F.K. 2005. Subcloning, localization, and expression of the rat intestinal sodium-hydrogen exchanger isoform 8. *Am. J. Physiol. Gastrointest. Liver Physiol.* 289: 36-41.

CHROMOSOMAL LOCATION

Genetic locus: SLC9A3 (human) mapping to 5p15.33; Slc9a3 (mouse) mapping to 13 C1.

SOURCE

NHE-3 (D-18) is an affinity purified goat polyclonal antibody raised against a peptide mapping within a C-terminal cytoplasmic domain of NHE-3 of human origin.

PRODUCT

Each vial contains 200 µg IgG in 1.0 ml of PBS with < 0.1% sodium azide and 0.1% gelatin.

Blocking peptide available for competition studies, sc-31384 P, (100 µg peptide in 0.5 ml PBS containing < 0.1% sodium azide and 0.2% BSA).

APPLICATIONS

NHE-3 (D-18) is recommended for detection of NHE-3 of human and, to a lesser extent, mouse and rat origin by Western Blotting (starting dilution 1:200, dilution range 1:100-1:1000), immunofluorescence (starting dilution 1:50, dilution range 1:50-1:500) and solid phase ELISA (starting dilution 1:30, dilution range 1:30-1:3000).

Suitable for use as control antibody for NHE-3 siRNA (h): sc-36059, NHE-3 siRNA (m): sc-36060, NHE-3 shRNA Plasmid (h): sc-36059-SH, NHE-3 shRNA Plasmid (m): sc-36060-SH, NHE-3 shRNA (h) Lentiviral Particles: sc-36059-V and NHE-3 shRNA (m) Lentiviral Particles: sc-36060-V.

Molecular Weight of NHE-3 (D-18): 93 kDa.

Positive Controls: MDCK cell lysate: sc-2252 or rat kidney extract: sc-2394.

RECOMMENDED SECONDARY REAGENTS

To ensure optimal results, the following support (secondary) reagents are recommended: 1) Western Blotting: use donkey anti-goat IgG-HRP: sc-2020 (dilution range: 1:2000-1:100,000) or Cruz Marker™ compatible donkey anti-goat IgG-HRP: sc-2033 (dilution range: 1:2000-1:5000), Cruz Marker™ Molecular Weight Standards: sc-2035, TBS Blotto A Blocking Reagent: sc-2333 and Western Blotting Luminol Reagent: sc-2048. 2) Immunofluorescence: use donkey anti-goat IgG-FITC: sc-2024 (dilution range: 1:100-1:400) or donkey anti-goat IgG-TR: sc-2783 (dilution range: 1:100-1:400) with UltraCruz™ Mounting Medium: sc-24941.

STORAGE

Store at 4° C, **DO NOT FREEZE**. Stable for one year from the date of shipment. Non-hazardous. No MSDS required.

RESEARCH USE

For research use only, not for use in diagnostic procedures.

PROTOCOLS

See our web site at www.scbt.com or our catalog for detailed protocols and support products.



MONOS
Satisfation
Guaranteed

Try **NHE-3 (53): sc-136368**, our highly recommended monoclonal alternative to NHE-3 (D-18).