


ABCD1 (A-15): sc-31875



BACKGROUND

The peroxisomal membrane contains several ATP-binding cassette (ABC) transporters, ABCD1-4 that are known to be present in the human peroxisome membrane. All four proteins are ABC half-transporters, which dimerize to form an active transporter. A mutation in the ABCD1 causes X-linked adrenoleukodystrophy (X-ALD), a peroxisomal disorder which affects lipid storage. ABCD2 in mouse, is expressed at high levels in the brain and adrenal organs, which are adversely affected in X-ALD. The peroxisomal membrane comprises 2 quantitatively major proteins, PMP22 and ABCD3. ABCD3 is associated with irregularly shaped vesicles which may be defective peroxisomes or peroxisome precursors. ABCD4 localizes to peroxisomes. The genes which encode ABCD1-4 map to human chromosome Xq28, 12q11-q12, 1p22-p21 and 14q24.3, respectively. ABCB7 is a half-transporter involved in the transport of heme from the mitochondria to the cytosol and maps to human chromosome Xq13.1-q13.3.

REFERENCES

1. Gartner, J., Kearns, W., Pearson, P. and Valle, D. 1992. Characterization and localization of the human 70 kDa peroxisomal membrane protein (PMP70) gene. *Am. J. Hum. Gen.* 51: 168.
2. Lombard-Platet, G., Savary, S., Sarde, C.O., Mandel, J.L. and Chimini, G. 1996. A close relative of the adrenoleukodystrophy (ALD) gene codes for a peroxisomal protein with a specific expression pattern. *Proc. Natl. Acad. Sci. USA* 93: 1265-1269.
3. Shani, N., Jimenez-Sanchez, G., Steel, G., Dean, M. and Valle, D. 1997. Identification of a fourth half ABC transporter in the human peroxisomal membrane. *Hum. Mol. Genet.* 6: 1925-1931.
4. Moser, H.W. 1997. Adrenoleukodystrophy: phenotype, genetics, pathogenesis and therapy. *Brain* 120: 1485-1508.
5. Savary, S., Troffer-Charlier, N., Gyapay, G., Mattei, M.G. and Chimini, G. 1997. Chromosomal localization of the adrenoleukodystrophy-related gene in man and mice. *Eur. J. Hum. Genet.* 5: 99-101.
6. Holzinger, A., Roscher, A.A., Landgraf, P., Lichtner, P. and Kammerer, S. 1998. Genomic organization and chromosomal localization of the human peroxisomal membrane protein-1-like protein (PXMP1-L) gene encoding a peroxisomal ABC transporter. *FEBS Lett.* 426: 238-242.
7. Shimada, Y., Okuno, S., Kawai, A., Shinomiya, H., Saito, A., Suzuki, M., Omori, Y., Nishino, N., Kanemoto, N., Fujiwara, T., Horie, M. and Takahashi, E. 1998. Cloning and chromosomal mapping of a novel ABC transporter gene (hABC7), a candidate for X-linked sideroblastic anemia with spinocerebellar ataxia. *J. Hum. Genet.* 43: 115-122.
8. LocusLink Report (LocusID:300100). <http://www.ncbi.nlm.nih.gov/LocusLink/>

CHROMOSOMAL LOCATION

Genetic locus: ABCD1 (human) mapping to Xq28; Abcd1 (mouse) mapping to X A7.3.

SOURCE

ABCD1 (A-15) is an affinity purified goat polyclonal antibody raised against a peptide mapping within an internal region of ABCD1 of human origin.

PRODUCT

Each vial contains 200 µg IgG in 1.0 ml of PBS with < 0.1% sodium azide and 0.1% gelatin.

Blocking peptide available for competition studies, sc-31875 P, (100 µg peptide in 0.5 ml PBS containing < 0.1% sodium azide and 0.2% BSA).

APPLICATIONS

ABCD1 (A-15) is recommended for detection of ABCD1 of mouse, rat and human origin by Western Blotting (starting dilution 1:200, dilution range 1:100-1:1000), immunofluorescence (starting dilution 1:50, dilution range 1:50-1:500) and solid phase ELISA (starting dilution 1:30, dilution range 1:30-1:3000).

ABCD1 (A-15) is also recommended for detection of ABCD1 in additional species, including equine, canine, bovine and porcine.

Suitable for use as control antibody for ABCD1 siRNA (h): sc-41143, ABCD1 siRNA (m): sc-41144, ABCD1 shRNA Plasmid (h): sc-41143-SH, ABCD1 shRNA Plasmid (m): sc-41144-SH, ABCD1 shRNA (h) Lentiviral Particles: sc-41143-V and ABCD1 shRNA (m) Lentiviral Particles: sc-41144-V.

Molecular Weight of ABCD1: 75 kDa.

Positive Controls: HeLa whole cell lysate: sc-2200 or mouse brain extract: sc-2253.

RECOMMENDED SECONDARY REAGENTS

To ensure optimal results, the following support (secondary) reagents are recommended: 1) Western Blotting: use donkey anti-goat IgG-HRP: sc-2020 (dilution range: 1:2000-1:100,000) or Cruz Marker™ compatible donkey anti-goat IgG-HRP: sc-2033 (dilution range: 1:2000-1:5000), Cruz Marker™ Molecular Weight Standards: sc-2035, TBS Blotto A Blocking Reagent: sc-2333 and Western Blotting Luminol Reagent: sc-2048. 2) Immunofluorescence: use donkey anti-goat IgG-FITC: sc-2024 (dilution range: 1:100-1:400) or donkey anti-goat IgG-TR: sc-2783 (dilution range: 1:100-1:400) with UltraCruz™ Mounting Medium: sc-24941.

STORAGE

Store at 4° C, **DO NOT FREEZE**. Stable for one year from the date of shipment. Non-hazardous. No MSDS required.

RESEARCH USE

For research use only, not for use in diagnostic procedures.

PROTOCOLS

See our web site at www.scbt.com or our catalog for detailed protocols and support products.