

IL-6 (Y-15): sc-323976

BACKGROUND

Interleukin-6, or IL-6, is a multifunctional protein, 212 amino acids in length, that plays critical roles in host defense, immune response and hematopoiesis. IL-6 is constitutively expressed by epidermal Langerhans cells and its expression is induced in stimulated keratinocytes. IL-6, IL-1 β and TNF α act as endogenous pyrogens, regulating the fever response to bacterial invasion. The IL-6 receptor is a trimeric complex composed of an IL-6-specific α chain and a homodimer of the gp130 glycoprotein common to the IL-6, IL-11, CNTF, OSM and LIF receptors. Stimulation with IL-6 leads to gp130 homodimerization and the activation of associated kinases JAK1 and JAK2. Once activated, JAK1 and JAK2 phosphorylate Stat3, causing its nuclear translocation and transcription of Stat3-responsive genes. IL-6 has also been shown to activate the Ras/MAP kinase pathway, which regulates NFIL6 transcription.

REFERENCES

- Hirano, T., et al. 1986. Complementary DNA for a novel human interleukin (BSF-2) that induces B lymphocytes to produce immunoglobulin. *Nature* 324: 73-76.
- Nakajima, T., et al. 1993. Phosphorylation at threonine 235 by a Ras-dependent mitogen-activated protein kinase cascade is essential for transcription factor NFIL6. *Proc. Natl. Acad. Sci. USA* 90: 2207-2211.
- Zhong, Z., et al. 1994. Stat3: A Stat family member activated by tyrosine phosphorylation in response to epidermal growth factor and interleukin-6. *Science* 264: 95-98.
- Wijdenes, J., et al. 1995. Interleukin-6 signal transducer gp130 has specific binding sites for different cytokines as determined by antagonistic and agonistic anti-gp130 monoclonal antibodies. *Eur. J. Immunol.* 25: 3474-3481.
- Wang, Y., et al. 1995. Interleukin-6 and ciliary neurotrophic factor trigger Janus kinase activation and early gene response in rat hepatocytes. *Gene* 162: 285-289.
- Holliday, M.R., et al. 1996. Stimulation by oxazolone of increased IL-6, but not IL-10, in the skin of mice. *Toxicology* 106: 237-242.
- Chai, Z., et al. 1996. Interleukin (IL)-6 gene expression in the central nervous system is necessary for fever response to lipopolysaccharide or IL-1 β : a study on IL-6-deficient mice. *J. Exp. Med.* 183: 311-316.

CHROMOSOMAL LOCATION

Genetic locus: IL6 (human) mapping to 7p15.3.

SOURCE

IL-6 (Y-15) is an affinity purified goat polyclonal antibody raised against a peptide mapping within an internal region of IL-6 of human origin.

PRODUCT

Each vial contains 200 μ g IgG in 1.0 ml of PBS with < 0.1% sodium azide and 0.1% gelatin.

Blocking peptide available for competition studies, sc-323976 P, (100 μ g peptide in 0.5 ml PBS containing < 0.1% sodium azide and 0.2% BSA).

APPLICATIONS

IL-6 (Y-15) is recommended for detection of IL-6 of human origin by Western Blotting (starting dilution 1:200, dilution range 1:100-1:1000), immunofluorescence (starting dilution 1:50, dilution range 1:50-1:500) and solid phase ELISA (starting dilution 1:30, dilution range 1:30-1:3000); non cross-reactive with other Interleukin proteins.

Suitable for use as control antibody for IL-6 siRNA (h): sc-39627, IL-6 shRNA Plasmid (h): sc-39627-SH and IL-6 shRNA (h) Lentiviral Particles: sc-39627-V.

Molecular Weight of IL-6: 21 kDa.

RECOMMENDED SECONDARY REAGENTS

To ensure optimal results, the following support (secondary) reagents are recommended: 1) Western Blotting: use donkey anti-goat IgG-HRP: sc-2020 (dilution range: 1:2000-1:100,000) or Cruz Marker[™] compatible donkey anti-goat IgG-HRP: sc-2033 (dilution range: 1:2000-1:5000), Cruz Marker[™] Molecular Weight Standards: sc-2035, TBS Blotto A Blocking Reagent: sc-2333 and Western Blotting Luminol Reagent: sc-2048. 2) Immunofluorescence: use donkey anti-goat IgG-FITC: sc-2024 (dilution range: 1:100-1:400) or donkey anti-goat IgG-TR: sc-2783 (dilution range: 1:100-1:400) with UltraCruz[™] Mounting Medium: sc-24941.

STORAGE

Store at 4° C, ****DO NOT FREEZE****. Stable for one year from the date of shipment. Non-hazardous. No MSDS required.

RESEARCH USE

For research use only, not for use in diagnostic procedures.

PROTOCOLS

See our web site at www.scbt.com or our catalog for detailed protocols and support products.