

nephrin (N-16): sc-32532

BACKGROUND

Nephrin is a member of the immunoglobulin family of cell adhesion molecules that localizes to opposing sites of the secondary foot processes formed by podocytes, a specialized epithelial cell that ensures size- and charge-selective ultrafiltration. The human nephrin gene maps to chromosome 19q13.12 and encodes a 1,241 amino acid protein that contains a transmembrane domain, 8 Ig-like modules, and one Fibronectin III-like module. Nephrin is expressed in embryonic and adult kidneys and localizes to glomerular podocytes and the glomerular slit diaphragm. Nephrin stimulates mitogen-activated protein kinases and is enhanced by Podocin, which binds to the cytoplasmic tail of nephrin. A293 cells treated with Phorbol-12-myristate-13-acetate can upregulate nephrin, suggesting that protein kinase C is part of an intracellular signalling system, which regulates nephrin.

CHROMOSOMAL LOCATION

Genetic locus: NPHS1 (human) mapping to 19q13.12; Nphs1 (mouse) mapping to 7 B1.

SOURCE

nephrin (N-16) is an affinity purified goat polyclonal antibody raised against a peptide mapping near the N-terminus of nephrin of human origin.

PRODUCT

Each vial contains 200 µg IgG in 1.0 ml of PBS with < 0.1% sodium azide and 0.1% gelatin.

Blocking peptide available for competition studies, sc-32532 P, (100 µg peptide in 0.5 ml PBS containing < 0.1% sodium azide and 0.2% BSA).

STORAGE

Store at 4° C, ****DO NOT FREEZE****. Stable for one year from the date of shipment. Non-hazardous. No MSDS required.

APPLICATIONS

nephrin (N-16) is recommended for detection of nephrin and nephrin precursor of mouse, rat and human origin by Western Blotting (starting dilution 1:200, dilution range 1:100-1:1000), immunoprecipitation [1-2 µg per 100-500 µg of total protein (1 ml of cell lysate)], immunofluorescence (starting dilution 1:50, dilution range 1:50-1:500) and solid phase ELISA (starting dilution 1:30, dilution range 1:30-1:3000).

nephrin (N-16) is also recommended for detection of nephrin and nephrin precursor in additional species, including canine, bovine and porcine.

Suitable for use as control antibody for nephrin siRNA (h): sc-36030, nephrin siRNA (m): sc-36031, nephrin shRNA Plasmid (h): sc-36030-SH, nephrin shRNA Plasmid (m): sc-36031-SH, nephrin shRNA (h) Lentiviral Particles: sc-36030-V and nephrin shRNA (m) Lentiviral Particles: sc-36031-V.


Molecular Weight of nephrin: 185 kDa.

Positive Controls: KNRK whole cell lysate: sc-2214 or human kidney extract: sc-363764.


RECOMMENDED SECONDARY REAGENTS

To ensure optimal results, the following support (secondary) reagents are recommended: 1) Western Blotting: use donkey anti-goat IgG-HRP: sc-2020 (dilution range: 1:2000-1:100,000) or Cruz Marker™ compatible donkey anti-goat IgG-HRP: sc-2033 (dilution range: 1:2000-1:5000), Cruz Marker™ Molecular Weight Standards: sc-2035, TBS Blotto A Blocking Reagent: sc-2333 and Western Blotting Luminol Reagent: sc-2048. 2) Immunoprecipitation: use Protein A/G PLUS-Agarose: sc-2003 (0.5 ml agarose/2.0 ml). 3) Immunofluorescence: use donkey anti-goat IgG-FITC: sc-2024 (dilution range: 1:100-1:400) or donkey anti-goat IgG-TR: sc-2783 (dilution range: 1:100-1:400) with UltraCruz™ Mounting Medium: sc-24941.

DATA



nephrin (N-16): sc-32532. Western blot analysis of nephrin expression in KNRK whole cell lysate.



nephrin (N-16): sc-32532. Immunofluorescence staining of formalin-fixed HeLa cells showing cytoplasmic localization.

RESEARCH USE

For research use only, not for use in diagnostic procedures.

PROTOCOLS

See our web site at www.scbt.com or our catalog for detailed protocols and support products.



Try **nephrin (G-8): sc-376522** or **nephrin (B-12): sc-377246**, our highly recommended monoclonal alternatives to nephrin (N-16). Also, for AC, HRP, FITC, PE, Alexa Fluor® 488 and Alexa Fluor® 647 conjugates, see **nephrin (G-8): sc-376522**.