

## SR-2A (L-20): sc-32537

### BACKGROUND

Serotonin (also designated 5-hydroxytryptamine or 5-HT) is a molecule that functions as a neurotransmitter, a hormone and a mitogen, and it is predominantly expressed in the gut, platelets and central nervous system (CNS). In the CNS, serotonin modulates several processes, including anxiety, sleep, appetite, behavior and drug abuse. In platelets and gut, serotonin plays a major role in cardiovascular function and motility of the gastrointestinal tract, respectively. Serotonin mediates its effects through several of G protein-coupled receptors, designated 5-HT receptors or alternatively SR receptors. The SR-2 receptors are comprised of three subtypes, SR-2A, SR-2B and SR-2C, which activate phospholipase C and release intracellular stores of calcium in response to serotonin. SR-2A has a specific role in tracheal smooth muscle contraction, bronchoconstriction and mediating aldosterone production, and it is also thought to play a role in several psychiatric disorders, including depression and schizophrenia. SR-2B is expressed in embryonic and adult cardiovascular tissues, gut and brain, and plays an important role in the pathology of cardiac disorders. SR-2C is thought to mediate the effects of atypical antipsychotic drugs.

### REFERENCES

1. Watts, S.W., et al. 1994. Contractile serotonin-2A receptor signal transduction in guinea pig trachea: importance of protein kinase C and extracellular and intracellular calcium but not phosphoinositide hydrolysis. *J. Pharmacol. Exp. Ther.* 271: 832-844.
2. Goppelt-Strube, M., et al. 1998. Signaling pathways mediating induction of the early response genes prostaglandin G/H synthase-2 and Egr-1 by serotonin via 5-HT<sub>2A</sub> receptors. *J. Cell. Physiol.* 175: 341-347.
3. Nebigil, C.G., et al. 2000. Serotonin-2B receptor is required for heart development. *Proc. Natl. Acad. Sci. USA* 97: 9508-9513.
4. Contesse, V., et al. 2000. Role of 5-HT in the regulation of the brain-pituitary-adrenal axis: effects of 5-HT on adrenocortical cells. *Can. J. Physiol. Pharmacol.* 78: 967-983.
5. Xu, T., et al. 2000. Cellular localization of serotonin-2A (5-HT<sub>2A</sub>) receptors in the rat brain. *Brain Res. Bull.* 51: 499-505.

### CHROMOSOMAL LOCATION

Genetic locus: HTR2A (human) mapping to 13q14.2; Htr2a (mouse) mapping to 14 D3.

### SOURCE

SR-2A (L-20) is an affinity purified goat polyclonal antibody raised against a peptide mapping at the N-terminus of SR-2A of human origin.

### PRODUCT

Each vial contains 200 µg IgG in 1.0 ml of PBS with < 0.1% sodium azide and 0.1% gelatin.

Blocking peptide available for competition studies, sc-32537 P, (100 µg peptide in 0.5 ml PBS containing < 0.1% sodium azide and 0.2% BSA).

### APPLICATIONS

SR-2A (L-20) is recommended for detection of serotonin 2A receptor (5-HT<sub>2</sub>) of human and, to a lesser extent, mouse and rat origin by Western Blotting (starting dilution 1:200, dilution range 1:100-1:1000), immunoprecipitation [1-2 µg per 100-500 µg of total protein (1 ml of cell lysate)], immunofluorescence (starting dilution 1:50, dilution range 1:50-1:500) and solid phase ELISA (starting dilution 1:30, dilution range 1:30-1:3000).


SR-2A (L-20) is also recommended for detection of serotonin 2A receptor (5-HT<sub>2</sub>) in additional species, including porcine.

Suitable for use as control antibody for SR-2A siRNA (h): sc-42231, SR-2A siRNA (m): sc-42232, SR-2A shRNA Plasmid (h): sc-42231-SH, SR-2A shRNA Plasmid (m): sc-42232-SH, SR-2A shRNA (h) Lentiviral Particles: sc-42231-V and SR-2A shRNA (m) Lentiviral Particles: sc-42232-V.

Molecular Weight of SR-2A: 55 kDa.

Positive Controls: HEK293 whole cell lysate: sc-45136.

### DATA



SR-2A (L-20): sc-32537. Western blot analysis of human recombinant SR-2A fusion protein.

### STORAGE

Store at 4° C, **\*\*DO NOT FREEZE\*\***. Stable for one year from the date of shipment. Non-hazardous. No MSDS required.

### RESEARCH USE

For research use only, not for use in diagnostic procedures.

### PROTOCOLS

See our web site at [www.scbt.com](http://www.scbt.com) or our catalog for detailed protocols and support products.



Try **SR-2A (A-4): sc-166775**, our highly recommended monoclonal alternative to SR-2A (L-20).