



GOA-1 (cT-14): sc-32684

BACKGROUND

GOA-1, the $G_{\alpha 0}$ subunit of a heterotrimeric G protein complex, is involved in asymmetric cell division in *C. elegans*. Specifically, GOA-1 assists in the positioning of spindle fibers in one-celled *C. elegans* embryos. Proper spindle fiber positioning also depends upon the interaction of GOA-1 in its GDP bound form to the receptor independent G protein activators GPR-1 and GPR-2, as well as the protein RIC-8, a guanine nucleotide exchange factor (GEF). GOA-1 is localized to the cell cortex and is found in most neurons. Additionally, muscles involved with male mating and egg laying are sites of GOA-1 expression in *C. elegans*. Mutations in the gene coding for GOA-1 are associated with a variety of deficits mimicking serotonin deficiency in *C. elegans*, including male impotence, premature egg laying and hyperactive movement. EGL-10 is a regulator of G protein signaling (RGS) protein that can selectively inhibit GOA-1 function through activation of G_{α} GTPase.

REFERENCES

1. Segalat, L., et al. 1995. Modulation of serotonin-controlled behaviors by G_0 in *Caenorhabditis elegans*. *Science* 267: 1648-1651.
2. Mendel, J.E., et al. 1995. Participation of the protein G_0 in multiple aspects of behavior in *C. elegans*. *Science* 267: 1652-1655.
3. Miller, K.G., et al. 1999. $G_{\alpha\alpha}$ and diacylglycerol kinase negatively regulate the $G_{\alpha\alpha}$ pathway in *C. elegans*. *Neuron* 24: 323-333.
4. Miller, K.G., et al. 2000. A role for RIC-8 (Synembryn) and GOA-1 ($G_{\alpha\alpha}$) in regulating a subset of centrosome movements during early embryogenesis in *Caenorhabditis elegans*. *Genetics* 156: 1649-1660.
5. van Swinderen, B., et al. 2001. $G_{\alpha\alpha}$ regulates volatile anesthetic action in *Caenorhabditis elegans*. *Genetics* 158: 643-655.
6. Robatzek, M., et al. 2001. EAT-11 encodes GPB-2, a $G_{\beta 5}$ ortholog that interacts with $G_{\alpha\alpha}$ and $G_{\alpha\alpha}$ to regulate *C. elegans* behavior. *Curr. Biol.* 11: 288-293.
7. Gotta, M., et al. 2001. Distinct roles for G_{α} and $G_{\beta\gamma}$ in regulating spindle position and orientation in *Caenorhabditis elegans* embryos. *Nat. Cell Biol.* 3: 297-300.
8. Patikoglou, G.A., et al. 2002. An N-terminal region of *Caenorhabditis elegans* RGS proteins EGL-10 and EAT-16 directs inhibition of $G_{\alpha\alpha}$ versus $G_{\alpha\alpha}$ signaling. *J. Biol. Chem.* 277: 47004-47013.
9. Manning, D.R., et al. 2003. Evidence mounts for receptor-independent activation of heterotrimeric G proteins normally *in vivo*: positioning of the mitotic spindle in *C. elegans*. *Sci. STKE*. 2003: 35.

SOURCE

GOA-1 (cT-14) is an affinity purified goat polyclonal antibody raised against a peptide mapping at the N-terminus of GOA-1 of *C. elegans* origin.

STORAGE

Store at 4° C, ****DO NOT FREEZE****. Stable for one year from the date of shipment. Non-hazardous. No MSDS required.

PRODUCT

Each vial contains 200 μ g IgG in 1.0 ml of PBS with < 0.1% sodium azide and 0.1% gelatin.

Blocking peptide available for competition studies, sc-32684 P, (100 μ g peptide in 0.5 ml PBS containing < 0.1% sodium azide and 0.2% BSA).

APPLICATIONS

GOA-1 (cT-14) is recommended for detection of GOA-1 of *Caenorhabditis elegans* origin by Western Blotting (starting dilution 1:200, dilution range 1:100-1:1000), immunofluorescence (starting dilution 1:50, dilution range 1:50-1:500) and solid phase ELISA (starting dilution 1:30, dilution range 1:30-1:3000).

Molecular Weight of GOA-1: 40 kDa.

RECOMMENDED SECONDARY REAGENTS

To ensure optimal results, the following support (secondary) reagents are recommended: 1) Western Blotting: use donkey anti-goat IgG-HRP: sc-2020 (dilution range: 1:2000-1:100,000) or Cruz Marker™ compatible donkey anti-goat IgG-HRP: sc-2033 (dilution range: 1:2000-1:5000), Cruz Marker™ Molecular Weight Standards: sc-2035, TBS Blotto A Blocking Reagent: sc-2333 and Western Blotting Luminol Reagent: sc-2048. 2) Immunofluorescence: use donkey anti-goat IgG-FITC: sc-2024 (dilution range: 1:100-1:400) or donkey anti-goat IgG-TR: sc-2783 (dilution range: 1:100-1:400) with UltraCruz™ Mounting Medium: sc-24941.

RESEARCH USE

For research use only, not for use in diagnostic procedures.

PROTOCOLS

See our web site at www.scbt.com or our catalog for detailed protocols and support products.