

Autotaxin (N-13): sc-33075

BACKGROUND

Autotaxin (ATX), also designated ectonucleotide pyrophosphatase/phosphodiesterase 2 (E-NPP 2), is a membrane-bound glycoprotein that cleaves diester bonds for a broad range of substrates. Originally isolated from the human melanoma cell line (A2058), Autotaxin is predominantly expressed in brain, placenta, ovary and small intestine. Autotaxin has significant homology to the cell membrane differentiation antigen PC-1, and is a stimulator of tumor cell motility. It also functions as a catalyst by hydrolytically removing 5'-nucleotides from the 3'-hydroxy termini of 3'-hydroxy-terminated oligonucleotides.

CHROMOSOMAL LOCATION

Genetic locus: ENPP2 (human) mapping to 8q24.12; Enpp2 (mouse) mapping to 15 D1.

SOURCE

Autotaxin (N-13) is an affinity purified goat polyclonal antibody raised against a peptide mapping within an N-terminal extracellular domain of Autotaxin of human origin.

PRODUCT

Each vial contains 200 µg IgG in 1.0 ml of PBS with < 0.1% sodium azide and 0.1% gelatin.

Blocking peptide available for competition studies, sc-33075 P, (100 µg peptide in 0.5 ml PBS containing < 0.1% sodium azide and 0.2% BSA).

STORAGE

Store at 4° C, ****DO NOT FREEZE****. Stable for one year from the date of shipment. Non-hazardous. No MSDS required.

RESEARCH USE

For research use only, not for use in diagnostic procedures.

APPLICATIONS

Autotaxin (N-13) is recommended for detection of Autotaxin of mouse, rat and human origin by Western Blotting (starting dilution 1:200, dilution range 1:100-1:1000), immunoprecipitation [1-2 µg per 100-500 µg of total protein (1 ml of cell lysate)], immunofluorescence (starting dilution 1:50, dilution range 1:50-1:500) and solid phase ELISA (starting dilution 1:30, dilution range 1:30-1:3000).

Autotaxin (N-13) is also recommended for detection of Autotaxin in additional species, including canine, bovine and porcine.

Suitable for use as control antibody for Autotaxin siRNA (h): sc-44906, Autotaxin siRNA (m): sc-44907, Autotaxin shRNA Plasmid (h): sc-44906-SH, Autotaxin shRNA Plasmid (m): sc-44907-SH, Autotaxin shRNA (h) Lentiviral Particles: sc-44906-V and Autotaxin shRNA (m) Lentiviral Particles: sc-44907-V.


Molecular Weight of Autotaxin: 125 kDa.

Positive Controls: Hep G2 cell lysate: sc-2227 or MOLT-4 cell lysate: sc-2233.

RECOMMENDED SECONDARY REAGENTS

To ensure optimal results, the following support (secondary) reagents are recommended: 1) Western Blotting: use donkey anti-goat IgG-HRP: sc-2020 (dilution range: 1:2000-1:100,000) or Cruz Marker™ compatible donkey anti-goat IgG-HRP: sc-2033 (dilution range: 1:2000-1:5000), Cruz Marker™ Molecular Weight Standards: sc-2035, TBS Blotto A Blocking Reagent: sc-2333 and Western Blotting Luminol Reagent: sc-2048. 2) Immunoprecipitation: use Protein A/G PLUS-Agarose: sc-2003 (0.5 ml agarose/2.0 ml). 3) Immunofluorescence: use donkey anti-goat IgG-FITC: sc-2024 (dilution range: 1:100-1:400) or donkey anti-goat IgG-TR: sc-2783 (dilution range: 1:100-1:400) with UltraCruz™ Mounting Medium: sc-24941.

DATA



Autotaxin (N-13): sc-33075. Western blot analysis of Autotaxin expression in Hep G2 (A) and MOLT-4 (B) whole cell lysates.

PROTOCOLS

See our web site at www.scbt.com or our catalog for detailed protocols and support products.


 MONOS
Satisfaction
Guaranteed

Try **Autotaxin (E-12): sc-374222**, our highly recommended monoclonal alternative to Autotaxin (N-13).