## SANTA CRUZ BIOTECHNOLOGY, INC.

# Bicoid (dD-20): sc-33288



## BACKGROUND

*Drosophila melanogaster* is a proven and effective model for studying developmental and cellular processes common to higher eukaryotes. Approximately 13,600 genes have been elucidated from more than 120 megabases of euchromatin, and they are organized among the chromosomes 2, 3, 4, X and Y, with the Y chromosome being predominately heterochromatic. *Drosophila* genes can be categorized based on the type of protein for which they encode and are represented by six major classifications, which include intracellular signaling proteins, transmembrane proteins, RNA binding proteins, secreted factors, transcription regulators (basic helix-loop-helix, homeodomain containing, zinc finger containing and chromatin associated) or other functional proteins. The Bicoid protein is critical for patterning the anterior-posterior body plan in *Drosophila*. Bicoid, a concentration-dependent transcriptional activator in the *Drosophila* embryo, regulates the expression of the hunchback gene, which shows a step-like-function expression pattern.

#### REFERENCES

- Adams, M.D., et al. 2000. The genome sequence of *Drosophila* melanogaster. Science 287: 2185-2195.
- Fu, D., et al. 2005. Interplay between positive and negative activities that influence the role of Bicoid in transcription. Nucleic Acids Res. 33: 3985-3993.
- Macdonald, P.M., et al. 2005. Translational repression by Bicoid: competition for the cap. Cell 121: 321-322.
- Lopes, F.J., et al. 2005. A kinetic mechanism for *Drosophila* Bicoid cooperative binding. J. Theor. Biol. 235: 185-198.
- Ochoa-Espinosa, A., et al. 2005. The role of binding site cluster strength in Bicoid-dependent patterning in *Drosophila*. Proc. Natl. Acad. Sci. USA 102: 4960-4965.
- Snee, M.J., et al. 2005. Recognition of the bcd mRNA localization signal in *Drosophila* embryos and ovaries. Mol. Cell. Biol. 25: 1501-1510.
- 7. The Interactive Fly. http://sdb.bio.purdue.edu/fly/aimain/1aahome.htm. http://sdb.bio.purdue.edu/fly/aimain/6biochem.htm

### SOURCE

Bicoid (dD-20) is an affinity purified goat polyclonal antibody raised against a peptide mapping at the C-terminus of Bicoid of *Drosophila melanogaster* origin.

## PRODUCT

Each vial contains 200  $\mu g$  lgG in 1.0 ml of PBS with < 0.1% sodium azide and 0.1% gelatin.

Blocking peptide available for competition studies, sc-33288 P, (100  $\mu$ g peptide in 0.5 ml PBS containing < 0.1% sodium azide and 0.2% BSA).

#### **STORAGE**

Store at 4° C, \*\*DO NOT FREEZE\*\*. Stable for one year from the date of shipment. Non-hazardous. No MSDS required.

#### APPLICATIONS

Bicoid (dD-20) is recommended for detection of Bicoid of *Drosophila melanogaster* origin by Western Blotting (starting dilution 1:200, dilution range 1:100-1:1000), immunofluorescence (starting dilution 1:50, dilution range 1:50-1:500) and solid phase ELISA (starting dilution 1:30, dilution range 1:30-1:3000).

#### **RECOMMENDED SECONDARY REAGENTS**

To ensure optimal results, the following support (secondary) reagents are recommended: 1) Western Blotting: use donkey anti-goat IgG-HRP: sc-2020 (dilution range: 1:2000-1:100,000) or Cruz Marker™ compatible donkey anti-goat IgG-HRP: sc-2033 (dilution range: 1:2000-1:5000), Cruz Marker™ Molecular Weight Standards: sc-2035, TBS Blotto A Blocking Reagent: sc-2333 and Western Blotting Luminol Reagent: sc-2048. 2) Immunofluores-cence: use donkey anti-goat IgG-FITC: sc-2024 (dilution range: 1:100-1:400) or donkey anti-goat IgG-TR: sc-2783 (dilution range: 1:100-1:400) with UltraCruz™ Mounting Medium: sc-24941.

## **RESEARCH USE**

For research use only, not for use in diagnostic procedures.

## PROTOCOLS

See our web site at www.scbt.com or our catalog for detailed protocols and support products.